MICROSTAAR® INJECTORS

MICROSTAAR® Injection System MODELS MSI-TF, MSI-PF, SFC-45, FTP

FOR USE WITH STAAR IMPLANTABLE COLLAMER® LENS
(ICL AND TICL)

CAUTION: Federal (USA) law restricts the sale of this device to, or upon the order of, a physician.

DIRECTIONS FOR USE

TABLE OF CONTENTS


EN (English)	1	IT (Italian)	41
BG (Bulgarian)	5	LT (Lithuanian)	44
CS (Czech)	8	LV (Latvian)	47
DA (Danish)	11	NL (Dutch)	50
DE (German)	14	NO (Norwegian)	53
EL (Greek)	17	PL (Polish)	56
ES (Spanish)	20	PT-PT (Portuguese)	59
ET (Estonian)	23	PT-BR (Brazilian Portuguese)	62
FI (Finnish)	26	RO (Romanian)	65
FR (French)	29	SK (Slovak)	68
HR (Croatian)	32	SL (Slovenian)	71
HU (Hungarian)	35	SV (Swedish)	74
IS (Icelandic)	38	TR (Turkish)	77

DEVICE DESCRIPTION

The MICROSTAAR Injection System, models MSI-TF and MSI-PF (Injectors), SFC-45 Cartridge (Cartridge), Foam Tip Plunger (FTP) and Foam Tip Plunger holder (FTP holder), are devices used to implant the STAAR Implantable Collamer® Lenses (ICLs and/or TICLs). The injection system provides a tubular pathway through an incision over the iris, to the center of the pupil. The ICL or TICL may then be gently injected from the tube into the eye.

The MICROSTAAR injection system consists of four parts:

1. The sterile reusable Injector Threaded (Fig. 1), or Plunger (Fig. 2).
2. Sterile, single-use, disposable SFC-45 Cartridge (Cartridge) (Fig. 3; do not resterilize).
3. Sterile, single-use, disposable Foam Tip Plunger (FTP) (Fig. 4; do not resterilize).
4. Sterile, single-use, disposable Foam Tip Plunger holder (FTP holder) (Fig. 5; do not resterilize).

INDICATION FOR USE

The MICROSTAAR Injection System is indicated for folding and inserting STAAR ICLs and TICLs for surgical placement in the eye.

CONTRAINDICATIONS


Use of the MICROSTAAR injection system components with other intraocular implants, Cartridges, or FTP has not been tested or approved for use by STAAR Surgical Company.

INSTRUCTIONS FOR USE

CAUTION: Do not use if unit package has been opened or damaged.

Prepare the FTP and the Injector

1. Open the pouch for the Injector and the outer pouch for the FTP and transfer contents to the sterile field.
2. Open the inner pouch and remove the FTP and holder. The base of the FTP protrudes out the back of the holder.
3. Insert the FTP into the Injector, base first, (the vertical tab of holder is not intended to be snap locked into notch of Injector). While holding the tab in place, advance the Injector cap until the ball end of the FTP interlocks with the Injector. A click can be felt and heard when the plunger is properly secured. (Fig. 6)
4. Retract the Injector plunger fully. The FTP will remain locked in place.
5. Remove the FTP holder by sliding it back out of the front of the Injector. (Fig. 7)
6. Place the assembled FTP and Injector into Balanced Salt Solution (BSS) such that the tip of the Injector is submerged up to the FTP. (Fig. 8)

CAUTION: Hydration of the foam tip is extremely important as insufficient hydration may cause trapping of the ICL or TICL between the plunger and Cartridge upon ICL or TICL delivery.

Load ICL/TICL into the Cartridge

CAUTION: Load the ICL or TICL immediately prior to injection. ICLs and TICLs should not be left folded in Cartridge for more than 1-2 minutes.


1. Open the pouch of the Cartridge package and transfer the tray to the sterile field. When ready to use, open the tray and remove the Cartridge.

CAUTION: Do not handle the Cartridge at the tip.

CAUTION: The loading of the ICLs and TICLs into the cartridge is a critical component of the overall procedure and should be performed under an operating microscope. The delivery of the implant into the anterior chamber is largely dependent on

the precise and careful loading of the ICLs and TICLs into the cartridge.


2. Fill the Cartridge first with BSS followed by a partial fill with approximately 0.3 mL methylcellulose type viscoelastic material, creating a trail out the back of the loading area. (Fig. 9)

3. Gently remove the ICL or TICL from the vial.

CAUTION: Do not use sharp or rough edge tools; they could damage the ICL or TICL. The ICL and TICL must remain hydrated.

4. Hold the ICL or TICL with forceps and inspect it under the operating microscope to identify and verify the correct orientation of the ICL or TICL. (Fig. 10)


It is important that the long jaws of the forceps avoid the optic area of the ICL or TICL to avoid damage to this thin and critical area of the ICL or TICL. The footplate marks on the leading right and trailing left haptic indicate that the correct (convex) side of the ICL or TICL is facing up. The other center marks are used to guide alignment during loading of the ICL or TICL into the Cartridge.

5. With the Cartridge held in one hand and the forceps holding the ICL or TICL in the other, the ICL or TICL is placed in the Cartridge loading chamber. Place the ICL or TICL into the Cartridge loading chamber so that the long axis of the ICL or TICL is positioned into the groove, under each side-rail of the Cartridge. This usually requires starting one long edge of the ICL or TICL under one rail and rolling the wrist to position the opposite side of the ICL or TICL under the opposing rail.

NOTE: The ICL or TICL must be convex in the Cartridge with both edges below the side rails of the Cartridge. (Fig. 11)

6. Close the jaws of the ICL and TICL loading forceps and insert them into the barrel, from the front of the Cartridge. (Fig. 12)
7. Advance the forceps through Cartridge until the jaws are about to contact the leading edge of the ICL or TICL. (Fig. 13)
8. Open the jaws of the forceps and grasp the footplate of the ICL or TICL, so that the ICL or TICL positioning mark is aligned with the jaws. (Fig. 14)
9. Slowly pull the ICL or TICL into the barrel while moving the Cartridge in the opposite direction. Observe the ICL or TICL positioning marks on either side of the optic to confirm alignment as you advance the ICL or TICL. Continue this process until the ICL or TICL is positioned within the Cartridge so that its leading edge is adjacent to the end of the Cartridge. Release the ICL or TICL and remove the forceps. (Fig. 15)


10. Any air bubbles should be evacuated from the loaded Cartridge with BSS injected into the tip of the Cartridge using the cannula to backfill the Cartridge.


Loading Cartridge into Injector and delivery of ICLs and TICLs

CAUTION: The ICLs and TICLs should be injected into the eye within 1-2 minutes after loading into the injector. Viscoelastic materials tend to lose their lubricity if exposed to the air too long.


1. Slide the Cartridge into the front of the Injector and snap lock the vertical tab into position. (Fig. 16 & 17)

2. Advance the plunger/FTP until it is in contact with the ICL or TICL. The final ICL or TICL position should be within approximately 2mm of the end of the Cartridge.
3. Inspect the orientation of the lens under the operating microscope. The clear funnel of the Cartridge enables identification of the center marks on each side of the optic of the lens. These marks should be visible at the 12 o'clock position and be in straight alignment down the shaft of the cartridge. (If there is not alignment or proper orientation, the ICL or TICL may be twisted and should be injected into the Cartridge tray and the loading process repeated). (Fig. 18)

- The fully assembled injector is placed tip down into a container of BSS to maintain lens hydration. (Fig.19)

CAUTION: The maximum recommended time for the loaded ICLs or TICLs to remain in the Injector prior to surgery is 1-2 minutes.

- Press (for PF Injector) or twist (for TF Injector) the cap/push-rod, and carefully advance the plunger/FTP forward in a smooth and slow motion until the initial contact with the ICL or TICL. Retract the plunger slightly and advance again. This step prevents the ICL or TICL from being overridden.
- Insert the tip of the Cartridge through the incision and into the Anterior Chamber of the eye.
- Rotate the injector to keep the center marks on the lens up and to allow the footplate of the ICL or TICL to unfold.
- If the ICL or TICL does not open, retract the plunger, which will allow the ICL or TICL to open.
- Inspect the orientation of the ICL or TICL after injection.
- When the ICL or TICL is completely released, remove the Injector and Cartridge from the incision.
- Manipulate the ICL or TICL into position behind the iris (Follow instructions as presented in the ICL or TICL training course).

Removal of the Cartridge and the Foam Tip Plunger

- Use the snap lock to unlock the Cartridge and then remove the Cartridge from the injector.
- Push or twist the cap to advance the foam tip plunger until it stops.
- Hold the foam tip plunger tightly. Retract the cap/push rod until the foam tip plunger separates from the injector.

NOTE: MICROSTAAR Cartridges, FTPs, and FTP holders are packaged and sterilized for single use only. MICROSTAAR Injectors are multi-use. The Injectors have been validated for cleaning and re-sterilization. (Refer to cleaning and sterilization instruction section.)

CLEANING INSTRUCTIONS for MICROSTAAR Injectors ONLY

Immediately after use, the Injector should be safely transported to the designated processing area and cleaned per the following instructions after removal of the cartridge and FTP.

- Wipe the nose cone with a non-linting wipe for a minimum of 15 seconds. Rinse the nose cone with deionized water at ambient temperature for a minimum of 15 seconds. Repeat wiping with a new non-linting wipe for a minimum of 15 seconds.
- For the MSI-TF only, loosen the plunger assembly by unthreading the plunger cap completely from the injector body prior to proceeding with the remaining cleaning.

- Fully immerse the Injector, nose cone downward, into a clean glass container of deionized water at ambient temperature for 1 minute. While immersed, use a soft bristled brush to gently clean the nose cone and plunger for the full immersion time.
- Remove the Injector allowing it to drain; then fully immerse the Injector, nose cone down, into a clean glass container of deionized water at ambient temperature. Sonicate the container at ambient temperature for a minimum of five minutes.
- Vertically remove the Injector and allow it to thoroughly drain, then fully immerse it, nose cone down, into a clean glass container of 70% isopropyl alcohol (IPA) and sonicate for a minimum of five minutes.
- Vertically remove the Injector allowing it to drain; then fully immerse the Injector, nose cone down, into a clean glass container of deionized water at ambient temperature and agitate the Injector while immersed for a minimum of 1 minute.
- Dry the Injector using non-linting wipes until visually dry. Use a syringe filled with air to dry the inside of the Injector until no moisture is visible.
- Visually inspect the Injector with adequate magnification and light to verify any soil has been removed. Re-clean if soil is observed. Failure to thoroughly clean the device could lead to inadequate sterilization. Verify that all parts are moving smoothly. Inspect the Injector for any cracks or other damage that may affect the functionality of the device. Discard Injector if any damage is detected.
- The MICROSTAAR Injector Models MSI-TF and MSI-PF can be cleaned and sterilized up to 20 times before disposal. The number of reprocessing cycles should be tracked and is dependent upon full compliance with these directions for use.

STERILIZATION INSTRUCTIONS for MICROSTAAR Injectors ONLY

After the Injector has been properly cleaned and inspected, the Injector should be placed into an appropriately sized, FDA cleared steam sterilization pouch and sterilized using either of the following steam sterilization cycles:

Cycle Type	Prevacuum	Gravity
Pulses	4	Not applicable
Set Point Temperature	132°C	132°C
Full Cycle Dwell Time	4 minutes	15 minutes
Dry Time	20 minutes	20 minutes

NOTE: The MICROSTAAR Injector Models MSI-TF and MSI-PF can be cleaned and sterilized up to 20 times before disposal. The number of reprocessing cycles should be tracked and is dependent upon full compliance with these directions for use.

WARNING

STAAR Surgical Cartridges, FTPs, and FTP holders are packaged and sterilized for single use only. The Cartridge, FTP, and FTP holder should not be resterilized, or reused. If one of these devices were reused after cleaning and/or resterilization, it is highly probable that it would be contaminated, and the contamination could result in endophthalmitis and inflammation.

Performance characteristics of Cartridges, FTPs, and FTP holders will degrade if they are reused.

SURGICAL PROCEDURE

The MICROSTAAR injection system is to be used by trained health care professionals for the surgical implantation of ICLs and TICLs into patients' eyes. Proper surgical technique is the responsibility of the individual surgeon. The surgeon must determine the suitability of any particular procedure based upon his or her medical training and experience. Dispose of any single use accessories that may have become contaminated with bodily fluids during the procedure as biohazardous waste according to standard surgical biohazard waste disposal procedure.

WARRANTY AND LIMITATION OF LIABILITY

STAAR Surgical Company warrants that reasonable care was taken in making this product.

STAAR Surgical Company shall not be responsible for any incidental or consequential loss, damage, or expense which arises directly or indirectly from the use of this product. Any liability shall be limited to the replacement of any MICROSTAAR Injector, Cartridge or FTP which is returned to and found to be defective by STAAR Surgical Company.

This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including but not limited to, any implied warranties or merchantability or fitness for use.

EXPIRATION DATE

The expiration date on the original Injector, Cartridge, and FTP packages is the sterility expiration date. The sterility is assured if the outer pouch seals are not punctured or damaged until the expiration date.

REPORTING

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as implant and Injector/Cartridge/FTP related and that were not previously expected in nature, severity and degree of incidence should be reported to STAAR Surgical. For surgeons/patient located in the EU, the competent authority should also be notified in the EU member state where the surgeon/patient is established.

USA/Canada Phone: +1 (800) 352-7842
International Phone: +(41) 32 332 8888

RETURN POLICY FOR STAAR INJECTION SYSTEM COMPONENTS

Contact STAAR Surgical.

STORAGE

Store the STAAR Injector, Cartridge, and FTP at room/ambient temperature.

HOW SUPPLIED

The MICROSTAAR Injectors, Cartridges, FTPs, and FTP holders are supplied in sealed sterile trays and/or pouch packaging. The packaging is sterilized by appropriate method and should be opened only under sterile conditions. MSI-TF and MSI-PF injectors are steam sterilized, SFC-45 Cartridges are sterilized using ethylene oxide, and FTP and FTP holders are sterilized using irradiation.

STAAR Surgical AG Spain
Avenida Meridiana 216,
Despacho 4
Barcelona 08027
Spain

STAAR Surgical AG
Hauptstrasse 104
CH-2560 Nidau
Switzerland
Tel: +(41) 32 332 8888

CE
0344

SYMBOL GLOSSARY

edfu.staar.com
(800) 352 7842
+41 32 332 88 88

Consult electronic instructions for use

Sterilized using steam

Sterilized using ethylene oxide

Sterilized using irradiation

Single sterile barrier system

Single sterile barrier system with
protective packaging inside

Catalog Number

Batch Code

Unique device identifier

Use by date

Date of manufacture

Do not re-use

Caution

Temperature limit – Store at
room/ambient temperature

Do not use if the product sterile barrier
system or its packaging is compromised

Manufacturer

Country of manufacture – United States

Authorized representative
in the European Community

Importer into the European Union

U.S. (Federal) law restricts this device
to sale by or on the order of a physician

Medical device

CE conformity marking per European
Council Directive 93/42/EEC or European
Council Regulation (EU) 2017/745

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