

EVO*
Visian ICL™

Evolution in Visual Freedom™

*EVO Visian ICL Family of lenses include EVO Visian ICL, EVO+ Visian ICL, EVO Visian Toric ICL, EVO+ Visian Toric ICL

PHAKIC INTRAOCULAR LENS

TABLE OF CONTENTS

PRODUCT INFORMATION	1	POST-APPROVAL STUDY OF THE EFFECT OF THE VISIAN MICL LENS ON AXIAL LENGTH MEASUREMENT	50
DEVICE DESCRIPTION	1	POST APPROVAL ADVERSE EVENT STUDY - VISIAN MICL LENS FOR MYOPIA	50
INDICATIONS	3	INSTRUCTIONS FOR USE	51
MODE OF ACTION	3	SURGICAL PRECAUTIONS/INFORMATION	52
CONTRAINDICATIONS	3	SURGICAL PROCEDURE	54
WARNINGS	4	PATIENT IMPLANT CARD	54
PRECAUTIONS	4	REPORTING	54
ADVERSE EVENTS	5	HOW SUPPLIED	54
CLINICAL TRIALS AND RESULTS	6	EXPIRATION DATE	54
POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT	15	RETURN POLICY FOR STAAR EVO/EVO+ ICL/TICL LENSES	54
POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL	22	LENS SPECIFIC RECOMMENDATION	54
PRE-APPROVAL VISIAN TICL LENS CLINICAL TRIAL AND RESULTS	27	WARRANTY AND LIMITATION OF LIABILITY	54
POST-APPROVAL STUDY OF THE VISIAN TICL LENS	32	STORAGE	54
ADDITIONAL CLINICAL OUTCOMES	39	SYMBOL GLOSSARY	55
PRE-APPROVAL PMA CLINICAL TRIAL AND RESULTS - VISIAN MICL LENS FOR MYOPIA	40		
POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL	44		



EVO|EVO+ VISIAN Implantable Collamer™ Lens (EVO|EVO+ ICL™) for Myopia

For the correction/reduction of moderate to high myopia

AND

EVO|EVO+ VISIAN TORIC Implantable Collamer™ Lens (EVO|EVO+ TICL™) for Myopia

For the correction/reduction of moderate to high myopic astigmatism

DIRECTIONS FOR USE

Manufactured and Distributed by

STAAR Surgical Company
1911 Walker Avenue
Monrovia, CA 91016
USA
Tel: (800) 352-7842
Fax: (800) 952-4923

CAUTION: U.S. (Federal) law restricts this device to sale by or on the order of a physician.

PRODUCT INFORMATION

Please review this product information completely before performing your initial clinical procedure. All physicians must complete the STAAR Surgical EVO/EVO+ ICL/TICL Lens Physician Certification Program prior to use.

DEVICE DESCRIPTION

The EVO ICL and EVO TICL lens (Implantable Collamer Lens) is an intraocular implant manufactured from Collamer, a proprietary hydroxyethyl methacrylate (HEMA)/porcine collagen containing biocompatible polymer material. The EVO/EVO+ ICL/TICL lens contains a UV absorber made from a UV absorbing material. The lens features a plate-haptic design with a central convex/concave optical zone and a 0.36 mm diameter central port; the lens incorporates a forward vault to minimize contact of the EVO/EVO+ ICL/TICL with the central anterior capsule.

While the parent devices (non-EVO/non-central port Visian MICL lens and Visian TICL lens) require preoperative peripheral iridotomies (PIs) to facilitate aqueous flow, the EVO/EVO+ ICL/TICL lenses include a central port that allows the flow of aqueous humor through the lens, thus eliminating the need for PIs prior to implantation.


The EVO/EVO+ ICL/TICL lenses feature an optic diameter that varies with the dioptric power; the smallest optic diameter being 4.9 mm and the largest 6.1 mm. The EVO/EVO+ ICL/TICL lenses are capable of being folded and inserted into the posterior chamber through an incision of 3.5 mm or less. The EVO/EVO+ ICL/TICL lenses are intended to be placed entirely within the posterior chamber directly behind the iris and in front of the anterior capsule of the human crystalline lens. When correctly positioned, the EVO/EVO+ ICL/TICL lenses function as a refractive element to optically reduce moderate to high myopia with or without astigmatism.

Table 1: EVO/EVO+ ICL Models

Brand Name	Model Name	Spherical Power (D)	Overall Diameter (mm)	Optic Diameter (mm)	Haptic Design
EVO	VICMO 12.1	-3.0 to -16.0	12.1	4.9 to 5.8	Flat, plate
EVO	VICMO 12.6	-3.0 to -16.0	12.6	4.9 to 5.8	Flat, plate
EVO	VICMO 13.2	-3.0 to -16.0	13.2	4.9 to 5.8	Flat, plate
EVO	VICMO 13.7	-3.0 to -16.0	13.7	4.9 to 5.8	Flat, plate
EVO+	VICM5 12.1	-3.0 to -16.0	12.1	5.0 to 6.1	Flat, plate
EVO+	VICM5 12.6	-3.0 to -16.0	12.6	5.0 to 6.1	Flat, plate
EVO+	VICM5 13.2	-3.0 to -16.0	13.2	5.0 to 6.1	Flat, plate
EVO+	VICM5 13.7	-3.0 to -16.0	13.7	5.0 to 6.1	Flat, plate

Table 2: EVO/EVO+ TICL Models

Brand Name	Model Name	Spherical Equivalent (D)	Cylindrical Power (D)	Overall Diameter (mm)	Optic Diameter (mm)	Haptic Design
EVO	VTICMO 12.1	-3.0 to -16.0	+1.0 to +4.0	12.1	4.9 to 5.8	Flat, plate
EVO	VTICMO 12.6	-3.0 to -16.0	+1.0 to +4.0	12.6	4.9 to 5.8	Flat, plate
EVO	VTICMO 13.2	-3.0 to -16.0	+1.0 to +4.0	13.2	4.9 to 5.8	Flat, plate
EVO	VTICMO 13.7	-3.0 to -16.0	+1.0 to +4.0	13.7	4.9 to 5.8	Flat, plate
EVO+	VTICM5 12.1	-3.0 to -16.0	+1.0 to +4.0	12.1	5.0 to 6.1	Flat, plate
EVO+	VTICM5 12.6	-3.0 to -16.0	+1.0 to +4.0	12.6	5.0 to 6.1	Flat, plate
EVO+	VTICM5 13.2	-3.0 to -16.0	+1.0 to +4.0	13.2	5.0 to 6.1	Flat, plate
EVO+	VTICM5 13.7	-3.0 to -16.0	+1.0 to +4.0	13.7	5.0 to 6.1	Flat, plate



*Artigas J. M., Felipe A., Navea A., et al. Spectral Transmission of the Human Crystalline Lens in Adult and Elderly Persons: Color and Total Transmission of Visible Light. Investigative Ophthalmology & Visual Science. 2012; 53 (7):4076- 4084.

Figure 1: Collamer ICL UV/Visible Spectrum




Figure 2: EVO/EVO+ ICL Lens Diagram

The EVO/EVO+ ICL lens has orientation markings on the footplates to ensure the lens is implanted right side up. When correctly oriented the orientation markings will be on the leading right/trailing left footplates.




Figure 3: EVO/EVO+ TICL Lens Diagram

The EVO/EVO+ TICL lens (**Figure 3**) is labeled using a plus cylinder axis format. The lenses are labeled to the nearest degree and as such lenses of any axis between 1° to 180° may be held in inventory. The EVO/EVO+ TICL lens is designed to be rotated up to 22.5° clockwise or counterclockwise in order to align the lens axis at the preoperative plus cylinder axis. The lens has two axis alignment markings, one on each side of the optic, these are to aid with the alignment of the lens. The markings indicate the meridian from which the cylinder axis is measured and do not indicate the cylinder axis of the lens.

The EVO/EVO+ TICL lens has orientation markings on the footplates to ensure the lens is implanted right side up. When correctly oriented the orientation markings will be on the leading right/trailing left footplates.

The sphere component of the EVO/EVO+ TICL lens label indicates the spherical power and not the spherical equivalent power.

INDICATIONS

The EVO/EVO+ ICL lens is indicated for use in patients 21-60 years of age:

1. for the correction of myopia with spherical equivalent ranging from -3.0 D to ≤ -15.0 D with less than or equal to 2.5 D of astigmatism at the spectacle plane;
2. for the reduction of myopia with spherical equivalent ranging from greater than -15.0 D to -20.0 D with less than or equal to 2.5 D of astigmatism at the spectacle plane;
3. with an anterior chamber depth (ACD) of 3.00 mm or greater, when measured from the corneal endothelium to the anterior surface of the crystalline lens, and a stable refractive history (within 0.5 D for 1 year prior to implantation).
4. The ICL lens is intended for placement in the posterior chamber (ciliary sulcus) of the phakic eye.

The EVO/EVO+ TICL lens is indicated for use in patients 21-60 years of age:

1. for the correction of myopic astigmatism with spherical equivalent ranging from -3.0 D to ≤ -15.0 D (in the spectacle plane) with cylinder (spectacle plane) of 1.0 D to 4.0 D.
2. for the reduction of myopic astigmatism with spherical equivalent ranging from greater than -15.0 D to -20.0 D (in the spectacle plane) with cylinder (spectacle plane) 1.0 D to 4.0 D.
3. with an anterior chamber depth (ACD) of 3.00 mm or greater, when measured from the corneal endothelium to the anterior surface of the crystalline lens and a stable refractive history (within 0.5 D for both spherical equivalent and cylinder for 1 year prior to implantation).
4. The TICL lens is intended for placement in the posterior chamber (ciliary sulcus) of the phakic eye.

MODE OF ACTION

The EVO/EVO+ ICL/TICL lenses function as a refractive element to optically reduce moderate to high myopia with or without astigmatism.

CONTRAINDICATIONS

The EVO/EVO+ ICL/TICL family of lenses is contraindicated in patients

1. with a true ACD of <3.00 mm*;
2. with anterior chamber angle less than Grade III as determined by gonioscopic examination;
3. who are pregnant or nursing;
4. less than 21 years of age;
5. who have moderate to severe glaucoma;
6. who do not meet the minimum endothelial cell density (ECD).

Table 3: Minimum Endothelial Cell Density for Age and True ACD*

Age	Minimum ECD ACD ≥ 3.0 mm	Minimum ECD ACD ≥ 3.2 mm	Minimum ECD ACD ≥ 3.5 mm
21-25	3875 cells/mm ²	3800 cells/mm ²	3250 cells/mm ²
26-30	3425 cells/mm ²	3375 cells/mm ²	2900 cells/mm ²
31-35	3025 cells/mm ²	2975 cells/mm ²	2625 cells/mm ²
36-40	2675 cells/mm ²	2625 cells/mm ²	2350 cells/mm ²
41-45	2350 cells/mm ²	2325 cells/mm ²	2100 cells/mm ²
>45	2075 cells/mm ²	2050 cells/mm ²	1900 cells/mm ²

* The true ACD is defined as the distance from the apex of the **posterior** corneal surface to the apex of the anterior crystalline lens surface. Many measuring devices provide an ACD measurement defined as the distance from the apex of the **anterior** corneal surface to the apex of the anterior crystalline lens surface. If the surgeon is using an instrument that measures from the anterior corneal surface, the thickness of the cornea must be subtracted to get the true ACD.

Table 3 indicates the minimum ECD per age group at time of implantation for three different ACD ranges. This data was developed as part of the STAAR ICL lens for Myopia Clinical Study (with the non-central port parent model ICL). This table was developed using rates of 2.47%, 2.44% and 2.15% (the upper 90% confidence interval of the average cell loss for eyes with the specified ACD) for the ≥ 3.0 mm, ≥ 3.2 mm, and ≥ 3.5 mm groups, respectively. It sets minimum ECD criteria as functions of age that should result in at least 1000 cells/mm² at 75 years of age. Specular microscopy should be performed preoperatively and ECD should be monitored postoperatively at intervals dictated by the physician's medical judgment.

WARNINGS

NOTE: All of these Warnings are applicable to the EVO/EVO+ ICL and EVO/EVO+ TICL

1. Some subjects in the STAAR Visian MICL lens for Myopia Clinical Study (with the non-central port parent model ICL) demonstrated endothelial cell loss >30% (range, 30.9% to 42.6%) at 5-7 years postoperatively. The long term effects (beyond 5 – 7 years) on the corneal endothelium have not been established. Additionally, some subjects in the EVO/EVO+ ICL/TICL Lens Study demonstrated endothelial cell loss > 30% at 3 years postoperatively. Patients should be advised about the potential risk of corneal edema, possibly requiring corneal transplantation. Patients' ECD should be monitored periodically as long as they remain implanted with the EVO/EVO+ ICL/TICL lens.
2. Secondary to implantation of the EVO/EVO+ ICL/TICL lens, patients have increased risk of development of cataract, including visually significant cataract that continues to increase with each year. The physician should monitor the patient for cataract periodically. The long term risk of visually significant cataract and related secondary surgery may be higher in older patients and those with higher myopia. The long-term rate (beyond 5-7 years) of cataract formation secondary to implantation, removal and/or replacement of the EVO/EVO+ ICL/TICL lens is unknown. In addition, prospective clinical trial data has not been collected in patients between 46-60 years of age. The surgeon should use clinical judgment to determine the benefit/risk ratio before implanting the ICL in older patients, with specific consideration to the potential development of cataracts.
3. Implantation of the EVO/EVO+ ICL/TICL lens is associated with an elevated risk of early postoperative increase in intraocular pressure (IOP). With the EVO/EVO+ ICL/TICL this is usually associated with incomplete removal of the OVD but could also be caused by angle closure (associated with pupillary block and/or excessive EVO/EVO+ ICL/TICL vault) that requires secondary surgical intervention. The risk of increased IOP due to incomplete removal of OVD can be mitigated by following the recommended OVD removal technique described briefly below (Intraoperative Information) and more fully in the EVO/EVO+ ICL/TICL Physician Certification Program. IOP should be initially checked 1 – 6 hours postoperatively, so that increased IOP can receive treatment as quickly as possible. The long-term risks of glaucoma, peripheral anterior synechiae and pigment dispersion are not well established.
4. Do not attempt to resterilize or repack the EVO/EVO+ ICL/TICL lens.
5. Do not autoclave the EVO/EVO+ ICL/TICL lens. Do not expose to temperature greater than 40°C. Do not freeze. If temperature requirements are not met, return the EVO/EVO+ ICL/TICL lens to STAAR Surgical.
6. The iridocorneal angle distance may decrease after implantation of the EVO/EVO+ ICL/TICL lens. Iridocorneal angle should be assessed 1 week after implantation and monitored if the angle is extremely narrow.
7. A patient with mesopic pupil size that is greater than the optic diameter of the EVO/EVO+ ICL/TICL lens may experience symptoms of glare and/or halos. Patients should be advised about this potential risk prior to EVO/EVO+ ICL/TICL lens implantation.
8. Complete removal of viscoelastic from the eye after completion of the surgical procedure is essential. STAAR Surgical recommends a low molecular weight 2% hydroxypropyl methylcellulose or dispersive, low viscosity ophthalmic viscosurgical device (OVD). Do not use short chain sodium hyaluronate acids (viscoelastics) due to increased risk of cataract formation related to trapped viscoelastic.

NOTE: The only viscoelastic used with the EVO/EVO+ ICL/TICL lens during the clinical trial was a low molecular weight 2% hydroxypropyl methylcellulose preparation.

PRECAUTIONS

Prior to surgery, the surgeon must provide prospective patients with a copy of the patient information booklet for this product and inform these patients of the possible benefits and complications associated with the use of this device.

NOTE: All of these Precautions are applicable to the EVO/EVO+ ICL and EVO/EVO+ TICL.

- Patients with higher degrees of myopia and/or myopic astigmatism experience lower efficacy and higher rates of adverse events (AEs) and complications.
- Inadequate flushing of the viscoelastic from the eye may lead to IOP spikes. IOP should be checked 1-6 hours postoperatively.
- The effectiveness of ultraviolet (UV) absorbing intraocular lenses (IOLs) in reducing the incidence of retinal disorders has not been established.
- The relationship between the EVO/EVO+ ICL/TICL lens and retinal detachment is undetermined.
- If a method of power calculation different from that used in the EVO/EVO+ ICL/TICL lens clinical study (i.e., lens power calculated by STAAR Surgical using STAAR's proprietary software) is used, the effectiveness of the EVO/EVO+ ICL/TICL lens for myopia with or without astigmatism may not be consistent with the results reported in the EVO/EVO+ ICL/TICL lens clinical study results section.
- The accuracy of ultra-sound based measurement of axial length in an eye with an EVO/EVO+ ICL/TICL lens is unknown. Axial length measurements based upon partial coherence laser interferometry appear to not be significantly affected by implantation of the ICL lens. See section on "Post-Approval Study of the Effect of the Visian MICL Lens on Axial Length Measurement."
- In the Visian TICL lens clinical study, surgeons were instructed to create one or two side port incisions, 60 – 90° away from the main incision, which should always be made at the horizontal temporal position. A 3.2 mm clear corneal tunnel incision was constructed parallel to the iris plane, with a tunnel length of 1.5 to 1.75 mm. If the surgeon uses a method of incision which is different from that used in the Visian TICL lens clinical study, the postoperative astigmatic results may not be consistent with the results reported for the Visian TICL lens clinical study, and the same precaution applies to implantation of EVO/EVO+ ICL and EVO/EVO+ TICL lenses. A temporal clear corneal tunnel incision of 3.5 mm or less constructed parallel to the iris plane, with a tunnel length of 1.5 to 1.75 mm, is recommended for implantation of EVO/EVO+ ICL and EVO/EVO+ TICL lenses.

The safety and effectiveness of the EVO/EVO+ ICL/TICL lens for the correction of moderate to high myopia has **NOT** been established in patients with

1. greater than 20.0 D of myopia;
2. greater than 2.5 D of astigmatism for the EVO/EVO+ ICL lens, or
3. astigmatism less than 1.0 D and greater than 4.0 D for the EVO/EVO+ TICL lens;
4. unstable or worsening myopia;
5. a diagnosis of ocular hypertension or glaucoma;
6. pseudoexfoliation;
7. pigment dispersion;
8. history or clinical signs of iritis/uveitis;
9. insulin-dependent diabetes or diabetic retinopathy;
10. history of previous ocular surgery;
11. progressive sight-threatening disease other than myopia;
12. serious (life-threatening) non-ophthalmic disease.

ADVERSE EVENTS

A list of adverse events associated with the EVO/EVO+ ICL and EVO/EVO+ TICL is provided below. Additionally, the location for specific adverse event data from the EVO/EVO+ ICL/TICL, Visian TICL and Visian MICL clinical studies is provided. For some events, the greatest detail is provided in the section that includes the adverse event data from the Visian MICL clinical studies (pre-approval study and extended follow-up post-approval study with the non-central port parent model ICL).

Table 4: Adverse Events

Adverse Event	For more information please refer to:
Implantation of the EVO/EVO+ ICL/TICL can be associated with insufficient EVO/EVO+ ICL/TICL vaulting over the crystalline lens, which can lead to anterior subcapsular opacities or clinically significant cataracts	<p>EVO/EVO+ ICL/TICL LENS PMA CLINICAL TRIAL AND RESULTS: Adverse Events</p> <p>POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events and Additional Safety Outcomes</p> <p>PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events</p> <p>PRE-APPROVAL PMA CLINICAL TRIAL AND RESULTS – VISIAN MICL LENS FOR MYOPIA: Adverse Events</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Adverse Events and Complications – Lens Opacity and Visually Significant Cataract Formation</p>
Implantation of the EVO/EVO+ ICL/TICL can be associated with excessive EVO/EVO+ ICL/TICL vaulting, which can cause a narrowing of the anterior chamber angle, possible pupillary block, increased intraocular pressure and glaucoma	<p>Adverse Events</p> <p>POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events and Additional Safety Outcomes</p> <p>PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL:</p> <ul style="list-style-type: none"> • Adverse Events • Surgical Reinterventions • Intraocular Pressure
<p>Implantation of the EVO/EVO+ ICL/TICL can be associated with substantial postoperative IOP increases/spikes, which can be associated with causes such as:</p> <ul style="list-style-type: none"> • retained OVD requiring interventions such as repeated irrigation, paracentesis / anterior chamber tap/aqueous tap, and hypotensive medication; • narrowed anterior chamber angle with or without pupillary block, requiring interventions such as iridotomy and/or hypotensive medication; and • pigmentary dispersion which can cause pigmentary glaucoma. 	<p>EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events</p> <p>POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events</p> <p>POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL</p> <p>PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events</p> <p>PRE-APPROVAL PMA CLINICAL TRIAL AND RESULTS – VISIAN MICL LENS FOR MYOPIA: Adverse Events</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Adverse Events and Complications</p> <p>POST-APPROVAL ADVERSE EVENT STUDY – VISIAN MICL LENS FOR MYOPIA</p>
Implantation of the EVO/EVO+ ICL/TICL is associated with an increased rate of chronic corneal endothelial cell loss, which may, over a period of time, lead to corneal edema and possibly the need for a corneal transplant	<p>EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events</p> <p>EVO/EVO+ ICL/TICL LENS Follow-Up Study: ECD Loss through Year 3 and Adverse Events</p> <p>PRE-APPROVAL VISIAN TICL CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL:</p> <ul style="list-style-type: none"> • Adverse Events • Endothelial Cell Density
EVO/EVO+ ICL/TICL may move out of its appropriate position	<p>EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events</p> <p>PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Vision TICL Related Additional Surgery</p> <p>POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Surgical Reinterventions</p>
There may be a need for secondary surgery for EVO/EVO+ ICL/TICL removal, replacement, or repositioning	<p>EVO/EVO+ ICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events</p> <p>POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL: Adverse Events</p> <p>PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Vision TICL Related Additional Surgery</p> <p>POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Surgical Reinterventions</p>
There may be a need for other types of secondary surgical intervention to treat some adverse events	<p>EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events</p> <p>POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL: Adverse Events</p> <p>PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications</p> <p>POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events</p> <p>POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Surgical Reinterventions</p>

Table 4: Adverse Events

Adverse Event	For more information please refer to:
There may be a loss of best spectacle-corrected visual acuity	EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events and Additional Safety Outcomes PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Best Corrected Distance Visual Acuity (CDVA) Loss
Implantation of the EVO/EVO+ ICL/TICL may cause an increase in refractive astigmatism	EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Refractive Cylinder (Target Variance) Distribution POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events
The EVO/EVO+ ICL/TICL may be associated with pigment dispersion and iris transillumination defects	EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Slit Lamp Findings
As with implantation of other types of intraocular lenses, potential adverse events can include, but are not limited to infection (endophthalmitis), hypopyon, corneal endothelial damage, IOL dislocation, cystoid macular edema, corneal edema, pupillary block, iritis, retinal detachment, retinal tear, transient or persistent glaucoma, vitritis, iris prolapse, secondary surgical intervention and increased visual symptoms related to the optical characteristics of the IOL including halos, glare and/or double vision	EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL: Adverse Events PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Optical Visual Symptoms PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Subjective Symptoms Stratified by Optic Diameter POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Surgical Reinterventions
Secondary surgical interventions may include, but are not limited to lens repositioning, lens replacement, vitreous aspiration, iridotomy/iridectomy for pupillary block, wound leak repair, retinal detachment repair and corneal transplantation	EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS: Adverse Events and Complications POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT: Adverse Events POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL: Adverse Events PRE-APPROVAL VISIAN TICL PMA CLINICAL TRIAL AND RESULTS: Vision TICL Related Additional Surgery POST-APPROVAL STUDY OF THE VISIAN TICL: Adverse Events POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL: Surgical Reinterventions; Other Complications

CLINICAL TRIALS AND RESULTS

Data from clinical studies of the EVO/EVO+ ICL/TICL lens and data from prior clinical studies of the parent Visian MICL lens and Visian TICL lens are included to support the safety and effectiveness of the EVO/EVO+ ICL/TICL lens. These include the following:

A PMA clinical study of the EVO/EVO+ ICL/TICL lenses was conducted to demonstrate the safety and effectiveness of the modifications of the previously approved parent Visian MICL and Visian TICL lenses, including the addition of a small central hole designed to permit aqueous flow and eliminate the requirement for laser peripheral iridotomy to prevent pupillary block. Post-approval follow-up of the PMA cohort through 3 years was required to evaluate the continued safety and effectiveness of the EVO/EVO+ ICL/TICL.

In the PMA clinical study of the EVO/EVO+ ICL/TICL lenses described above, 19.9% (125/629) of treated eyes experienced postoperative IOP spikes at 1-6 hours after surgery, believed to be due to incomplete OVD removal. Of note, no prophylactic IOP lowering medication was permitted in this study. A modification of the physician training program including instructions on appropriate OVD removal procedures was instituted to mitigate this risk. In addition, a new enrollment, post-approval study of the EVO/EVO+ ICL/TICL with two-week follow-up was performed to assess the effectiveness of this mitigation.

A prospective, nonrandomized clinical study of the parent Visian TICL lens in 210 eyes of 124 subjects was performed to demonstrate the safety and effectiveness of the modification of the previously approved Visian MICL lens model by the addition of a toric optic. A new enrollment, post-approval study of the parent Visian TICL lens with a two-year follow-up was performed to evaluate the long-term clinical performance of the lens with respect to rotational stability, refractive and visual outcomes, and ocular adverse event (AE) rates. Clinical studies of the parent Visian MICL lens including the primary safety and effectiveness study in 526 eyes of 294 subjects and three post-approval studies were performed: (1) extended follow-up of the pre-approval cohort to further characterize safety; (2) a new enrollment patient survey study to collect safety information from patients, and; (3) a post-approval study to assess the effect of the Visian MICL lens on axial length measurement.

The EVO/EVO+ ICL/TICL lens was evaluated in a prospective nonrandomized clinical study in 629 eyes of 327 subjects and two post-approval studies; extended follow-up of the pre-approval cohort to evaluate continued safety and effectiveness, and a prospective nonrandomized post-approval IOP study of 408 eyes of 205 subjects. The following sections provide details about each of these clinical studies.

PRE-APPROVAL EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS

This section includes clinical data on the EVO/EVO+ ICL/TICL lenses from a U.S. clinical study of these lenses.

The EVO and EVO+ sphere and toric Visian ICL (EVO/EVO+ ICL/TICL) lenses were evaluated through 6 months postoperative in a prospective nonrandomized single arm, three year study enrolling 629 eyes of 327 subjects. The primary analysis for the study was to occur after a minimum of 300 primary eyes completed the Month 6 Visit. Subject follow-up through 3 years was planned to obtain long-term data on clinical performance. The purpose of the study was to evaluate the safety, and to collect supportive data concerning the effectiveness of the EVO/EVO+ ICL/TICL lenses. Study subjects with moderate-to high myopia ranging from -3.00 to -20.00 D spherical equivalent (SE) in the spectacle plane or moderate to high myopic astigmatism with SE ranging from -3.00 to -20.00 D (in the spectacle plane) and cylinder ranging from 1.00 D to 4.00 D of cylinder (in the spectacle plane), with preoperative best spectacle corrected visual acuity (CDVA) of 20/40 or better and no pre-existing progressive sight-threatening ocular disorders other than pathological refractive error were eligible for the study.

The primary study (safety) endpoints were evaluated in primary (first eye treated) eyes only:

- Incidence of peripheral iridotomy (PI) required to treat elevated IOP caused by mechanical pupillary block through Month 6 Visit.
- Distribution of percent ECD losses and the percent of eyes that had ECD <1500 cells/mm² and ECD < 1000 cells/mm² through Month 6 visit (no prespecified performance target).
- Incidence of AEs through Month 6 Visit (no prespecified performance target).

Secondary (safety) endpoints were evaluated in all eyes (primary and fellow eyes) and had no prespecified performance targets:

- Incidence of PI required to treat elevated IOP caused by mechanical pupillary block through Month 6 Visit
- Distribution of percent ECD losses and the percent of eyes that had ECD <1500 cells/mm² and ECD <1000 cells/mm² through Month 6 Visit
- Incidence of AEs through Month 6 Visit

Effectiveness endpoints for this study had no prespecified performance targets:

- MRSE within ± 0.50 D and ± 1.00 D of target at Month 6 Visit
- UDVA of 20/40 or better at Month 6 Visit (for those eyes with CDVA 20/20 or better at Preoperative/Screening Visit)
- CDVA through Year 3 Visit (Day 1050 – 1170)

Demographics of the Study Cohort are presented in **Table 5**.

Table 5: Demographics

Demographics	Subjects (N=327) n (%)
Gender	
Male	114 (34.9)
Female	213 (65.1)
Race	
Caucasian	274 (83.8)
African American/Black	11 (3.4)
Asian	38 (11.6)
Native Hawaiian or Other Pacific Islander	3 (0.9)
American Indian or Alaska Native	1 (0.3)
Ethnicity	
Hispanic or Latino	34 (10.4)
Not Hispanic or Latino	293 (89.6)
Age (years)	
Mean (SD)	35.6 (5.1)
Median	36.0
Min, Max	22, 45

*Percentage calculated as (n/N)*100.

Accountability

A total of 327 patients (327 primary and 302 fellow eyes, 629 total eyes) were enrolled and underwent EVO/EVO+ ICL/TICL implantation in this study. One subject was discontinued from the study following lens explantation due to complaint of glare and halos. The interim analysis for PMA P030016/S035 included 303 primary eyes and 266 fellow eyes (569 total eyes) that completed the Month 6 visit. An update of safety data was submitted after all remaining treated eyes completed the Month 6 visit. Therefore, the safety data provided below includes all implanted eyes. Effectiveness data were not updated after all eyes completed the Month 6 visit; effectiveness data presented in this document are based on the 303 primary eyes (569 total eyes) that completed the Month 6 visit included in the interim analysis. Refer to Effectiveness Findings in POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT for all eyes that completed the Month 6 Visit.

Table 6 provides accountability for primary eyes and **Table 7** provides accountability for all eyes treated in the study.

Table 6: Accountability – Primary Eyes

Eye Status	Total #	Op Visit (Day 0) n (%)	Postop V1 (Day 1) n (%)	Postop V2 (Day 5-9) n (%)	Postop V3 (Day 21-35) n (%)	Postop V4 (Day 70-98) n (%)	Postop V5 (Day 147-182) n (%)	Postop V6 (Day 330-420) n (%)
All eyes treated (N)	327							
Available for analysis		327 (100.0)	327 (100.0)	325 (99.4)	325 (99.4)	324 (99.1)	321 (98.2)	42 (12.8)
Missing eye/data								
Discontinued		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	1 (0.3)	1 (0.3)	1 (0.3)
Missing at scheduled visit but seen earlier/later ¹		0 (0.0)	0 (0.0)	0 (0.0)	9 (2.8)	2 (0.6)	7 (2.1)	1 (0.3)
Missing but accounted for ²		0 (0.0)	0 (0.0)	2 (0.6)	2 (0.6)	1 (0.3)	3 (0.9)	12 (3.7)
Lost to follow-up		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	1 (0.3)	2 (0.6)	2 (0.6)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	270 (82.6)
% Accountability ⁴		327/327 (100.0)	327/327 (100.0)	325/327 (99.4)	325/327 (99.4)	324/326 (99.4)	321/326 (98.5)	42/56 (75.0)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Missing but accounted for: represents the total number of eyes that missed the visit but have not been discontinued/lost to follow-up.

³ Active: represents the total number of eyes that have not reached the time associated with the visit. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ % Accountability = [Available for Analysis/(Treated-Discontinued-Active)].

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Table 7: Accountability – All Eyes

Eye Status	Total #	Op Visit (Day 0) n (%)	Postop V1 (Day 1) n (%)	Postop V2 (Day 5-9) n (%)	Postop V3 (Day 21-35) n (%)	Postop V4 (Day 70-98) n (%)	Postop V5 (Day 147-182) n (%)	Postop V6 (Day 330-420) n (%)
All eyes treated (N)	629							
Available for analysis		629 (100.0)	628 (99.8)	624 (99.2)	626 (99.5)	624 (99.2)	619 (98.4)	81 (12.9)
Missing eye/data								
Discontinued		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	1 (0.2)	1 (0.2)	1 (0.2)
Missing at scheduled visit but seen earlier/later ¹		0 (0.0)	0 (0.0)	1 (0.2)	16 (2.5)	2 (0.3)	13 (2.1)	1 (0.2)
Missing but accounted for ²		0 (0.0)	1 (0.2)	5 (0.8)	3 (0.5)	2 (0.3)	5 (0.8)	0 (0.0)
Lost to follow-up		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	2 (0.3)	4 (0.6)	4 (0.6)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	543 (86.3)
% Accountability ⁴		629/629 (100.0)	628/629 (99.8)	624/629 (99.2)	626/629 (99.5)	624/628 (99.4)	619/628 (98.6)	81/85 (95.3)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Missing but accounted for: represents the total number of eyes that missed the visit but have not been discontinued/lost to follow-up.

³ Active: represents the total number of eyes that have not reached the time associated with the visit. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ % Accountability = [Available for Analysis/(Treated-Discontinued-Active)].

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Safety Outcomes

Incidence of Peripheral Iridotomy (PI) Required to Treat Elevated IOP Caused by Mechanical Pupillary Block

No primary eyes (0/327, 0.0%) and no fellow eyes (0/302, 0.0%) experienced pupillary block, and no PIs were performed through Month 6.

Table 8: Incidence of PI Required to Treat Elevated IOP Caused by Mechanical Pupillary Block

	Primary Eyes (N=327)		All Eyes (N=629)	
	No. Primary Eyes n (%)	No. Events	All Eyes n (%)	No. Events
Required PI to treat elevated IOP through Month 6	0 (0.0%)	0	0 (0.0%)	0

* Percentage calculated as (n/N)*100

ECD Losses through Month 6

No instances of ECD <1500 or <1000 cells/mm² through Month 6 were reported in the study, as shown in **Table 9**. Mean ECD loss (SD) from baseline was 2.4% (4.3%) in primary eyes and 2.3% (4.0%) for all eyes at Month 6. The range of change in ECD from baseline was +6.3% to -46.7%, with 97.3% (602/619) of all eyes experiencing ≤10% ECD loss from preoperative values. Three eyes of 3 subjects (3/619, 0.5%) reported ECD loss > 30% which was related to the surgical procedure.

Table 9: ECD Change from Baseline Through Month 6

Parameter	Primary Eyes (N=321*)		All Eyes (N=619*)	
	Value	95% CI	Value	95% CI
% ECD Change from Baseline				
N - Missing	319	-	614	-
Mean (SD)	-2.4 (4.3)	-2.860, -1.922	-2.26 (4.01)	-2.576, -1.941
Median	-1.8	-	-1.68	-
Min, Max	-46.7, 6.3	-	-46.7, 6.3	-
Distribution of % ECD Change from Baseline	n (%)	95% CI	n (%)	95% CI
Gain > 5%	2 (0.6)	0.08, 2.23	2 (0.3)	0.04, 1.16
Gain ≥ 2% to ≤ 5%	9 (2.8)	1.29, 5.26	22 (3.6)	2.24, 5.33
Gain < 2% to Loss < 2%	161 (50.2)	44.55, 55.76	320 (51.7)	47.68, 55.70
Loss ≥ 2% to ≤ 5%	101 (31.5)	26.42, 36.85	190 (30.7)	27.08, 34.49
Loss > 5% to ≤ 10%	41 (12.8)	9.32, 16.93	68 (11.0)	8.63, 13.72
Loss > 10% to ≤ 20%	3 (0.9)	0.19, 2.71	8 (1.3)	0.56, 2.53
Loss > 20% to ≤ 30%	0 (0.0)	0.00, 1.14	1 (0.2)	0.00, 0.90
Loss > 30%	2 (0.6)	0.08, 2.23	3 (0.5)	0.10, 1.41
Missing	2	-	5	-
ECD less than 1500 (n, %)	0 (0.0)	0.00, 1.14	0 (0.0)	0.00, 0.59
ECD less than 1000 (n, %)	0 (0.0)	0.00, 1.14	0 (0.0)	0.00, 0.59

* N is the number of eyes present at both the Preoperative and Month 6 Visits
Percentage calculated as (n/N)*100.

Adverse Events

All ocular AEs (only eyes implanted with study lenses) and all serious AEs (both ocular and nonocular) were to be reported in this study. Non-serious non-ocular AEs were not reported. All secondary surgical interventions (SSIs) and the events that caused these interventions were required to be reported as SAEs.

Experience with intraocular surgery and the implantation of IOLs has shown that some events can be considered normal or expected after these procedures. Early, low grade anterior chamber cell/flare, corneal edema, and increase in IOP can often be considered normal or expected after phakic IOL surgery and were not to be reported as AEs if they occurred prior to 1 week postoperatively and if they met the following criteria:

- AC cells or flare of ≤ grade 2 (using the SUN criteria) that require no change in standard postoperative medication regimen
- Corneal edema of ≤ grade 2 that does not reduce CDVA to 20/40 or worse and does not require any change in standard postoperative medication regimen
- Increased IOP that is <10 mmHg above baseline or is <25 mmHg and requires no change in standard postoperative medication regimen or any other special treatment
- Loss of CDVA ≥10 letters up to 1 week postoperatively

All other untoward events that occurred during the study, and all events with sequelae were to be reported as AEs, regardless of when they occurred.

Adverse Events – EVO/EVO+ ICL/TICL Clinical Trial

A total of 203 ocular AEs were reported for 25.8% (162/629) of all implanted eyes (**Table 10**).

Ocular AEs reported in the EVO/EVO+ ICL/TICL PMA study through the update of safety data that was submitted after all eyes completed the Month 6 visit are provided in **Table 10**. Details on the ocular AEs that were categorized as serious are provided in **Table 11**. The incidence of cumulative and persistent ocular AEs identified in the ISO 11979-7:2018 historical grid for Primary (n=327) and All (n=629) eyes are presented in **Table 12**. The results of AE analyses based on the consensus definitions as set forth by American Academy of Ophthalmology's (AAO) Task Force (Masket et al, 2017) are provided in **Table 14**.

The most frequent AE observed in the EVO/EVO+ ICL/TICL clinical trial was increased IOP caused by retained OVD (19.9%, 125/629), steroid response (2.4%, 15/629) or secondary surgical intervention (0.5%, 3/629). Increased IOP is discussed in more detail in the next section.

Three eyes (3/629, 0.5%) of 3 subjects reported ECD losses of > 30% from baseline at the 6 Month visit that was related to the surgical procedure. No instances of ECD less than 1500 or 1000 cells/mm² through Month 6 have been reported for any eye in this study.

No anterior subcapsular opacities or anterior subcapsular cataracts were reported through Month 6. There was a single report of a nuclear sclerotic cataract (0.16%, 1/629).

Three eyes of 2 subjects experienced retinal events, for an overall incidence of 0.5% (3/629). Surgical intervention (4 SSIs of retinal laser in 2 eyes of 1 subject and pars plana vitrectomy in 2 eyes of 2 subjects) was performed to treat each of these events.

Two eyes (2/629, 0.3%) of 2 subjects experienced anterior chamber angle narrowing that required secondary surgical intervention (SSI). Both of these events resolved following an initial repositioning of the lens and subsequent lens exchange. Neither event was associated with increased IOP. One subject complained of halo and glare in 1 eye (1/629, 0.2%) which resolved following explantation of the lens, and 1 subject complained of blurred vision related to residual astigmatism in 1 eye (1/629, 0.2%) which resolved following rotational repositioning of the toric lens.

No significant persistent loss of CDVA greater than or equal to 2 lines (10 letters) was reported through Month 6; only one eye (1/629, 0.2%) experienced a transient loss of 2 lines (10 letters), which resolved by the next study visit.

Table 10: Cumulative Ocular Adverse Events

Cumulative Ocular AEs	Primary Eyes (N=327)		All Eyes (N=629)	
	Eyes ¹ n (%) ²	Events n	Eyes ¹ n (%) ²	Events n
Eyes experienced any ocular AE	90 (27.5)	108	162 (25.8)	203
Intraocular pressure increased ³	75 (22.9)	77	136 (21.6)	143
Anterior chamber cell/flare ⁴	7 (2.1)	7	11 (1.7)	11
Corneal epithelial defect	3 (0.9)	3	6 (1.0)	6
Narrow anterior chamber angle ⁵	2 (0.6)	3	2 (0.3)	3
Corneal endothelial cell loss ⁶	2 (0.6)	2	4 (0.6)	4
Dry eye	2 (0.6)	2	4 (0.6)	4
Intraocular lens exchange	2 (0.6)	2	2 (0.3)	2
Intraocular lens repositioning	2 (0.6)	2	3 (0.5)	3
Retinal surgery	1 (0.3)	1	3 (0.5)	7
Retinal detachment ⁷	1 (0.3)	1	3 (0.5)	3
Glaucoma	1 (0.3)	1	2 (0.3)	2
Contact dermatitis	1 (0.3)	1	2 (0.3)	2
Intraocular lens removal	1 (0.3)	1	1 (0.2)	1
Cataract nuclear	1 (0.3)	1	1 (0.2)	1
Glare/Halo ⁸	1 (0.3)	1	1 (0.2)	1
Hordeolum	1 (0.3)	1	1 (0.2)	1
Iris incarceration	1 (0.3)	1	1 (0.2)	1
Visual acuity reduced ⁹	1 (0.3)	1	1 (0.2)	1
Retinal tear	0 (0.0)	0	1 (0.2)	2
Vitreous detachment	0 (0.0)	0	2 (0.3)	2
Astigmatism ¹⁰	0 (0.0)	0	1 (0.2)	1
Eye discharge	0 (0.0)	0	1 (0.2)	1
Punctate keratitis	0 (0.0)	0	1 (0.2)	1

¹ Percentage calculated as (n/N)*100.

² Only the first incidence of an event is counted for any given eye.

³ IOP ≥ 10 mmHg above baseline to a minimum of 25 mmHg or that required a change in the standard postoperative medication regimen or other special treatment was reported as an AE.

⁴ Anterior chamber cell/flare was reported as an AE if it met criteria for chronic anterior uveitis or was greater than grade 2 at Visit 2 (Day 5 – 9) or later.

⁵ Only those cases in which the investigator observed a reduction in anterior chamber angle and believed that a Secondary Surgical Intervention (SSI) was necessary. See **Table 17** for more information on gonioscopic evaluation.

⁶ Cases of endothelial cell loss that were counted as AEs included only cases of loss >30%. Refer to ECD Losses Through Month 6 section for additional information.

⁷ Refer to **Table 89** for additional information (rates of retinal detachment in original FDA study of the Visian MCL).

⁸ Only glare/halo leading to lens explantation was reported as an AE.

⁹ Loss of CDVA ≥10 letters at any time point > 1 week postoperatively was reported as an AE. Refer to Visual Acuity section and Additional Safety Outcomes section of POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT for more detail on loss of CDVA.

¹⁰ Residual astigmatism requiring second surgery of lens rotational repositioning.

Table 11: Ocular SAEs – All Eyes

Cumulative Ocular SAEs	All Eyes (N=629)	
	Eyes ¹ n (%) [*]	Events n
Eye experienced any ocular SAE	7 (1.1)	22
Eye disorders		
Glare/Halo	1 (0.2)	1
Narrow anterior chamber angle	2 (0.3)	3
Retinal detachment	3 (0.5)	3
Retinal tear	1 (0.2)	2
Surgical Reinterventions		
Intraocular lens exchange	2 (0.3)	2
Intraocular lens removal	1 (0.2)	1
Intraocular lens repositioning	3 (0.5)	3
Retinal surgery	3 (0.5)	7

¹ Only the first incidence of an event is counted for any given eye.

^{*} Percentage calculated as (n/N)*100.

Table 12: Cumulative and Persistent Ocular AEs¹

Adverse Event	Primary Eyes ²	All Eyes
Cumulative	N=327 n, %³	N=629 n, %³
Cystoid Macular Edema	0, 0%	0, 0%
Hypopyon	0, 0%	0, 0%
Endophthalmitis	0, 0%	0, 0%
IOL Dislocation	0, 0%	0, 0%
Pupillary Block	0, 0%	0, 0%
Retinal Detachment ⁴	1, 0.3%	3, 0.5%
Secondary Surgical Intervention	6, 1.8%	9, 2.8%
Persistent⁵	N=321 n, %³	N=619 n, %³
Corneal Stroma Edema	0, 0%	0, 0%
Cystoid Macular Edema	0, 0%	0, 0%
Iritis	0, 0%	0, 0%
Raised IOP Requiring Treatment	0, 0%	0, 0%

¹ Refer to Table B.2 in ISO 11979-7 2018: Ophthalmic implants - Intraocular lenses Part 7: Clinical investigations for AE categories included in table.

² Only the first incidence of an event is counted for any given eye.

³ Percentage calculated as (n/N)*100.

⁴ Comparison should be made to literature for retinal detachment rates for high myopia. Retinal detachment rates increase with increasing myopia. Refer to **Table 89** for additional information (rates of retinal detachment in original FDA study of the Visian MCL).

⁵ Persistent events are those that are present at the Month 6 visit. N is the number of eyes available at the Month 6 Visit (321 primary eyes and 619 total eyes).

Table 13: Secondary Surgical Reinterventions

Surgical Reinterventions	All Eyes (N=629)	
	Eyes ¹ n (%) ²	Events n
Intraocular lens exchange	2 (0.3)	2
Intraocular lens removal	1 (0.2)	1
Intraocular lens repositioning	3 (0.5)	3
Retinal surgery	3 (0.5)	7

¹ Only the first incidence of an event is counted for any given eye.

² Percentage calculated as (n/N)*100.

Table 14: Supportive Characterization of Ocular Adverse Events based on a Modified Version of AAO Consensus¹

Adverse Event	Primary Eyes N=327 n, % ²	All Eyes N=629 n, % ²
Chronic Anterior Uveitis	0, 0%	0, 0%
Clinically Significant Cystoid Macular Edema ≥ 1 month	0, 0%	0, 0%
Corneal Edema ≥ 1 week	0, 0%	0, 0%
Endophthalmitis	0, 0%	0, 0%
Mechanical Pupillary Block	0, 0%	0, 0%
Increased IOP	75, 22.9%	136, 21.6%
Retinal Detachment	1, 0.3%	3, 0.5%
Toxic anterior segment syndrome	0, 0%	0, 0%
Hypopyon	0, 0%	0, 0%
IOL Dislocation	0, 0%	0, 0%
Secondary IOL intervention - Exchange	2, 0.6%	2, 0.3%
Secondary IOL intervention - Removal	1, 0.3%	1, 0.2%
Secondary IOL intervention - Reposition	2, 0.6%	3, 0.5%

¹ Masket S, Rorer E, Stark W, Holladay J, MacRae S, Tarver ME, Glasser A, Calogero D, Hilmantel G, Nguyen T, Eydelman M. Special Report: The American Academy of Ophthalmology Task Force Consensus Statement on Adverse Events with Intraocular Lenses. Ophthalmology. 2017;124: 142-144.

² Percentage calculated as (n/N)*100.

Increased Intraocular Pressure

Increased IOP was the most frequently reported AE in the study through at least Month 6. No instances of increased IOP were attributed by investigators to pupillary block, anterior chamber angle narrowing, pigment dispersion or intraocular inflammation. No prophylactic IOP lowering medications were allowed during the study. These AEs commonly occurred either at PO visit 0 (1 – 6 hours) due to retained OVD or 6 to 31 days postoperative due to steroid response. An increase in IOP with onset 1 – 6 hours postoperatively was reported for 19.9% (125/629) of treated eyes. These AEs, related to incomplete removal of the dispersive OVD at the end of the surgical procedure, were managed either without treatment or with aqueous tap and/or ocular hypotensive medication and all resolved without sequelae by the first postoperative day. **Table 15** provides the distribution of maximum IOP in these cases, and **Table 16** provides the numbers of eyes treated with aqueous tap and/or medication.

Table 15: Maximum IOP Among Incidences of Elevated IOP with Onset on Day 0

Adverse Event - Elevated IOP	Primary Eyes (N=327) n (%) ¹	All Eyes (N=629) n (%) ¹
Number of elevated IOP events	67 (20.5)	125 (19.9)
Maximum IOP (mmHg)		
< 30	17 (5.2)	40 (6.4)
≥ 30	50 (15.3)	85 (13.5)
≥ 40	23 (7.0)	38 (6.0)
≥ 50	13 (4.0)	24 (3.8)
≥ 60	6 (1.8)	11 (1.7)
≥ 70	0 (0.0)	1 (1.6)

¹ Percentage calculated as (n/N)*100.

Table 16: Elevated IOP Requiring Treatment with Onset on Day 0 (All Treated Eyes)

Number of elevated IOP events requiring treatment	Primary Eyes (N=55) n (%) ¹	All Eyes (N=97) n (%) ¹
Events treated with concomitant medication(s)	53 (96.4)	94 (96.9)
Events treated with paracentesis/ AC tap*	39 (70.9)	70 (72.2)

¹ Percentage calculated as (n/N)*100.

Note: "paracentesis/AC tap" refers to burping the existing corneal incision to release aqueous; in no case was a needle paracentesis performed.

Investigators were previously certified ICL surgeons (through required training) and had experience implanting the U.S.-approved Visian MICL/TICLs. The OVD used in the study was hydroxypropylmethylcellulose 2% (HPMC), the OVD recommended by STAAR, and training and labeling pointed out the importance of thorough removal of the OVD to reduce the risk of postoperative increases in IOP. Investigators provided responses to a questionnaire regarding their surgical techniques of OVD removal, including the thoroughness of removal and the volume of balanced salt solution (BSS) used for irrigation. Comparison of the questionnaire responses with the incidence of elevated IOP at the 1 – 6 hour postoperative visit demonstrated that the 2 surgeons who practiced the least thorough methods of OVD removal and used the least volume of BSS for irrigation accounted for all of the events of increased IOP ≥ 40 mmHg, and the 4 surgeons who practiced the least thorough methods of OVD removal accounted for 84.0% (105/125) of events of elevated IOP but only 38% of enrolled eyes. Conversely, the 7 surgeons reporting the most thorough methods of OVD removal accounted for 55% of enrolled eyes but only 13.6% (17/125) of events of elevated IOP. These results support that the thoroughness of OVD removal is related to the incidence of elevated IOP at the 1 – 6 hour postoperative visit.

An additional 15 events (15/629, 2.4%) of increased IOP with onset from 6 to 31 days postoperative were related to the use of a topical corticosteroid and resolved with continued steroid taper and/or topical ocular hypotensive medication. Increased IOP as a result of secondary surgical intervention was reported for an additional three eyes (3/629, 0.5%). None of these events was attributed by investigators to the EVO/EVO+ ICL/TICL lens, nor was any event attributed by investigators to either blockage of the flow of aqueous through the central port or narrowing of the anterior chamber angle.

Other Safety Outcomes

NOTE: For other safety outcomes (Gonioscopy, Loss of CDVA from baseline, and Vault), N is 569 eyes at the Month 6 Visit as these data are based on the treated eyes available for the interim analysis, prior to the safety update.

Gonioscopy

Table 17 provides the results of gonioscopy at baseline and Month 6. A total of 60 eyes (60/569, 10.5%) demonstrated a narrower angle at Month 6 than at the preoperative visit.

Table 17: Gonioscopy by Visit in All Eyes (Safety Population)

Gonioscopy	Preoperative Visit (N=629) n (%) ¹	Month 6 Visit (N=569) n (%) ¹
Angle grade		
0	0 (0.0)	0 (0.0)
1	0 (0.0)	2 (0.4)
2	0 (0.0)	9 (1.6)
3	66 (10.5)	87 (15.3)
4	563 (89.5)	469 (82.4)
Missing	0	2
Pigmentation grade		
0	497 (79.0)	430 (75.6)
1	102 (16.2)	110 (19.3)
2	13 (2.1)	11 (1.9)
3	17 (2.7)	16 (2.8)
4	0 (0.0)	0 (0.0)
Missing	0	2
Peripheral anterior synechiae		
Absent	628 (99.8)	566 (99.5)
Present (specify clock hours)	1 (0.2)	1 (0.2)
0.5-2.0	1 (0.2)	1 (0.2)
2.5-4.0	0 (0.0)	0 (0.0)
4.5-6.0	0 (0.0)	0 (0.0)
6.5-8.0	0 (0.0)	0 (0.0)
8.5-10.0	0 (0.0)	0 (0.0)
10.5-12.0	0 (0.0)	0 (0.0)
Missing	0	2

¹ Percentage calculated as (n/N)*100.

Loss of CDVA from Baseline

No significant persistent loss of CDVA ≥ 2 lines (10 letters) was reported in this study; only 1 eye experienced a transient loss of 2 lines at Week 1, which resolved by the next visit. Overall, 91.7% (522/569) of all eyes reported unchanged or increased CDVA at Month 6 compared with the preoperative visit.

Vault

Table 18 provides the number and percent of eyes with vault measurements <250 microns and >900 microns, as well as mean vault and quartiles for vault at the Month 6 visit. The preoperative factors showing the greatest correlation to achieved vault were crystalline lens rise above the ATA (angle to angle) plane and lens diameter (Figure 4 and Figure 5). Crystalline lens rise was defined as the distance between the anterior pole of the crystalline lens and a plane joining the following anatomical landmarks: the iridocorneal recess (angle to angle plane), the iris pigment end (pigment to pigment plane), or the ciliary sulcus (sulcus to sulcus plane).

Table 18: Lens Vault at Month 6 Visit (Interim Analysis)

Parameter	Primary Eyes	All Eyes
Number of eyes with vault measurement (N)	301	566
Number (%) of eyes measured with vault < 250 μ	33 (11.0)	69 (12.2)
Number (%) of eyes measured with vault > 900 μ	16 (5.3)	30 (5.3)
Mean vault (μ)	503.2	496.8
0 th percentile for measured vault (μ)	10.0	10.0
25 th percentile for measured vault (μ)	350.0	346.0
50 th percentile for measured vault (μ)	475.0	470.0
75 th percentile for measured vault (μ)	637.0	634.0
100 th percentile for measured vault (μ)	1240.0	1240.0

¹ Percentage calculated as (n/N)*100




Figure 4: Vault at Month 6 by Crystalline Lens Rise




Figure 5: Vault at Month 6 by Lens Diameter

Effectiveness Outcomes

Accuracy of Refractive Outcome

MRSE by visit is provided in Table 19. As shown in Table 20, 89.4% (271/303) and 98.3% (298/303) primary eyes and 90.5% (563/569) and 98.9% (563/569) of all eyes achieved MRSE within ± 0.5 D and ± 1.0 D from target at the 6 month examination, respectively.

Table 19: MRSE by Visit

MRSE (D)	PreOp	Month 1	Month 3	Month 6
Primary Eyes				
N	327	325	324	321
Mean (SD)	-7.63 (2.80)	-0.11 (0.29)	-0.05 (0.31)	-0.09 (0.38)
Median	-7.38	-0.120	0.000	0.000
Min, Max	-15.62, -3.00	-1.25, 1.00	-1.62, 1.12	-3.88, 1.12
Missing	0	0	0	0
All Eyes				
N	629	626	624	619
Mean (SD)	-7.62 (2.76)	-0.11 (0.30)	-0.03 (0.31)	-0.08 (0.33)
Median	-7.38	-0.120	0.000	0.000
Min, Max	-15.62, -3.00	-1.25, 1.00	-1.62, 1.12	-3.88, 1.12
Missing	0	0	0	0

Table 20: MRSE Within ± 0.50 D and ± 1.00 D of Target at Month 6

	Primary Eyes (N=321)		All Eyes (N=619)	
	n	Proportion (95% CI)	n	Proportion (95% CI)
± 0.50 D	286	0.891 (0.8516 - 0.9229)	560	0.905 (0.8788 - 0.9266)
± 1.0 D	316	0.984 (0.9640 - 0.9949)	613	0.990 (0.9790 - 0.9964)

Visual Acuity

The 6 Month postoperative results provided in **Table 21** and **Table 22** – **Table 23** demonstrate that the EVO/EVO+ ICL/TICL lens provides accurate refractive correction and levels of uncorrected distance visual acuity (UDVA) consistent with the non-central port Visian MICL and TICL parent lenses.

Table 21: UDVA at 6 Months (Where emmetropia was the goal (± 0.50 D) and Preoperative Best Corrected Visual Acuity (CDVA) was 20/20 or better)

	All Eyes
N (506)	n, % ¹
20/20 or better	407, 80.4%
20/40 or better	503, 99.4%

¹ Percentage calculated as (n/N)*100.

Table 22: Best Corrected Distance Visual Acuity (CDVA) at 6 Months (Eyes with Preoperative CDVA 20/20 or better)

	6 Months
N (506)	n, % ¹
20/20 or better	500, 98.8%
20/40 or better	506, 100%

¹ Percentage calculated as (n/N)*100.

Table 23: Best Corrected Distance Visual Acuity (CDVA) at 6 Months (All Eyes)

	6 Months
N (619)	n, % ¹
20/20 or better	599, 96.8%
20/40 or better	619, 100%

¹ Percentage calculated as (n/N)*100.

POST-APPROVAL FOLLOW-UP OF THE EVO/EVO+ ICL/TICL PMA COHORT

Study Objective

As a condition of FDA approval, extended follow-up of the PMA cohort through 3 years was required. The objective of this post-approval study was to evaluate the continued safety and effectiveness of the EVO/EVO+ ICL/TICL in the original PMA clinical study cohort through 3 years of follow-up after surgery.

Study Design

The PMA study cohort was seen for 3 scheduled visits at Year 1 (Days 330–420), Year 2 (Days 690–810), and Year 3 (Days 1050–1170) after ICL surgery. Assessments performed at these visits included UDVA, CDVA, manifest refraction, gonioscopy, slit lamp examination, crystalline lens status, specular microscopy (corneal endothelial cell count), IOP, dilated fundus examination, optical coherence tomography (OCT; lens vault), and collection of AEs.

Study Endpoints

The following co-primary endpoints were evaluated in all eyes and had no prespecified performance targets.

- Distribution of percent ECD losses and the percent of eyes that had ECD <1500 cells/mm² and ECD <1000 cells/mm² through Year 3 Visit (Day 1050–1170)
- Incidence of AEs through Year 3 Visit (Day 1050–1170)

Additional safety and effectiveness parameters assessed in the original PMA study were evaluated through the Year 3 Visit. Key findings from these assessments are provided in this section.

Study Population and Data Source

This study enrolled 324 subjects of the 327 subjects who completed the Month 6 visit in the PMA study. At the Year 3 Visit, the 629 eyes of 327 subjects enrolled at 14 US sites in the original PMA study were accounted for, and 579 eyes (93.4% Accountability) were available for analysis (**Table 24**).

Table 24: Accountability – Post-Approval Follow-Up of the EVO/EVO+ ICL/TICL PMA Cohort

Eye Status	Total N	Year 1 (Days 330–420) n (%)	Year 2 (Days 690–810) n (%)	Year 3 (Day 1050–1170) n (%)
All eyes treated (N)	629			
Available for analysis		615 (97.8)	585 (93.0)	579 (92.1)
Missing eye/data				
Discontinued		1 (0.2)	4 (0.6)	9 (1.4)
Missing at scheduled visit but seen earlier/ later ¹		8 (1.3)	23 (3.7)	3 (0.5)
Missing but accounted for ²		0 (0)	7 (1.1)	0 (0.0)
Lost to follow-up		13 (2.1)	33 (5.2)	41 (6.5)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)
% Accountability ⁴		615/628 (97.9)	585/625 (93.6)	579/620 (93.4)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Missing but accounted for: represents the total number of eyes that missed the visit but have not been discontinued/lost to follow-up.

³ Active: represents the total number of eyes that have not reached the time associated with the visit. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ % Accountability = [Available for Analysis/(Treated-Discontinued-Active)].

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Safety Findings

ECD Loss through Year 3

No instances of ECD <1500 or <1000 cells/mm² were reported in the study, as shown in **Table 25**. The mean (SD) ECD loss from the preoperative visit at Year 3 for all eyes was 6.7% (5.4%). All instances of ECD losses of ≥20% compared to preoperative baseline values were reported as AEs. ECD loss > 30% from baseline was reported in 4 eyes of 4 subjects during the study (4/579; 0.7%). In three eyes, the decrease in ECD occurred in the first 6 months after ICL surgery, was attributed to the surgical procedure and not the study device and stabilized following the initial postoperative decline. In the fourth eye, the decrease was reported at the Year 3 Visit. A >20% to ≤30% ECD loss from baseline was reported for 9 eyes of 9 subjects at Year 3 (9/579; 1.6%) at Year 3.

Per the study protocol, the investigator was required to consult with the Medical Monitor to determine if implant removal might be warranted. In each case of ECD loss >20%, the investigator decided that implant removal was not warranted, and the subject would continue to be followed closely with serial specular microscopy per the investigator's standard of care.

Table 25: ECD Change from Baseline at Each Visit

Parameter	Year 1 (N=615) ¹		Year 2 (N=585) ¹		Year 3 (N=579) ¹	
	Value	95% CI	Value	95% CI	Value	95% CI
% ECD change from baseline						
n	615	–	585	–	575	–
Mean (SD)	-3.0 (4.0)	-3.355 – -2.716	-4.0 (4.1)	-4.374 – -3.709	-6.7 (5.4)	-7.130 – -6.254
Median	-2.5	–	-3.5	–	-6.0	–
Minimum, maximum	-43.3, 8.6	–	-36.4, 10.5	–	-41.7, 11.5	–
Distribution of % ECD change from baseline, n (%)²						
Gain >5%	2 (0.3)	0.04 – 1.17	3 (0.5)	0.11 – 1.49	1 (0.2)	0.00 – 0.96
Gain ≥2% to ≤5%	13 (2.1)	1.13 – 3.59	13 (2.2)	1.19 – 3.77	5 (0.9)	0.28 – 2.00
Gain <2% to loss <2%	236 (38.4)	34.51 – 42.35	152 (26.0)	22.47 – 29.74	82 (14.2)	11.42 – 17.27
Loss ≥2% to ≤5%	242 (39.3)	35.47 – 43.34	232 (39.7)	35.67 – 43.75	151 (26.1)	22.55 – 29.86
Loss >5% to ≤10%	108 (17.6)	14.63 – 20.80	149 (25.5)	21.99 – 29.20	217 (37.5)	33.52 – 41.56
Loss >10% to ≤20%	11 (1.8)	0.90 – 3.18	33 (5.6)	3.91 – 7.83	107 (18.5)	15.40 – 21.89
Loss >20% to ≤30%	0	0.00 – 0.60	1 (0.2)	0.00 – 0.95	9 (1.6)	0.71 – 2.93
Loss >30%	3 (0.5)	0.10 – 1.42	2 (0.3)	0.04 – 1.23	3 (0.5) ³	0.11 – 1.51
Missing	0	–	0	–	4	–
ECD <1500 cells/mm²	0	0.00 – 0.60	0	0.00 – 0.63	0	0.00 – 0.64
ECD <1000 cells/mm²	0	0.00 – 0.60	0	0.00 – 0.63	0	0.00 – 0.64

¹ N is the number of eyes present at both the Preoperative and Year 1, Year 2, and Year 3 Visits.

² Percentage calculated as (n/N)*100.

³ One eye reported at >30% loss at the 6 Month, Year 1 and Year 2 Visits. At the Year 3 Visit, ECD loss was 28.7%.

Adverse Events through Year 3

A total of 265 ocular AEs were reported for 29.7% (187/629) eyes of 134 subjects through Year 3 (**Table 26**). No unanticipated AEs or deaths caused by or associated with the use of the device were reported in the study. Refer to Adverse Events section of EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS for more information on events reported in the first 6 months of follow-up.

The most frequently reported AE was increase in IOP due to retained OVD (19.9%; 125/629) in the 1–6 hr postoperative period, steroid response (2.4%, 15/629) from 6 to 31 days postoperative due to steroid response or as a result of SSIs (0.5%, 4/629). The next most frequently reported AEs were corneal endothelial cell loss (2.1%, 13/629; 3 events in the first 6 months), anterior chamber cell/flame (1.9%, 12/629; 11 events in the first 6 months), dry eye (1.3%, 8/629; 4 events in the first 6 months), and retinal surgery (0.6%, 4/629; 3 events in the first 6 months).

Table 26: Cumulative Ocular Adverse Events

Cumulative Ocular AEs	Subjects ¹ (N=327)	Eyes ¹ (N =629)	Events ² N
Eye experienced any ocular AE	134 (41.0)	187 (29.7)	265
Intraocular pressure increased ³	103 (31.5)	137 (21.8)	144
Corneal endothelial cell loss ⁴	13 (4.0)	13 (2.1)	13
Anterior chamber cell/flare ⁵	10 (3.1)	12 (1.9)	12
Retinal surgery	2 (0.6)	4 (0.6)	8
Dry eye ⁶	4 (1.2)	8 (1.3)	8
Vitreous detachment	6 (1.8)	7 (1.1)	7
Corneal epithelial defect	5 (1.5)	6 (1.0)	6
Chalazion	3 (0.9)	4 (0.6)	5
Iris transillumination defect	4 (1.2)	5 (0.8)	5
Intraocular lens repositioning	4 (1.2)	4 (0.6)	4
Visual acuity reduced ⁷	3 (0.9)	3 (0.5)	4
Intraocular lens exchange	3 (0.9)	3 (0.5)	3
Glare/Halo ⁸	2 (0.6)	3 (0.5)	3
Lenticular opacities ⁹	3 (0.9)	3 (0.5)	3
Narrow anterior chamber angle ¹⁰	2 (0.6)	2 (0.3)	3
Retinal detachment ¹¹	2 (0.6)	3 (0.5)	3
Retinal tear	2 (0.6)	2 (0.3)	3
Vitreous floaters	2 (0.6)	3 (0.5)	3
Contact dermatitis	1 (0.3)	2 (0.3)	2
Cataract surgery	2 (0.6)	2 (0.3)	2
Lens rotation	1 (0.3)	2 (0.3)	2
Astigmatism ¹²	1 (0.3)	1 (0.2)	2
Blepharitis	1 (0.3)	2 (0.3)	2
Glaucoma	1 (0.3)	2 (0.3)	2
Punctate keratitis	2 (0.6)	2 (0.3)	2
Vitreous degeneration	2 (0.6)	2 (0.3)	2
Hordeolum	1 (0.3)	1 (0.2)	1
Traumatic anterior subcapsular opacity	1 (0.3)	1 (0.2)	1
Periorbital cellulitis	1 (0.3)	1 (0.2)	1
Intraocular lens removal	1 (0.3)	1 (0.2)	1
Nuclear cataract	1 (0.3)	1 (0.2)	1
Anterior subcapsular cataract	1 (0.3)	1 (0.2)	1
Keratomileusis	1 (0.3)	1 (0.2)	1
Corneal epithelium disorder	1 (0.3)	1 (0.2)	1
Epiretinal membrane	1 (0.3)	1 (0.2)	1
Eye discharge	1 (0.3)	1 (0.2)	1
Eye pain	1 (0.3)	1 (0.2)	1
Iris incarceration	1 (0.3)	1 (0.2)	1

¹ Percentage calculated as (n/N)*100.² Only the first incidence of an event is counted for any given subject or eye.³ IOP 10 mmHg above baseline to a minimum of 25 mmHg or that required a change in the standard postoperative medication regimen or other special treatment was reported as an AE.⁴ Cases of endothelial cell loss that were counted as AEs included only cases of loss >20%.⁵ Anterior chamber cell/flare was reported as an AE if it met criteria for chronic anterior uveitis or was greater than Grade 2 at Visit 2 (Day 5 – 9) or later.⁶ All cases classified by the investigator as not related to the study device.⁷ Loss of CDVA ≥10 letters at any time point >1 week postoperatively was reported as an AE.⁸ Only glare/halo leading to lens explantation was reported as an AE in the PMA study. In the long-term follow-up of the PMA cohort, one subject reported glare/halo OU that was classified as unrelated to the study device by the Investigator.⁹ Refer to **Table 29** for additional information.¹⁰ Only those cases in which the investigator observed a reduction in anterior chamber angle and believed that a SSI was necessary.¹¹ Refer to **Table 89** for rates of retinal detachment in original FDA study of the MICL.¹² Residual astigmatism requiring second surgery of lens rotational repositioning or LASIK enhancement.

A total of 34 ocular SAEs were reported in 1.9% eyes (12/629) through Year 3 (**Table 27**). Retinal surgery was the most common ocular SAE (0.6%, 4/629) with 8 events, followed by 4 cases of lens repositioning (0.6%, 4/629) and 3 cases each of narrow anterior chamber angle, intraocular lens exchange, retinal detachment, and retinal tear. Of the 34 ocular SAEs reported, 14 were classified as related to study device or study device and study procedures. These events include 2 eyes of 2 subjects treated for narrowing of the anterior chamber angle with lens repositioning and subsequent lens exchange and one eye treated for halo and glare with lens explantation, all of which occurred in the first 6 months after ICL surgery; one eye of one subject who developed ASC cataract and underwent cataract surgery with IOL implantation; one eye of one subject treated with lens repositioning for residual astigmatism and one eye of another subject that underwent lens repositioning and subsequent exchange following rotation of the TICL.

Table 27: Ocular SAEs Through Year 3

Cumulative Ocular SAEs ¹	All Eyes (N=629)	
	Eyes ² n (n/N%)	Events n
Eye experienced any ocular SAE	12 (1.9)	34
Retinal surgery	4 (0.6)	8
Intraocular lens repositioning	4 (0.6)	4
Narrow anterior chamber angle	2 (0.3)	3
Intraocular lens exchange	3 (0.5)	3
Retinal detachment	3 (0.5)	3
Retinal tear	2 (0.3)	3
Cataract surgery	2 (0.3)	2
Lens rotation	2 (0.3)	2
Cataract nuclear	1 (0.2)	1
Glare/Halo	1 (0.2)	1
Intraocular lens removal	1 (0.2)	1
Astigmatism	1 (0.2)	1
Anterior subcapsular cataract	1 (0.2)	1
Keratomileusis	1 (0.2)	1

¹ Refer to EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS for discussion of SAEs reported through Month 6. Percentage calculated as (n/N)*100.

² Only the first incidence of an event is counted for any given eye.

Additional Safety Outcomes

Loss of CDVA From Baseline

No clinically significant persistent loss of CDVA ≥ 2 lines (10 letters) was reported through Year 3 (**Table 28**). A transient loss of 2 lines (10 letters) was reported at the Week 1 visit for one (0.2%, 1/629) eye, which resolved at the next visit. Between Month 6 and Year 3, a 2-line loss of CDVA from baseline was reported in 3 (0.5%, 3/629) eyes of 3 subjects. In each of these cases, the change of CDVA was considered not clinically significant by the investigator, as CDVA remained equal to or better than 20/20. Overall, 90.5% of eyes (524/579) reported unchanged or increased CDVA at Year 3 compared with the preoperative visit.

Table 28: Within-Eye Change in CDVA From Preoperative Visit Through Year 3 (All Eyes)

Within Eye Change in CDVA	Week 1 (N=624)	Month 1 (N=626)	Month 3 (N=624)	Month 6 (N=619)	Year 1 (N=615)	Year 2 (N=585)	Year 3 (N=579)
Change in Letters Read from Preop Visit (Letters Read)							
n	624	626	624	619	615	581	579
Mean (SD)	3.6 (3.91)	4.3 (3.86)	4.5 (4.07)	4.7 (4.13)	4.6 (4.27)	4.2 (4.02)	4.4 (4.09)
Median	4.0	4.0	4.0	5.0	5.0	4.0	4.0
Min, Max	-10, 16	-8, 16	-7, 16	-8, 17	-10, 15	-9, 15	-12, 18
Change in Lines from Preoperative, n (n/N%)							
≥ 2 Lines (≥ 10 Letters) Increase	49 (7.9)	54 (8.6)	82 (13.1)	82 (13.2)	87 (14.1)	60 (10.3)	67 (11.6)
1 to < 2 Lines (5 to <10 Letters) Increase	202 (32.4)	236 (37.7)	225 (36.1)	242 (39.1)	224 (36.4)	218 (37.3)	216 (37.3)
< 1 Line (<5 Letters) Increase	233 (37.3)	237 (37.9)	215 (34.5)	196 (31.7)	189 (30.7)	195 (33.3)	212 (36.6)
No Change	63 (10.1)	46 (7.3)	53 (8.5)	41 (6.6)	44 (7.2)	45 (7.7)	29 (5.0)
< 1 Line (<5 Letters) Decrease	66 (10.6)	42 (6.7)	39 (6.3)	49 (7.9)	62 (10.1)	55 (9.4)	44 (7.6)
1 to < 2 Lines (5 to <10 Letters) Decrease	10 (1.6)	11 (1.8)	10 (1.6)	9 (1.5)	8 (1.3)	8 (1.4)	9 (1.6)
≥ 2 Lines (≥ 10 Letters) Decrease	1 (0.2)	0 (0.0)	0 (0.0)	0 (0.0)	1 (0.2)	0 (0.0)	2 (0.3)
Missing	0	0	0	0	0	4	0

Crystalline Lens Status

Crystalline lens status (Lens Opacity Classification System [LOCS] III grade) reported during the study is provided in **Table 29**. One case of nuclear cataract, resulting in cataract surgery and classified as not related to the EVO ICL, was reported at Month 6. One case of clinically significant ASC cataract, classified as related to the EVO ICL by the investigator, was reported at the Year 2 Visit and subsequently resulted in cataract surgery due to cataract progression and subject complaints of night driving difficulty and glare during the day.

Four cases of asymptomatic lenticular opacity were reported in the study. A case of traumatic ASC opacity caused by blunt trauma in one eye, and a case of asymptomatic cortical opacity reported in another eye were each classified as unrelated to the EVO ICL by the investigators. The remaining 2 cases of asymptomatic ASC opacity reported for one eye in each of two subjects at Year 2, were classified as related to the EVO ICL by the investigator.

Table 29: Crystalline Lens Status by Visit (All Eyes)

Crystalline Lens Status	Preop Visit (N=629) n (n/N %)	Month 1 (N=626) n (n/N %)	Month 3 (N=624) n (n/N %)	Month 6 (N=619) n (n/N %)	Year 1 (N=615) n (n/N %)	Year 2 (N=585) n (n/N %)	Year 3 (N=579) n (n/N %)
Opacity Present?							
Absent	629 (100.0)	625 (99.8)	624 (100.0)	618 (99.8)	612 (99.5)	581 (99.3)	573 (99.0)
Present	0 (0.0)	1 (0.2)	0 (0.0)	1 (0.2)	3 (0.5)	4 (0.7)	6 (1.0)
Missing	0	0	0	0	0	0	0
If Opacity present¹:							
Nuclear Color							
N	0	0	0	1	1	3	4
Mean (SD)	-	-	-	0.50 (.)	0.50 (.)	1.77 (0.577)	1.80 (0.476)
Median	-	-	-	0.50	0.50	2.10	2.00
Min, Max	-	-	-	0.5, 0.5	0.5, 0.5	1.1, 2.1	1.1, 2.1
Nuclear Opalescence							
N	0	0	0	1	1	3	4
Mean (SD)	-	-	-	1.80 (.)	2.60 (.)	1.80 (0.000)	1.78 (0.050)
Median	-	-	-	1.80	2.60	1.80	1.80
Min, Max	-	-	-	1.8, 1.8	2.6, 2.6	1.8, 1.8	1.7, 1.8
Cortical Cataract							
N	0	1	0	0	1	1	3
Mean (SD)	-	0.50 (.)	-	-	0.50 (.)	0.50 (.)	1.43 (1.210)
Median	-	0.50	-	-	0.50	0.50	1.00
Min, Max	-	0.5, 0.5	-	-	0.5, 0.5	0.5, 0.5	0.5, 2.8
Posterior Subcapsular Cataract							
N	0	0	0	0	0	0	0
Mean (SD)	-	-	-	-	-	-	-
Median	-	-	-	-	-	-	-
Min, Max	-	-	-	-	-	-	-
Anterior Subcapsular Cataract							
N	0	0	0	0	1	3	3
Mean (SD)	-	-	-	-	2.30 (.)	2.93 (0.551)	3.53 (1.124)
Median	-	-	-	-	2.30	3.20	3.80
Min, Max	-	-	-	-	2.3, 2.3	2.3, 3.3	2.3, 4.5

¹ Lens Opacity Classification System [LOCS] III grade used in study

Note: The Lens Opacities Classification System III (LOCS III)* provides standardized grading of nuclear color (NC), nuclear opalescence (NO), cortical cataract (C), and posterior subcapsular (P) cataract. In this study, the grading scale for P was also used to grade anterior subcapsular cataract. Cataract severity was graded on a decimal scale, and the standards have regularly spaced intervals on a decimal scale. The scale ranges from 0.1 (clear or colorless) to 5.9 (very opaque [in cases of C and P] or 6.9 (very opaque or brunescent [in cases of NO and NC])).

* Reference: Chylack LT Jr, Wolfe JK, Singer DM, Leske MC, Bullimore MA, Bailey IL, Friend J, McCarthy D, Wu SY. The Lens Opacities Classification System III. The Longitudinal Study of Cataract Study Group. Arch Ophthalmol. 1993 Jun;111(6):831-6.

Intraocular Pressure

Mean (SD) IOP values were similar and stable between Month 3 and Year 3 Visits (**Table 30**). Refer to EVO/EVO+ ICL/TICL PMA CLINICAL TRIAL AND RESULTS for discussion of postoperative increases in IOP due to incomplete removal of OVD in the 1 – 6 hour early postoperative period or steroid response in the first month after surgery.

Table 30: Intraocular Pressure by Visit (All Eyes)

IOP (mmHg)	Preop (N=629)	Day 1 (N=628)	Week 1 (N=624)	Month 1 (N=626)	Month 3 (N=624)	Month 6 (N=619)	Year 1 (N=613)	Year 2 (N=585)	Year 3 (N=579)
n	629	628	624	626	624	619	613	585	579
Mean (SD)	15.9 (2.83)	15.1 (3.30)	16.0 (3.12)	16.2 (3.08)	15.3 (2.56)	15.3 (2.42)	15.2 (2.55)	15.2 (2.59)	15.5 (2.64)
Median	16.0	15.0	16.0	16.0	16.0	15.0	15.0	15.0	16.0
Min, Max	9, 22	4, 24	7, 34	8, 33	8, 22	8, 21	7, 32	9, 22	8, 22

N=number of eyes that completed the visit

Gonioscopy

Overall, data show the degree of angle narrowing that is expected to occur following implantation of EVO ICL lenses. No significant increase in pigmentation of the trabecular meshwork and no new peripheral anterior synechiae were observed.

Table 31: Gonioscopy Through Year 3 in All Eyes

Gonioscopy	Preoperative Visit (N=629) Eyes n (n/N%)	Month 6 (N=619) Eyes n (n/N %)	Year 1 (N=615) Eyes n (n/N %)	Year 2 (N=585) Eyes n (n/N %)	Year 3 (N=579) Eyes n (n/N %)
Angle Grade					
0	0	0	0	0	0
1	0	2 (0.3) ²	0	0	0
2	0	10 (1.6)	3 (0.5)	2 (0.3)	1 (0.2)
3	66 (10.5)	93 (15.0)	93 (15.1)	83 (14.2)	74 (12.8)
4	563 (89.5)	512 (82.7)	517 (84.1)	494 (84.4)	504 (87.0)
Missing	0	2	2	6	0
Pigmentation Grade					
0	497 (79.0)	471 (76.1)	445 (72.4)	412 (70.4)	413 (71.3)
1	102 (16.2)	117 (18.9)	131 (21.3)	134 (22.9)	126 (21.8)
2	13 (2.1)	13 (2.1)	19 (3.1)	17 (2.9)	22 (3.8)
3	17 (2.7)	16 (2.6)	16 (2.6)	16 (2.7)	18 (3.1)
4	0	0	2 (0.3)	0	0
Missing	0	2	2	6	0
Peripheral Anterior Synechiae					
Absent	628 (99.8)	616 (99.5)	610 (99.2)	578 (98.8)	578 (99.8)
Present, specify clock hours	1 (0.2)	1 (0.2)	3 (0.5)	1 (0.2)	1 (0.2)
0.5 – 2.0	1 (0.2)	1 (0.2)	3 (0.5)	1 (0.2)	1 (0.2)
2.5 – 4.0	0	0	0	0	0
4.5 – 6.0	0	0	0	0	0
6.5 – 8.0	0	0	0	0	0
8.5 – 10.0	0	0	0	0	0
10.5 – 12.0	0	0	0	0	0
Missing	0	2	2	6	0

N=number of eyes that completed the visit

Achieved Vault and Preoperative Biometrics

Achieved vault, measured using OCT, from Month 1 through Year 3 is shown in **Table 32**, and **Table 33** provides the number and percent of eyes with vault measurements <250 µm and >900 µm, as well as mean vault and quartiles for vault at the Month 6, Year 1, Year 2, and Year 3 Visits.

Table 32: Achieved Vault (All Eyes)

Lens Vault, µm	Month 1 (N=626)	Month 3 (N=624)	Month 6 (N=619)	Year 1 (N=615)	Year 2 (N=585)	Year 3 (N=579)
n	620	622	615	615	581	577
Mean (SD)	517.7 (231.7)	509.7 (231.9)	491.9 (227.0)	458.5 (214.6)	422.4 (199.0)	412.1 (196.7)
Median	494.5	480.0	466.0	425.0	390.0	390.0
Min, Max	22, 1399	3, 1390	10, 1240	10, 1158	16, 1100	19, 1059

N=number of eyes that completed the visit

Table 33: Lens Vault Through Year 3 (All Eyes)

Parameter	Month 6 (N=619)	Year 1 (N=615)	Year 2 (N=585)	Year 3 (N=579)
Eyes with vault measurement, n	615	615	581	577
Eyes measured with vault <250 µm, n (n/N%)	79 (13)	92 (15)	106 (18)	120 (21)
Eyes measured with vault >900 µm, n (n/N %)	33 (5)	24 (4)	12 (2)	11 (2)
Mean vault, µm	492	459	422	412
0 th percentile for measured vault, µm	10.0	10.0	16.0	19.0
25 th percentile for measured vault, µm	330.0	306.0	280.0	284.0
50 th percentile for measured vault, µm	466.0	425.0	390.0	390.0
75 th percentile for measured vault, µm	630.0	589.0	538.0	530.0
100 th percentile for measured vault, µm	1240.0	1158.0	1100.0	1059.0

Results of correlation analyses between achieved vault and preoperative biometry through Year 3 are provided in **Table 34**. No statistically significant correlations were observed between achieved vault and keratometry at any time point. Significant correlations were observed between achieved vault and axial length, true ACD, and white-to-white at Month 6, Year 1, Year 2, and Year 3. For UBM, significant correlations were observed between achieved vault and both horizontal and vertical sulcus-to-sulcus at Month 6, Year 1, Year 2, and Year 3. For anterior segment OCT, only 1 significant correlation was observed, between achieved vault and iris pigment end-to-pigment end at Year 3. Crystalline lens rise was defined as the distance between the anterior pole of the crystalline lens and a plane joining the following anatomical landmarks: the iridocorneal recess (angle to angle plane), the iris pigment end (pigment to pigment plane), or the ciliary sulcus (sulcus to sulcus plane). Significant correlations were observed between achieved vault and angle-to-angle and sulcus-to-sulcus at Month 6, Year 1, Year 2, and Year 3.

Table 34 also provides the coefficients of determination (R^2) and p-values for correlations between preoperative parameters and postoperative vault through Year 3. Clinically significant correlations were found for lens diameter, crystalline lens rise as measured from the sulcus-to-sulcus (STS) plane and crystalline lens rise as measured from the angle-to-angle (ATA) plane. For lens diameter, vault tended to increase with increasing diameter (**Figure 5** – see Pre-Approval EVO/EVO+ ICL/TICL PMA Clinical Trial and Results: Vault). Crystalline lens rise was defined as the distance between the anterior pole of the crystalline lens and a plane joining the following anatomical landmarks: the iridocorneal recess (angle to angle plane), the iris pigment end (pigment to pigment plane), or the ciliary sulcus (sulcus to sulcus plane). Vault tended to decrease with increasing crystalline lens rise measured from the STS and ATA planes.

Table 34: Significance of Associations Between Achieved Vault and Preoperative Biometrics

Significance of Associations Between Achieved Vault and Preoperative Biometrics	Month 6 (N=619)		Year 1 (N=615)	Year 2 (N=585)	Year 3 (N=579)
	R^2	p-value	p-value	p-value	p-value
Keratometry					
Steep diopter	<0.01	0.2971	0.4812	0.4585	0.3029
Steep axis	NR	0.8041	0.6586	0.7475	0.9741
Flat diopter	<0.01	0.4520	0.7009	0.5612	0.2047
Flat axis	NR	0.4679	0.4375	0.7619	0.8929
Biometry					
Axial length	<0.01	0.0120	0.0015	0.0006	<0.0001
True ACD ¹	0.07	<0.0001	<0.0001	<0.0001	<0.0001
White-to-white	0.07	<0.0001	<0.0001	<0.0001	<0.0001
Pupil size		0.3328	0.3146	0.2380	0.1615
UBM, Substudy Only²	n=228		n=226	n=207	n=199
Horizontal sulcus-to-sulcus	0.05	0.0010	0.0006	0.0019	<0.0001
Vertical sulcus-to-sulcus	0.05	0.0007	0.0002	0.0005	0.0001
Anterior Segment OCT, Substudy Only²	n=203		n=201	n=182	n=174
True ACD ²	NR	0.3134	0.4002	0.1083	0.1068
Central corneal thickness	0.01	0.2458	0.2453	0.5228	0.9319
Angle-to-angle distance	<0.01	0.5916	0.5822	0.5121	0.1081
Iris pigment end-to-pigment end, distance	0.02	0.0663	0.0510	0.6247	0.0473
Crystalline Lens Rise³					
Angle-to-angle	n=37		n=37	n=37	n=37
	0.40	0.0001	<0.0001	<0.0001	0.0001
Pigment-to-pigment	n=22		n=22	n=22	n=22
	<0.01	0.9061	0.8067	0.8764	0.4327
Sulcus-to-sulcus	n=144		n=142	n=123	n=115
	0.22	<0.0001	<0.0001	<0.0001	<0.0001

¹ True ACD is defined as the distance from the apex of the posterior corneal surface to the apex of the anterior crystalline lens surface.

² Only sites with appropriate equipment and experience as required for the study were included in this sub study.

³ Methodology for calculating crystalline lens rise is site specific. Measurements were taken from the pigment-to-pigment plane, the sulcus-to-sulcus plane, or from the angle-to-angle plane.
NR = Not Reported

Effectiveness Findings

Accuracy of Refractive Outcome

Mean (SD) preoperative MRSE was -7.62 (2.76) D in all eyes. Postoperative mean (SD) MRSE at the Month 6, and Year 1, Year 2 and Year 3 Visits is provided in **Table 35**.

As shown in **Table 36**, 99% of eyes were within ± 1.00 D of target MRSE at each time point from Month 6 to Year 3; at least 90% of eyes were within ± 0.50 D of target MRSE during that time frame.

Table 35: MRSE by Visit

MRSE (D)	Preop	Month 6	Year 1	Year 2	Year 3
All Eyes (N)	629	619	615	585	579
Mean (SD)	-7.62 (2.76)	-0.08 (0.33)	-0.10 (0.34)	-0.11 (0.29)	-0.12 (0.30)
Median	-7.38	0.000	0.00	-0.12	-0.12
Min, Max	-15.62, -3.00	-3.88, 1.12	-4.25, 1.00	-1.75, 0.75	-2.00, 1.00

Table 36: MRSE Within $\pm 0.50D$ and $\pm 1.00D$ of Target Through Year 3 (All Eyes)

MRSE	n	Month 6 (N=619) Proportion (95% CI)	n	Year 1 (N=615) Proportion (95% CI)	n	Year 2 (N=585) Proportion (95% CI)	n	Year 3 (N=579) Proportion (95% CI)
$\pm 0.50 D$	560	0.905 (0.8788 – 0.9266)	559	0.909 (0.8834 – 0.9305)	529	0.904 (0.8775 – 0.9269)	525	0.907 (0.8801 – 0.9292)
$\pm 1.00 D$	613	0.990 (0.9790 – 0.9964)	609	0.990 (0.9789 – 0.9964)	580	0.991 (0.9802 – 0.9972)	573	0.990 (0.9776 – 0.9962)

Visual Acuity

Table 37: UDVA 20/40 or Better Through Year 3 (Eyes with Preoperative CDVA 20/20 or better)

	Month 6		Year 1		Year 2		Year 3	
	n/N	(95% CI)	n/N	(95% CI)	n/N	(95% CI)	n/N	(95% CI)
20/20 or better	407/506 (80.4)	76.71 – 83.80	411/503 (81.7)	78.05 – 84.99	381/480 (79.4)	75.48 – 82.91	375/470 (79.8)	75.87 – 83.33
20/40 or better	503/506 (99.4)	98.28 – 99.88	502/503 (99.8)	98.90 – 99.99	473/480 (98.5)	97.02 – 99.41	467/470 (99.4)	98.15 – 99.87
Missing UDVA	0	-	0	-	4	-	0	-

Table 38: CDVA 20/40 or Better Through Year 3 (Eyes with Preoperative CDVA 20/20 or better)

	Month 6		Year 1		Year 2		Year 3	
	n/N (%)	95% CI	n/N	(95% CI)	n/N	(95% CI)	n/N	(95% CI)
20/20 or better	500/506 (98.8)	97.44 – 99.56	496/503 (98.6)	97.15 – 99.44	471/480 (98.1)	96.47 – 99.14	462/470 (98.3)	96.67 – 99.26
20/40 or better	506/506 (100.0)	99.27 – 100.00	503/503 (100.0)	99.27 – 100.00	476/480 (99.2)	97.88 – 99.77	470/470 (100.0)	99.22 – 100.00
Missing CDVA	0	-	0	-	4	-	0	-

Study Strengths and Weaknesses

Strengths

1. This study provides long term follow-up of subjects implanted with EVO/EVO+ ICL/TICL lenses, including assessment of endothelial cell density.
2. Sufficient retention of subjects over the full course of study provides valid scientific evidence of long term safety and effectiveness.
3. Adherence to protocol requirements including low numbers of protocol deviations insures the validity of the data.

Weaknesses

1. The absence of randomized control group, e.g., subjects implanted with the parent lenses, represents a limitation of the study.
2. Variation in enrollment among sites may limit wider applicability of the findings.

POST-APPROVAL IOP STUDY OF THE EVO/EVO+ ICL/TICL

The STAAR EVO/EVO+ ICL/TICL lens models are modifications of the previously approved Visian MICL and TICL phakic IOL models. A central 360 μm hole was incorporated into the lens to eliminate the requirement for laser peripheral iridotomy to prevent pupillary block. The PMA study of the EVO/EVO+ ICL/TICL included experienced ICL surgeons and was successful in demonstrating that the central hole modification eliminated the need for preoperative iridotomies (which were needed for the parent Visian MICL/TICL lenses). However, at a 1 – 6 hour scheduled postoperative visit, 125 of 629 implanted eyes (19.9%) experienced elevated IOP (no prophylactic IOP lowering medication was permitted in this study). Of these 125 eyes, 85 had IOP ≥ 30 mmHg, 38 had IOP ≥ 40 mmHg, 24 had IOP ≥ 50 mmHg, and 11 had IOP ≥ 60 mmHg. Survey of study investigators revealed that an increased rate of IOP spikes was associated with less thorough removal of HPMC (hydroxypropylmethylcellulose 2%) OVD after lens implantation. To mitigate this risk, STAAR modified their required physician certification training program to specifically include instructions regarding the most effective methods of OVD removal. As a condition of approval, the FDA required that STAAR conduct a short post-approval clinical study to assess the effectiveness of this training program in lowering the rate of early IOP spikes. This study and its results are described in this section.

Study Objective

The objective of this study was to assess the rate of early intraocular pressure (IOP) increases following EVO/EVO+ ICL/TICL surgery by surgeons who were trained and certified in the EVO/EVO+ ICL/TICL Physician Certification Program compared with rates observed in the original PMA study.

Study Design

The EVO/EVO+ ICL/TICL lenses were evaluated in a prospective, multicenter, single-arm, 2-week follow-up post-approval study. Subjects who met study eligibility criteria received bilateral implantation and were treated and seen for 4 scheduled study visits (per eye) at 1–6 hours (Visit 1), 1 day (Visit 2), 1 week (Visit 3, Day 5–9) and 2 weeks (Visit 4, Day 10–18) after surgery.

Study Population

A total of 408 eyes of 205 subjects were enrolled at 8 sites in the U.S.

Table 39: Demographics

Demographics	Subjects (N=205)
Gender, n (%)	
Male	79 (38.5%)
Female	126 (61.5%)
Race, n (%)	
Caucasian	176 (85.9%)
Black/African American	1 (0.5%)
Asian	24 (11.7%)
Native Hawaiian/Pacific Islander	0
American Indian/Alaska Native	2 (1.0%)
Other	2 (1.0%)
Ethnicity, n (%)	
Hispanic or Latino	33 (16.1%)
Not Hispanic or Latino	172 (83.9%)
Age, years	
n	205
Mean (SD)	36.1 (4.38)
Median	36.0
Minimum, maximum	26, 45

* Percentage calculated as (n/N)*100.

Study Endpoints

Primary endpoints, evaluated in primary eyes only, were:

- the proportion of primary eyes with IOP \geq 30 mmHg at 1 – 6 hours postoperative,
- the proportion of primary eyes with IOP \geq 40 mmHg at 1 – 6 hours postoperative.

Secondary endpoints, evaluated in fellow eyes only, were:

- the proportion of fellow eyes with IOP \geq 30 mmHg at 1 – 6 hours postoperative,
- the proportion of fellow eyes with IOP \geq 40 mmHg at 1 – 6 hours postoperative.

Other endpoints were evaluated in primary and all (primary + fellow) eyes for the entire postoperative follow-up period (including the 1–6 hour postop Visit through Visit 4, postop Day 10 –18):

- Rates of increased IOP due to retained ophthalmic viscoelastic device (OVD),
- Rates of increased IOP due to other causes (e.g., pupillary block, steroid response, etc.),
- Rates of all categories of AEs.

Accountability

Table 40: Accountability – Primary Eyes

Eye Status	Total #	Op Visit (Day 0) n (%)	Postop V1 (1–6 hr) n (%)	Postop V2 (Day 1) n (%)	Postop V3 (Day 5–9) n (%)	Postop V4 (Day 10–18) n (%)
Eyes treated (N)	205					
Available for Analysis		205 (100%)	205 (100%)	205 (100%)	199 (97.1%)	202 (98.5%)
Missing Eye/Data						
Discontinued		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Missing at scheduled visit but seen earlier/later ¹		0 (0.0)	0 (0.0)	0 (0.0)	5 (2.4%)	3 (1.5%)
Not seen but accounted for ²		0 (0.0)	0 (0.0)	0 (0.0)	1 (0.5%)	0 (0.0)
Lost to follow-up		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Accountability ⁴		205/205 (100%)	205/205 (100%)	205/205 (100%)	199/205 (97.1%)	202/205 (98.5%)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Not seen but accounted for=The total number of eyes that missed the visit but have not been discontinued/lost to follow-up

³ Active: represents the total number of eyes that have not reached the time associated with the visit. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ %Accountability = [Available for Analysis / (N - [Discontinued + Active])

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Table 41: Accountability – All Eyes

Eye Status	Total #	Op Visit (Day 0) n (%)	Postop V1 (1-6 hr) n (%)	Postop V2 (Day 1) n (%)	Postop V3 (Day 5-9) n (%)	Postop V4 (Day 10-18) n (%)
Eyes treated (N)	408					
Available for Analysis		408 (100%)	408 (100%)	407 (99.8%)	401 (98.3%)	402 (98.5%)
Missing Eye/Data						
Discontinued		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Missing at scheduled visit but seen earlier/later ¹		0 (0.0)	0 (0.0)	1 (0.2%)	5 (1.2%)	6 (1.5%)
Missing but accounted for		0 (0.0)	0 (0.0)	0 (0.0)	2 (0.5%)	0 (0.0)
Lost to follow-up ²		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
Accountability ⁴		408/408 (100%)	408/408 (100%)	407/408 (99.8%)	401/408 (98.3%)	402/408 (98.5%)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Not seen but accounted for=The total number of eyes that missed the visit but have not been discontinued/lost to follow-up

³ Active: represents the total number of eyes that have not reached the time associated with the visit. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ %Accountability = [Available for Analysis / (N - [Discontinued + Active])]

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Safety Outcomes

Incidence of IOP ≥ 30 mmHg and ≥40 mmHg

The analysis of the primary study endpoint was based on the superiority of the proportion of primary eyes with IOP ≥ 30 mmHg and IOP ≥ 40 mmHg at 1–6 hours postoperatively compared with the outcomes of the original PMA study of the EVO/EVO+ ICL/TICL, i.e., 15.3% (50/327) primary eyes with IOP ≥ 30 mmHg, and 7.0% (23/327) primary eyes with IOP ≥ 40 mmHg. In this PAS, 7.3% (15/205) and 2.0% (4/205) of primary eyes presented with IOP ≥ 30 mmHg and IOP ≥ 40 mmHg, respectively, at the 1–6 hour postoperative visit (**Table 42**).

Table 42: Incidence of Primary Eyes with IOP ≥ 30 mmHg and ≥40 mmHg at 1-6 hr Postoperative

Endpoint:	IOP ≥ 30 mmHg		
Eyes with IOP spike/ Implanted Eyes	n/N	Percent (%)	95% CI(%) ¹
EVO PAS (new enrollment study)	15/205	7.3	4.2, 11.8
EVO PMA Study (for approval)	50/327	15.3	-
Difference (%)	-	8.0	-
p-value of difference ²	p=0.0004		
Endpoint:	IOP ≥ 40 mmHg		
Eyes with IOP spike/ Implanted Eyes	n/N	Percent	95% CI(%) ¹
EVO PAS (new enrollment study)	4/205	2.0	0.5, 4.9
EVO PMA Study (for approval)	23/327	7.0	-
Difference (%)	-	5.1	-
p-value of difference ²	p=0.0010		

¹ 95% Clopper-Pearson confidence interval for proportion of eyes meeting the criterion

² Statistical significance of the difference between the rates of increased IOP in the EVO/EVO+ ICL/TICL PAS study and the original EVO/EVO+ ICL/TICL PMA study.

A p-value < 0.05 indicates a statistically significant reduction in the rate of increased IOP in the EVO/EVO+ ICL/TICL PAS study.

The secondary study endpoint was the proportion of secondary (fellow) eyes that had IOP ≥ 30 mmHg and ≥ 40 mmHg at 1–6 hours postoperative. Per protocol, analysis of fellow eye rates was used to support the findings of the two primary eye analyses. As shown in **Table 43**, 7.9% (16/203) and 3.0% (6/203) of fellow eyes presented with IOP ≥ 30 mmHg and IOP ≥ 40 mmHg, respectively. In comparison, 15.3% (50/327) and 11.9% (36/302) of fellow eyes in the original PMA study of the EVO/EVO+ ICL/TICL presented with IOP ≥ 30 mmHg and IOP ≥ 40 mmHg, respectively at the 1–6 hour postoperative visit.

Based on these results, the incidence of IOP ≥ 30 mmHg and IOP ≥ 40 mmHg in primary and fellow eyes was significantly lower than that reported in the original PMA study.

Table 43: Incidence of Fellow Eyes with IOP ≥ 30 mmHg and ≥40 mmHg at 1-6 hr Postoperative

Endpoint: Eyes with IOP spike/ Implanted Eyes	n/N	IOP ≥ 30 mmHg Percent (%)	95% CI(%) ¹
EVO PAS (new enrollment study)	16/203	7.9	4.6, 12.5
EVO PMA Study (for approval)	50/327	15.3	-
Difference (%)	-	8.0	-
p-value of difference ²	p=0.0012		
Endpoint: Eyes with IOP spike/ Implanted Eyes	n/N	IOP ≥ 40 mmHg Percent	95% CI(%) ¹
EVO PAS (new enrollment study)	6/203	3.0	1.1, 6.3
EVO PMA Study (for approval)	36/302	11.9	-
Difference (%)	-	5.1	-
p-value of difference ²	p=0.0010		

¹ 95% Clopper-Pearson confidence interval for proportion of eyes meeting the criterion

² Statistical significance of the difference between the rates of increased IOP in the EVO/EVO+ ICL/TICL PAS study and the original EVO/EVO+ ICL/TICL PMA study.

A p-value < 0.05 indicates a statistically significant reduction in the rate of increased IOP in the EVO/EVO+ ICL/TICL PAS study.

Rates of increased IOP due to retained OVD and rates of increased IOP due to other causes (e.g., pupillary block, steroid response, etc.) over the entire postoperative period were evaluated for all treated eyes. Any increased IOP ≥ 10 mmHg above baseline to a minimum of 25 mmHg or any increase in IOP that was treated in any way, was reported as an AE.

A total of 49 AEs were reported for increased IOP through Visit 4 (**Table 44**). **Table 45** provides the distribution of maximum IOP in these cases. Forty-two (42) of the 49 events occurred at the 1–6 hour postoperative IOP check (Postop Visit 1). In all cases, gonioscopy and/or AS-OCT performed prior to any treatment of IOP demonstrated that the anterior chamber angle was open. No cases of angle closure or pupillary block were reported. For all 42 events, increased IOP either resolved without treatment or was treated with hypotensive medication alone, aqueous tap alone or a combination of hypotensive medication and aqueous tap (**Table 46**). No other procedures to reduce IOP were performed in this study. These events were classified as related to the surgical procedure, specifically to the use of OVD, and not related to the study lens. The remaining 7 of 49 AEs of elevated IOP occurred at postoperative Visit 2 (Day 1), postoperative Visit 3 (days 5–9) or Visit 4 (day 10–18). In all but 1 case, elevated IOP was attributed to a postoperative topical corticosteroid response and either resolved with discontinuation/tapering of steroid (5 events) or initiation of hypotensive medication (1 event). In the 7th case, an AE of elevated IOP was reported at Visit 2 that was attributed to retained OVD; no intervention was taken, and the event resolved without sequelae at Visit 4.

Table 44: Other Endpoint – Cumulative Elevated IOP AEs by Visit for all Implanted (Primary + Fellow) Eyes

N= 408	Postop V1 (1-6 hr)		Postop V2 (Day 1)		Postop V3 (1 Week, Day 5-9)		Postop V4 (2 Week, Day 10-18)	
	n (%) ¹ Eyes	n Events	n (%) ¹ Eyes	n Events	n (%) ¹ Eyes	n Events	n (%) ¹ Eyes	n Events
Increased IOP	42 (10.2%)	42	43 (10.5%)	43	44 (10.8%)	44	45 (11.0%)	49
Attributed to:								
Retained OVD	42 (10.2%)	42	43 (10.5%)	43	43 (10.5%)	43	43 (10.5%)	43
Steroid Response	0	0	0	0	1 (0.2%)	1	6 (1.5%)	6

¹ Percentage calculated as (n/N)*100.

Table 45: Maximum IOP for Adverse Events of Elevated IOP

Adverse Event - Elevated IOP	All Eyes (N=408) n (%) ¹
Number of elevated IOP events	49 (12.0)
Maximum IOP (mmHg)	
< 30	16 (3.9)
≥ 30	33 (8.1)
≥ 40	10 (2.5)
≥ 50	7 (1.7)
≥ 60	2 (0.5)
≥ 70	1 (0.2)

¹ Percentage calculated as (n/N)*100

Table 46: Treatment for Adverse Events of Elevated IOP

N=408	Postop V1 (1-6 hr) n (%) ¹	Postop V2 (Day 1) n (%) ¹	Postop V3 (1 Week, Day 5-9) n (%) ¹	Postop V4 (2 Week, Day 10-18) n (%) ¹
Events Treated with Concomitant Medication(s) and Aqueous Tap	15 (3.7%)	0 (0%)	0 (0%)	0 (0%)
Events Treated with Concomitant Medication(s) Only	13 (3.2%)	0 (0%)	1 (0.2%)	0 (0%)
Events Treated with Aqueous Tap Only	3 (0.7%)	0 (0%)	0 (0%)	0 (0%)
Events Not Treated	11 (2.7%)	1 (0.2%)	0 (0.0%)	5 (1.2%)

¹ Percentage calculated as (n/N)*100

Adverse Events

A total of 67 ocular AEs were reported for 57 eyes (14.0%) of 43 subjects in this study (**Table 47**). The most frequently reported AE was increased IOP with a total of 49 events occurring in 45 eyes (11%) of 33 subjects (refer to previous section for additional information). All other ocular AEs that occurred during this clinical trial were previously anticipated in nature, severity, and frequency based on prior clinical studies as well as the published literature regarding the ICL family of lenses.

Table 47: Cumulative Ocular Adverse Events

Cumulative Ocular AEs	All Eyes (N= 408)	
	Eyes n (%) ¹	Events n
Eyes experienced any ocular AE ²	57 (14.0%)	67
Intraocular pressure increased	45 (11.0%)	49
Anterior chamber cell/flare	9 (2.2%)	9
Superficial punctate keratitis secondary to cosmetic product	2 (0.5%)	2
Corneal abrasion	1 (0.2%)	1
Epithelial defect	1 (0.2%)	1
Eyelid contact dermatitis	1 (0.2%)	1
Incorrect lens power implantation	1 (0.2%)	1
Intraocular lens exchange	1 (0.2%)	1
Intraocular cilium	1 (0.2%)	1
Removal of intraocular cilium	1 (0.2%)	1

¹ Percentage calculated as (n/N)*100.

² Only the first incidence of an event is counted for any given subject.

A total of 4 SAEs (including 2 secondary surgical interventions (SSIs) reported, per protocol, as SAEs) were reported for 2 primary eyes of 2 subjects. In one subject, an SAE of incorrect lens power implantation was reported requiring an SSI of intraocular lens exchange. In the other subject, an intraocular cilium was noted at the 1–6 hour postoperative visit requiring an SSI of removal of the cilium. All 4 SAEs were reported as related to the study protocol and unrelated to the EVO/EVO+ ICL/ TICL.

Effectiveness Outcomes

The mean MRSE improved from -7.87 D preoperative to -0.11 D at the 2 Week postoperative visit (**Table 48**).

Table 48: MRSE by Visit

	Preop	Postop V3 (1 Week, Day 5-9)	Postop V4 (2 Week, Day 10-18)
All Eyes (N)	408	406	408
MRSE (D)			
Mean (SD)	-7.87 (2.47)	-0.13 (0.32)	-0.11 (0.31)
Median	-7.88	-0.13	0.00
Min, Max	-14.88, -3.00	-1.38, 0.75	-1.38, 0.75
Missing	0	2	0

At the 2 Week Visit, 84.6% (345/408) of all eyes had an MRSE within ± 0.50 D and 98.5% (402/408) had an MRSE within ± 1.00 D of target. In addition, 99.7% (396/397) of all eyes with preoperative CDVA of 20/20 or better reported a postoperative UDVA of 20/40 or better (**Table 49**).

Table 49: MRSE Within ± 0.50 D and ± 1.00 D of Target and UDVA 20/40 or Better

	All Eyes (N=408)	
	n	Proportion (95% CI) ¹
± 0.50 D	345	0.846 (0.8068, 0.8793)
± 1.00 D	402	0.985 (0.9683, 0.9946)
UDVA 20/40 or better ²	396	0.997 (0.9860, 0.9999)

¹ 95% Clopper-Pearson confidence interval for the proportion of eyes meeting the criterion.

² Only eyes with CDVA 20/20 or better at the Preoperative Visit are included in summary statistics for UDVA 20/40 or better.

Visual acuities at the Preoperative Visit and at the final 2 Week study visit are provided in the following tables.

Table 50: UDVA

N (408)	Preop	Postop V4 (2 Week, Day 10-18)
	n, % ¹	n, % ¹
20/12.5 or better	0, 0%	8, 2.0%
20/16 or better	0, 0%	246, 60.3%
20/20 or better	0, 0%	371, 90.9%
20/40 or better	0, 0%	407, 99.8%

¹ Percentage calculated as (n/N)*100.

Table 51: CDVA

N (408)	Preop	Postop V4 (2 Week Visit, Day 10-18)
	n, % ¹	n, % ¹
20/12.5 or better	1 (0.2)	16 (3.9)
20/16 or better	155 (38.0)	364 (89.2)
20/20 or better	397(97.3)	408 (100)
20/25 or better	408 (100)	408 (100)

¹ Percentage calculated as (n/N)*100.

Study Strengths and Limitations

Strengths of this study include the following observations:

1. This study demonstrated that training of surgeons in the thorough removal of OVD following EVO/EVO+ ICL/TICL implantation significantly reduced the incidence of elevated IOP at the 1 – 6 hour postoperative visit compared with the outcomes of the original PMA study of the EVO/EVO+ ICL/TICL.
2. Gonioscopy or anterior segment imaging performed during events of elevated IOP demonstrated that the anterior chamber angle was open in all cases, confirming that retained OVD was the sole etiology of increased IOP at the 1 – 6 hour postoperative visit.
3. The resolution of all adverse events of increased IOP confirmed that measurement and treatment of IOP as indicated at 1 – 6 hours postoperative represent important mitigations of potential sequelae.

Limitations of this study include the following observations:

1. The short term follow-up period, while sufficient to address this study's objectives, did not allow collection of additional longer term data on effectiveness and safety.
2. This study incorporated a historical control group, i.e., subjects in the original PMA study, rather than a concurrent randomized control group.

PRE-APPROVAL VISIAN TICL LENS CLINICAL TRIAL AND RESULTS

The Visian TICL lens was evaluated in a prospective nonrandomized study of 210 eyes of 124 subjects, 194 eyes of which were followed for 12 months. Study Cohort demographics are as follows:

Table 52: Demographics – Visian TICL Study

N=124 (Subjects)	
Age	
Mean (SD)	35.0 (6.8) yrs
Range	21 to 45 yrs
Race	
	n, % ¹
Caucasian	102, 82.3%
Hispanic	10, 8.1%
Black	6, 4.8%
Other	6, 4.8%
Gender	
Female	69, 55.6%
Male	55, 44.4%

¹ Percentage calculated as (n/N)*100.

Adverse Events and Complications

A total of 210 eyes of 124 subjects were evaluated in the clinical trial of the Visian TICL lens. Anterior subcapsular opacities, not all clinically significant, were observed postoperatively in six eyes (6/210, 2.9%). Two of these 6 cases (2/210, 1.0%) had a clinically significant cataract. The remaining 4 cases were asymptomatic with 20/16 or better CDVA and 20/25 or better UCVA at their last reported visit. There were no cases of greater than trace nuclear color, nuclear opalescence, cortical or posterior subcapsular changes preoperatively or at any postoperative visit.

A total of 3 eyes (3/210, 1.4%) reported a loss of ≥ 2 lines of CDVA between the preoperative and 12 month visit. A loss of > 2 lines of CDVA (20/25 to 20/50) occurred at the 12 month visit in one eye (1/210, 0.5%) due to anterior subcapsular cataract. There was no information regarding treatment or resolution at the time of study closure. A loss of 2 lines of CDVA was reported in two eyes (2/210, 1.0%). In one eye, the preoperative CDVA was 20/12.5 and at the 12 month visit the CDVA was 20/20. There were no lens opacities noted at any visit and the patient consistently rated her satisfaction with the procedure as very satisfied. The other eye was amblyopic with preoperative CDVA of 20/40 and postoperative CDVA of 20/60 at both the 6 and 12 month visits. This patient was subsequently seen 5 months after the 12 month visit and CDVA was within 1 line of preoperative CDVA. No eyes (0%) had CDVA worse than 20/40 (if preoperative CDVA 20/20 or better) between 1 and 12 months postoperative.

Corneal edema and iritis were not reported after the 1 week visit. There was 1 case (1/210, 0.5%) with a retinal detachment. One eye (1/210, 0.5%) had increased IOP at one day postoperative, which was related to a pupillary block and resolved with an additional Nd:YAG iridotomy. IOP at the one day follow up visit after Nd:YAG iridotomy was 12 mmHg. At the final 12 month post op visit, the BCVA was 20/25 and IOP was 14 mmHg. One eye (1/210, 0.5%) experienced an IOP > 25 mmHg at 6 months postoperative, which dropped to 17 mmHg at 12 months. Two eyes (2/210, 1.0%) of two subjects experienced an increase of > 10 mmHg over preoperative IOP during the 12 month follow-up period. These eyes experienced IOP increases from 8 mmHg to 21 mmHg and from 10 mmHg to 22 mmHg. No treatment was reported in any of these cases. No cases of endophthalmitis, corneal ulcer, ocular hypertension, corneal haze/edema (after 1 week), or corneal melting were reported during the study. The 8 cases (8/210, 3.8%) of surgical intervention all had improvement/no change in CDVA or no significant loss in CDVA (1 line in 1 case) at the last follow-up visit.

Incidence of key AEs/complications are provided in **Table 53**. For a benchmark, they are compared with the ISO historical rate for posterior chamber IOLs for aphakia, implanted in the capsular bag (from ISO 11979-7). Surgical reinterventions occurred in 3.8% (8/210) of eyes. Details concerning the types of surgical reinterventions are provided in **Table 54**.

Table 53: Incidence of Key Adverse Events and/or Complications – Visian TICL Study

Adverse Event	Cumulative N=210 Eyes n/210, %*	ISO ¹ Historical Rate %	Persistent (12 Months) n/194, %*	ISO Historical Rate %
Endophthalmitis	0, 0%	0.1%	0, 0%	---
Hyphema ⁴	0, 0%	---	0, 0%	---
Hypopyon	0, 0%	0.3%	0, 0%	---
IOL Dislocation	0, 0%	0.1%	0, 0%	---
Cystoid Macular Edema	0, 0%	3.0%	0, 0%	0.5%
Raised IOP Requiring Treatment ⁴	1, 0.5%	---	0, 0%	0.4%
Pupillary Block	1, 0.5%	0.1%	0, 0%	---
Retinal Detachment ²	1, 0.5%	0.3%	0, 0%	---
Surgical Reintervention ³	8, 3.8%	0.8%	0, 0%	---
CDVA loss ≥ 2 lines ⁴	3, 1.5%	---	3, 1.5%	---
Corneal Edema ⁴ (after 1 week)	0, 0%	---	0, 0%	0.3%
Iritis ⁴ (after 1 week)	0, 0%	---	0, 0%	0.3%
Anisocoria ⁴	1, 0.5%	---	0, 0%	---

¹ ISO-11979-7: Ophthalmic implants– Intraocular Lenses Part 7: Clinical Investigations

² Comparison should be made to literature for retinal detachment rates for high myopia. Retinal detachment rates increase with increasing myopia. The risk of retinal detachment within one year of implantation of this device is 0.5%. The risk of retinal detachment for high myopes following implantation with the Visian MICL lens¹ is more than 10 times the risk without surgery, i.e., greater than 10 fold the background rate of retinal detachment for high myopes (>-3 D) 5.0% in myopes >-6 D and 0.8% to 7.5% in pseudophakic eyes with high axial myopia.

³ Visian MICL Clinical Trial

⁴ Refer to table below for details on Surgical Reinterventions.

⁴ There is no ISO historical rate for cumulative hyphema, raised IOP requiring treatment, iritis (after 1 week), CDVA loss ≥ 2 lines, corneal edema (after 1 week) and anisocoria.

* Percentage calculated as (n/N)*100.

Table 54: Visian TICL Lens Related Additional Surgery

	n/210*	% ¹
Visian TICL Lens Repositioning	1	0.5%
Visian TICL Lens Replacement (too long)	1	0.5%
Visian TICL Lens Removal (no ICL lens or IOL replacement)	3	1.4%
YAG Iridotomy**	3	1.4%
TOTAL	8	3.8%

* Total Eye Cohort (N = 210)

¹ Percentage calculated as (n/N)*100.

** Three cases (3/210, 1.4%) underwent an additional iridotomy. One of these was performed on the day of surgery because the surgeon felt the previous YAG procedure was inadequate. The IOP was 14 mmHg or less at all postoperative visits. The second case had an additional YAG iridotomy performed at 5 days postoperative to deepen the anterior chamber which was successful. This case was not associated with an increase in IOP. In the third case, the procedure was performed at 1 day postoperative to enlarge the preoperative iridotomy which was occluded by retained viscoelastic material, resulting in elevated IOP. Subsequent to the YAG procedure, the IOP returned to normal and stayed normal for the remainder of the follow-up.

Anterior subcapsular opacities, not all clinically significant, were observed postoperatively in six eyes (6/210, 2.9%). Two of these six eyes (2/210, 1.0%) had clinically significant cataracts. The remaining 4 subjects were asymptomatic with 20/16 or better CDVA and 20/25 or better UDVA at their last reported visit.

Decrease in Refractive Myopia and Cylinder

Reduction in refractive myopia and cylinder (manifest refraction spherical equivalent [MRSE] and cylinder) were the primary efficacy outcomes for the study. **Table 55** and **Table 57** provide MRSE and cylinder over time, and **Table 56** and **Table 58** provide a comparison between preoperative and 12 month MRSE and cylinder for the consistent cohort. The mean MRSE improved from -9.34 D preoperative to 0.03 D at the 12 month follow-up visit. There was a highly significant (p< 0.001) 1.43 D mean decrease in cylinder from preoperative to 12 months postoperative (paired t-test).

Table 55: MRSE by Visit – Visian TICL Study

	Preop	1 Week	1 Month	3 Months	6 Months	12 Months
N (eyes)	210	205	200	191	182	194
Mean (D)	-9.38	0.02	0.13	0.13	0.11	0.03
SD	2.67	0.45	0.43	0.39	0.49	0.46
Range (D)	-19.50 to -2.38	-1.50 to 1.38	-1.63 to 1.75	-1.25 to 1.25	-1.75 to 2.63	-2.25 to ±1.00

Table 56: MRSE – Preoperative vs. 12 Months (consistent cohort) – Visian TICL Study

	Preop	12 Months
N (eyes)	194	194
Mean (D)	-9.34	0.03
SD	2.63	0.46
Range (D)	-19.50 to -2.38	-2.25 to ±1.00

Table 57: Manifest Refraction Cylinder by Visit

Spherical Equivalent	Preop	1 Week	1 Month	3 Months	6 Months	12 Months
N (eyes)	210	205	200	191	182	194
Mean (D)	1.95	0.50	0.50	0.52	0.45	0.52
SD	0.84	0.54	0.49	0.49	0.45	0.48
Range (D)	1.00 to 4.00	0.00 to 3.00	0.00 to 3.00	0.00 to 3.00	0.00 to 2.00	0.00 to 3.00

Table 58: Manifest Refraction Cylinder: Preoperative vs. 12 Months (consistent cohort) – Visian TICL Study

	Preop	12 Months
N (eyes)	194	194
Mean (D)	1.95	0.52
SD	0.85	0.48
Range (D)	1.00 to 4.00	0.00 to 3.00

Visual Acuity

The visual acuities at 6 and 12 months are described in Table 59 – Table 61.

Table 59: UDVA (Eyes with Preoperative CDVA 20/20 or Better) – Visian TICL Study

	Preoperative N=173 Eyes	6 Months N=155 Eyes	12 Months N=159 Eyes
UDVA	n/173, % ¹	n/155, % ¹	n/159, % ¹
20/12.5 or better	0, 0%	41, 26.5%	40, 25.2%
20/16 or better	0, 0%	117, 75.5%	101, 63.5%
20/20 or better	0, 0%	140, 90.3%	142, 89.3%
20/40 or better	0, 0%	155, 100%	159, 100%
20/50 or worse	173, 100%	0, 0.0%	0, 0.0%
20/200 or worse	173, 100%	0, 0.0%	0, 0.0%

¹ Percentage calculated as (n/N)*100.

Table 60: CDVA (Eyes with Preoperative CDVA 20/20 or better) – Visian TICL Study

	6 Months N=155 Eyes	12 Months N=159 Eyes
CDVA	n/155, % ¹	n/159, % ¹
20/12.5 or better	71, 45.8%	72, 45.3%
20/16 or better	141, 91.0%	143, 89.9%
20/20 or better	155, 100%	159, 100%
20/25 or better	155, 100%	159, 100%
20/40 or better	155, 100%	159, 100%

¹ Percentage calculated as (n/N)*100.

Table 61: Comparison of Preoperative CDVA to 12 Month Postoperative UDVA – Visian TICL Study

	Preop CDVA N=193 Eyes n/N, % ¹	12 Month UDVA N=193 Eyes n/N, % ¹
20/12.5 or better	7, 3.6%	40, 20.7%
20/16 or better	79, 40.9%	104, 53.9%
20/20 or better	159, 82.4%	158, 81.9%
20/25 or better	181, 93.8%	175, 90.7%
20/32 or better	190, 98.4%	180, 93.3%
20/40 or better	193, 100.0%	184, 95.3%
20/80 or better	193, 100.0%	191, 99.0%
20/200 or better	193, 100.0%	193, 100.0%
Worse than 20/200	0, 0%	0, 0%

¹ Percentage calculated as (n/N)*100.

Predictability of Refraction

The MRSE of the refraction was predictable with 97.4% (189/194) of eyes achieving within ± 1.0 D from target at the 12 month examination.

Table 62: Accuracy of MRSE to Target – Visian TICL Study

	N=194 Eyes n/194, % ¹
Within ± 0.50 D	149, 76.8%
Within ± 1.0 D	189, 97.4%

¹ Percentage calculated as (n/N)*100.

The manifest cylinder was predictable with 92.3% (179/194) of eyes achieving within ± 1.0 D from target at the 12 month examination.

Table 63: Accuracy of Manifest Cylinder to Target (at the corneal plane) – Visian TICL Study

	12 Months N=194 Eyes n/194, % ¹
Within ± 0.50 D	134, 69.1%
Within ± 1.0 D	179, 92.3%

¹ Percentage calculated as (n/N)*100.

The effect of a temporal corneal incision on corneal toricity was analyzed. On average, implantation of the Visian TICL lens contributes less than 0.5 D of “with-the-rule” astigmatism to the net corneal toricity.

Stability

MRSE was stable with 99.4% (176/177) of eyes achieving less than or equal to ± 1.0 D of shift between 6 and 12 months after surgery.

Table 64: MRSE Change between Visits – Visian TICL Study

	1 Month to 3 Months N=184 Eyes	3 Months to 6 Months N=172 Eyes	6 Months to 12 Months N=177 Eyes
Change	n/184, %¹	n/172, %¹	n/177, %¹
Within ± 0.25 D	136, 73.9%	129, 75.0%	139, 78.5%
Within ± 0.50 D	169, 91.8%	159, 92.4%	167, 94.4%
Within ± 1.0 D	184, 100%	170, 98.8%	176, 99.4%
> 1.0 D	0, 0%	2, 1.2%	1, 0.6%
Mean Change	0.010	-0.009	0.081
SD	0.311	0.330	0.360
95% CI of the Mean	-0.04 to 0.05	-0.06 to 0.04	0.03 to 0.13

¹ Percentage calculated as (n/N)*100.

Manifest cylinder was stable with 97.2% (172/177) -98.8% (165/167) of eyes achieving less than or equal to ± 1.0 D of shift between 6 to 12 months after surgery, depending on analysis method.

Table 65: Manifest Cylinder Change Between Visits – Visian TICL Study

Analysis Group	Exam Interval	N (Eyes)	Within ± 0.5 D n/N, % ¹	Within ± 1.0 D n/N, % ¹	Mean Change for Interval [95% Confidence Interval]
Vector Stability	1 to 3 Months	184	143/184, 77.7%	179/184, 97.3%	0.26 D [0.23 to 0.3]
	3 to 6 Months	172	145/172, 84.3%	167/172, 97.1%	0.23 D [0.19 to 0.26]
	6 to 12 Months	177	141/177, 79.7%	172/177, 97.2%	0.26 D [0.22 to 0.29]
Vector Stability Consistent cohort	1 to 3 Months	167	130/167, 77.8%	162/167, 97.0%	0.26 D [0.23 to 0.3]
	3 to 6 Months		140/167, 83.8%	162/167, 97.0%	0.23 D [0.19 to 0.27]
	6 to 12 Months		134/167, 80.2%	163/167, 97.6%	0.24 D [0.21 to 0.28]
Stability of Absolute Cylinder	1 to 3 Months	184	154/184, 83.7%	181/184, 98.4%	0.00 D [-0.05 to 0.05]
	3 to 6 Months	172	153/172, 89.0%	170/172, 98.8%	-0.03 D [-0.08 to 0.01]
	6 to 12 Months	177	151/177, 85.3%	174/177, 98.3%	0.04 D [0 to 0.09]
Stability of Absolute Cylinder Consistent Cohort	1 to 3 Months	167	140/167, 83.8%	164/167, 98.2%	0.00 D [-0.05 to 0.05]
	3 to 6 Months		148/167, 88.6%	165/167, 98.8%	-0.03 D [-0.08 to 0.01]
	6 to 12 Months		143/167, 85.6%	165/167, 98.8%	0.03 D [-0.02 to 0.07]

¹ Percentage calculated as (n/N)*100.

Study investigators were asked to examine the patient at the slit lamp and estimate the orientation of the long axis of the Visian TICL lens based upon the alignment markings or haptic edges if visible. The lens orientation was then recorded in clock hours. For instance, if the lens was oriented exactly horizontally it would be recorded as at either 3:00 or 9:00 (clock hour position). Rotation was evaluated based upon the change in clock hour orientation of the Visian TICL lens postoperatively. A change of a half clock hour would represent 15 degrees of rotation and a change of a quarter clock hour would represent 7.5 degrees of rotation.

Table 66: Rotation of the Visian TICL Lens Between Visits (from direct observation of Visian TICL Lens)

	1 Day – 1 Week	1 Week – 1 Month	1 Month – 3 Months	3 Months – 6 Months	6 Months – 12 Months
N (Eyes)	121	155	148	136	140
Rotation	n/121, % ¹	n/155, % ¹	n/148, % ¹	n/136, % ¹	n/140, % ¹
≤ 5°	118, 97.5%	148, 95.5%	141, 95.3%	133, 97.8%	132, 94.3%
≤ 10°	121, 100%	155, 100%	147, 99.3%	135, 99.3%	137, 97.9%

¹ Percentage calculated as (n/N)*100.

Optical Visual Symptoms

A standardized subjective patient questionnaire was administered across all investigative sites to all subjects in the Visian TICL Lens Study Cohort preoperatively and after Visian TICL lens implantation. Study subjects' subjective assessments of ocular symptoms of glare, halos, double vision, night vision and night driving difficulties were evaluated for each eye at the preoperative and at the 3 and 12 month postoperative follow-up visits. Subjects were asked to grade the level of the specific ocular symptom in one of five categories: Absent, Mild, Moderate, Marked or Severe.

Table 67: Eyes with Symptoms Worse at 12 Months compared to Preoperative – Visian TICL Study

Visual Symptom	Worse at 12 Months than Preoperative n/N, % ¹
Glare	28/185, 15.1%
Halos	33/185, 17.8%
Double Vision	3/185, 1.6%
Night Vision	22/184, 11.9%
Night Driving Difficulties	24/182, 13.2%

¹ Percentage calculated as (n/N)*100.

NOTE: The questionnaire and methodology used to evaluate these subjective symptoms were not considered by the FDA to be validated.

Additional Clinical Outcomes - Visian TICL Study

Table 68 provides predictability of intended refraction (within ± 0.50 D and ± 1.0 D) for all eyes and by the level of preoperative refraction.

Table 68: Accuracy of MRSE vs. Intended Target* by Preoperative MRSE – Visian TICL Study

Lens Group	Exam Interval	N	Within ± 0.50 D n/N, % ¹	Within ± 1.0 D n/N, % ¹	Within ± 2.0 D n/N, % ¹
Study Cohort	1 Week	201	149/201, 74.1%	194/201, 96.5%	201/201, 100%
	1 Month	198	155/198, 78.3%	189/198, 95.5%	198/198, 100%
	3 Months	190	142/190, 74.7%	185/190, 97.4%	190/190, 100%
	6 Months	181	122/181, 67.4%	174/181, 96.1%	180/181, 99.4%
	12 Months	194	149/194, 76.8%	189/194, 97.4%	194/194, 100%
≤ -7 D Cohort	12 Months	33	28/33, 84.8%	33/33, 100%	33/33, 100%
> -7 to -10 D Cohort	12 Months	93	76/93, 81.7%	92/93, 98.9%	93/93, 100%
> -10 D to -15 D Cohort	12 Months	62	42/62, 67.7%	59/62, 95.2%	62/62, 100%
>-15 D Cohort	12 Months	6	3/6, 50.0%	5/6, 83.3%	6/6, 100%

* All Study Cohort Eyes

¹ Percentage calculated as (n/N)*100.

Table 69: Accuracy of Manifest Cylinder vs. Intended Target By Visit – Visian TICL Study

Lens Group ¹	Exam Interval	N (Eyes)	Within 0.25 D n/N ² , % ³	Within 0.50 D n/N ² , % ³	Within 1.00 D n/N ² , % ³	Within 2.00 D n/N ² , % ³
Study Cohort	Preop	210	0/210, 0%	0/210, 0%	43/210, 20.5%	134/210, 63.8%
	1 Week	205	92/201, 45.8%	128/201, 63.7%	184/201, 91.5%	198/201, 98.5%
	1 Month	200	84/198, 42.4%	128/198, 64.6%	180/198, 90.9%	197/198, 99.5%
	3 Months	191	77/190, 40.5%	123/190, 64.7%	174/190, 91.6%	186/190, 97.9%
	6 Months	182	87/181, 48.1%	128/181, 70.7%	167/181, 92.3%	181/181, 100%
	12 Months	194	78/194, 40.2%	127/194, 65.5%	177/194, 91.2%	193/194, 99.5%

¹ All Study Cohort Eyes

² Eyes with non-missing data

³ Percentage calculated as (n/N)*100.

Table 70: Percent Reduction of Absolute (non-vector) Cylinder Attempted ‘vs’ Achieved at the Spectacle Plane – Visian TICL Study

Preoperative Cylinder	N=194* Eyes	Percent Reduction of Absolute Cylinder		
		Mean	Range	[% CI]
ALL	n/194, %	77.8	-62.7 to 151.9	[73.9 to 81.6]
> 0.5 D to ≤ 1.0 D	39, 20.1%	75.1	-26.4 to 125.2	[65.4 to 84.8]
> 1.0 D to ≤ 2.0 D	86, 44.3%	71.4	-62.7 to 137.3	[64.9 to 77.9]
> 2.0 D to ≤ 3.0 D	45, 23.2%	87.1	44.8 to 151.9	[82.2 to 91.9]
> 3.0 D to ≤ 4.0 D	24, 12.4%	87.6	29.0 to 125.6	[80.3 to 95]

* All Study Cohort Eyes
Percentage calculated as (n/N)*100.

Table 71 shows the UDVA by the level of preoperative refraction for all eyes implanted that had a CDVA of 20/20 or better preoperatively.

Table 71: UDVA* Over Time and by Preoperative MRSE – Visian TICL Study

MRSE Group	Exam Interval	N Eyes	20/20 or Better n/N, % ¹	20/40 or Better n/N, % ¹
Study Cohort	1 Week	171	131/171, 76.6%	170/171, 99.4%
	1 Month	166	139/166, 83.7%	164/166, 98.8%
	3 Months	161	140/161, 87.0%	161/161, 100%
	6 Months	155	140/155, 90.3%	155/155, 100%
	12 Months	159	142/159, 89.3%	155/155, 100%
≤ -7 D	12 Months	33	31/33, 93.9%	32/33, 97.0%
> -7 D to -10 D	12 Months	93	78/93, 83.9%	91/93, 97.8%
> -10 D to -15 D	12Months	61	47/61, 77.0%	59/61, 96.7%
> -15 D	12Months	6	2/6, 33.3%	2/6, 33.3%

* In eyes with preoperative CDVA of 20/20 or better

¹ Percentage calculated as (n/N)*100.

Subjective Quality of Vision

A standardized subjective patient questionnaire was administered across all investigative sites to all subjects in the Visian TICL Lens Study preoperatively and after Visian TICL lens implantation. Study subject's subjective assessments of their quality of vision were evaluated for each eye at the preoperative and at the 3 and 12 month postoperative follow-up visits. Subjects were asked to rate their level of quality of vision in one of five categories: Excellent, Very Good, Good, Poor or Very Poor.

Table 72: Subjective Quality of Vision (All Eyes) – Visian TICL Study

Quality of Vision Grading	Preoperative N=210 n/210, % ¹	12 Months N=184 n/184, % ¹
Excellent /Very Good	135/210, 64.3%	174/184, 94.6%
Good	53/210, 25.2%	10/184, 5.4%
Poor/Very Poor	22/210, 10.5%	0/184, 0%

¹ Percentage calculated as (n/N)*100.

NOTE: The questionnaire and methodology used to evaluate these subjective symptoms were not considered by the FDA to be validated.

POST-APPROVAL STUDY OF THE VISIAN TICL LENS

This study was conducted to evaluate the previously approved parent Visian TICL lens (without the central port). As a condition of approval, the FDA required that STAAR conduct a study to evaluate the long-term clinical performance of the TICL. The first subject was enrolled in February of 2020 and the last completed follow-up visit occurred in December of 2024. The study design and results are described in this section.

Study Objective

The objective of this study was to evaluate the long-term clinical performance of the TICL with respect to rotational stability, refractive and visual outcomes, and ocular adverse event (AE) rates.

Study Design

The TICL lenses were evaluated in a prospective, multicenter, single-arm, two-year follow-up post-approval study. Subjects who met study eligibility criteria received unilateral or bilateral TICL implantation and were seen for 9 scheduled postoperative study visits (per eye) at approximately 1-6 hours, 1-2 days, 1 week (5-9 days), 1 month (3-5 weeks), 3 months (10-14 weeks), 6 months (21-26 weeks), 12 months (11-14 months), 18 months (17-21 months), and 24 months (23-27 months).

Study Eligibility Criteria

The subjects were eligible for study enrollment if they met the Visian TICL indications and did not meet the Visian TICL contraindications, in accordance with the product labeling.

Other inclusion criteria were:

- Correctable (CDVA) to at least 20/40 in the eye to be treated and absent of ocular pathology (except that myopic degeneration was allowed)
- Able to achieve a minimum dilated pupil size of 7 mm (via mydriatic and/or cycloplegic eye drops) in each eye to be treated
- Able to read, understand and provide written informed consent
- Able and willing to return for scheduled follow-up examinations after surgery

Other exclusion criteria were:

- History or clinical signs of iritis/uveitis
- Diabetic retinopathy
- History of previous eye surgery
- Serious (i.e., life threatening), acute, chronic or systemic, non-ophthalmic disease or illness that would have increased the operative risk, confounded the outcome(s) of the study or which may have precluded study completion
- Progressive sight threatening disease or other previous or current ocular conditions, other than myopia, that may have predisposed for future complications
- Ocular hypertension
- Insulin-dependent diabetes
- Pseudoexfoliation
- Pigment dispersion
- Conditions associated with fluctuation of hormones that could have led to refractive changes
- Involvement in another clinical study

Demographics and Baseline Characteristics

A total of 232 eyes of 125 subjects were enrolled at 8 clinical sites in the U.S.

Table 73: Demographics

Demographics	Subjects (N = 125)
Age, yrs	
Mean (SD)	35.7 (5.26) yrs
Range	21 to 45 yrs
Gender	n, %¹
Male	48, 38.4%
Female	77, 61.6%
Race	n, %¹
Caucasian	81, 64.8%
African American/Black	10, 8.0%
Asian	32, 25.6%
American Indian or Alaska Native	1, 0.8%
Other	1, 0.8%
Ethnicity	n, %¹
Hispanic or Latino	12, 9.6%
Not Hispanic or Latino	113, 90.4%

¹ Percentage calculated as (n/N)*100.

The distribution of preoperative cylinder (in the spectacle plane) for all eyes is shown in **Table 74**. The mean (SD) preoperative manifest refractive cylinder was -2.06 (0.86) with a range of -1.00 D to -4.00 D.

Table 74: Preoperative Manifest Refraction Cylinder

Preoperative Manifest Refraction Cylinder	Primary Eyes (N=125) n (n/N%)	All Eyes (N=232) n (n/N%)
Cylinder¹ (D)		
N	125	232
Mean (SD)	-2.26 (0.91)	-2.06 (0.86)
Median	-2.250	-1.875
Min, Max	-4.00, -1.00	-4.00, -1.00
Distribution (D)		
1.00	11 (8.8)	30 (12.9)
1.25 to 2.00	51 (40.8)	105 (45.3)
2.25 to 3.00	36 (28.8)	61 (26.3)
3.25 to 4.00	27 (21.6)	36 (15.5)

Percentage calculated as (n/N)*100

D=Diopter, N=number of eyes in the analysis set, SD=standard deviation

Cylinder reported in the negative scale

¹Preoperative cylinder in this study ranged -1.00 to -4.00.

As shown in **Table 75**, all available lengths of TICL lenses were implanted in the study, the majority 12.6 mm (39.7%) and 13.2 mm (47.4%), respectively.

Table 75: Implanted Lens Diameters

Lens Diameters	Primary Eyes (N=125) n (n/N%)	All Eyes (N=232) n (n/N%)
Model/Length (mm)		
TICL 12.1	5 (4.0)	9 (3.9)
TICL 12.6	51 (40.8)	92 (39.7)
TICL 13.2	56 (44.8)	110 (47.4)
TICL 13.7	13 (10.4)	21 (9.1)

Percentage calculated as (n/N)*100.

Study Endpoints

The primary endpoint, evaluated in primary eyes only was:

- rotation of less than or equal to five degrees between 18 and 24 months postoperative.

Secondary effectiveness and safety endpoints, evaluated in all treated (primary and fellow) eyes, were:

- absolute rotation between visits,
- absolute rotation <5 degrees, <10 degrees, <20 degrees, and <30 degrees from the intended orientation at each visit,
- absolute rotation from the intended orientation at each visit,
- postoperative manifest refraction spherical equivalent (MRSE) and cylinder at each visit,
- ocular AE rates.

Accountability

Table 76: Accountability – Primary Eyes

Eye Status	Total #	Op Visit (Day 0) n (%)	PO V2 (1-2 Days) n (%)	PO V3 (Day 5–9) n (%)	PO V4 (Day 21–35) n (%)	PO V5 (Day 70–98) n (%)	PO V6 (Day 147–182) n (%)	PO V7 (Day 330–420) n (%)	PO V8 (Day 510–630) n (%)	PO V9 (Day 690–810) n (%)
Eyes Treated (N)	125									
Available for Analysis		125 (100.0%)	124 (99.2%)	123 (98.4%)	122 (97.6%)	119 (95.2%)	116 (92.8%)	118 (94.4%)	114 (91.2%)	111 (88.8%)
Missing Eye/Data										
Discontinued		0 (0.0)	0 (0.0)	2 (1.6%)	2 (1.6%)	4 (3.2%)	4 (3.2%)	4 (3.2%)	4 (3.2%)	7 (5.6%)
Missing at scheduled visit but seen earlier/later ¹		0 (0.0)	0 (0.0)	0 (0.0)	2 (1.6%)	4 (3.2%)	5 (4.0%)	2 (1.6%)	2 (1.6%)	2 (1.6%)
Not seen but accounted for ²		0 (0.0)	1 (0.8%)	0 (0.0)	1 (0.8%)	2 (1.6%)	3 (2.4%)	1 (0.8%)	3 (2.4%)	0 (0.0)
Lost to follow-up		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	2 (1.6%)	2 (1.6%)	4 (3.2%)	7 (5.6%)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
% Accountability ⁴		125/125 (100.0%)	124/125 (99.2%)	123/123 (100.0%)	122/123 (99.2%)	119/121 (98.3%)	116/121 (95.9%)	118/121 (97.5%)	114/121 (94.2%)	111/118 (94.1%)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Not seen but accounted for: represents the total number of eyes that missed the visit but have not been discontinued/lost to follow-up.

³ Active: represents the total number of eyes that have not reached the time associated with the visit or are currently within the window for completion. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ % Accountability = [Available for Analysis / (N – [Discontinued-Active])].

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Table 77: Accountability – All Eyes

Eye Status	Total #	Op Visit (Day 0) n (%)	PO V2 (1-2 Days) n (%)	PO V3 (Day 5–9) n (%)	PO V4 (Day 21–35) n (%)	PO V5 (Day 70–98) n (%)	PO V6 (Day 147–182) n (%)	PO V7 (Day 330–420) n (%)	PO V8 (Day 510–630) n (%)	PO V9 (Day 690–810) n (%)
Eyes Treated (N)	232									
Available for Analysis		232 (100.0%)	230 (99.1%)	228 (98.3%)	227 (97.8%)	223 (96.1%)	217 (93.5%)	220 (94.8%)	212 (91.4%)	208 (89.7%)
Missing Eye/Data										
Discontinued		0 (0.0)	0 (0.0)	4 (1.7%)	4 (1.7%)	6 (2.6%)	6 (2.6%)	6 (2.6%)	6 (2.6%)	10 (4.3%)
Missing at scheduled visit but seen earlier/later ¹		0 (0.0)	0 (0.0)	0 (0.0)	4 (1.7%)	7 (3.0%)	9 (3.9%)	4 (1.7%)	4 (1.7%)	4 (1.7%)
Not seen but accounted for ²		0 (0.0)	2 (0.9%)	0 (0.0)	1 (0.4%)	3 (1.3%)	5 (2.2%)	2 (0.9%)	6 (2.6%)	0 (0.0)
Lost to follow-up		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	4 (1.7%)	4 (1.7%)	8 (3.4%)	14 (6.0%)
Active ³		0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)	0 (0.0)
% Accountability ⁴		232/232 (100.0%)	230/232 (99.1%)	228/228 (100.0%)	227/228 (99.6%)	223/226 (98.7%)	217/226 (96.0%)	220/226 (97.3%)	212/226 (93.8%)	208/222 (93.7%)

¹ Missing at scheduled visit but seen earlier/later: represents the total number of eyes that were seen outside the time window associated with the visit.

² Not seen but accounted for: represents the total number of eyes that missed the visit but have not been discontinued/lost to follow-up.

³ Active: represents the total number of eyes that have not reached the time associated with the visit or are currently within the window for completion. The investigation at the visit is considered complete when the number of active eyes is zero.

⁴ % Accountability = [Available for Analysis / (N – [Discontinued-Active])].

The denominator for percentages is the number of treated eyes. Percentage calculated as (n/N)*100.

Primary Study Endpoint - Effectiveness

The primary study effectiveness endpoint was considered met if at least 90% of the treated primary eyes exhibited rotation ≤ 5 degrees between 18-and-24 Months postoperative. Of the 104/108 (96.3%) primary eyes available for analysis at 24 Months, all (104/104, 100%) exhibited rotation ≤ 5 degrees during this period (Table 78).

Table 78: Absolute Lens Rotation $\leq 5^\circ$ Between 18 and 24 Months Postoperative – Primary Eyes

Primary Endpoint:	Eyes with absolute rotation $\leq 5^\circ$ between 18 and 24 months	
Primary Eyes (N = 108)	Eyes n (%) ¹	(95% CI) ³
$\leq 5^\circ$	104 (96.3%)	0.963 (0.9079 - 0.9898)
Missing	4 (3.7%) ²	-

¹ Percentage calculated as (n/N)*100.

² Four primary eyes of four subjects were missing for this analysis.

³ 95% Clopper-Pearson confidence interval for proportion of eyes available for analysis.

Safety Outcomes

The rates of ocular AEs in all treated eyes was a secondary safety endpoint in the study. A total of 141 ocular AEs were reported for 83 (35.8%) treated eyes of 53 subjects (**Table 79**). Of the 141 AEs included in **Table 79**, 6 events in 4 (1.7%) eyes of 2 subjects occurred prior to ICL surgery. In two eyes of one subject, non-patent PIs required repeat YAG treatment prior to ICL surgery. These same eyes also experienced iritis prior to ICL surgery. In another subject, one eye required a repeat YAG PI prior to ICL surgery, and the other eye experienced increased IOP. The remaining 135 AEs occurred after TICL implantation.

The most frequently reported postoperative AE was increased IOP with a total of 48 events occurring in 40 (17.2%) eyes of 28 subjects (**Table 79** number of events = 49, one instance of increased IOP occurred prior to implant). Most of the AEs for increased IOP were associated with OVD retention. Six instances of increased IOP in 4 eyes of 2 subjects were reported as serious and were associated with SAEs of narrowing of the anterior chamber angle that resulted in explantation of the study lenses. The next most frequently reported AEs were dry eye (10/232, 4.3%), residual refractive error (9/232, 3.9%), narrowing of the anterior chamber (6/232, 2.6%), and pupillary block and TMICL misalignment (each with 5/232, 2.2%, respectively).

Loss of CDVA of 2 lines or greater was reported in 2/232 eyes (0.9%). In one eye, CDVA at one month was 20/32 compared to 20/20 preoperatively as a result of lens dislocation. Lens removal surgery was performed a week later and CDVA improved to 20/25 one month post lens removal surgery prior to study exit. In the other eye with 2 line loss, the CDVA loss was associated with misalignment of the lens. CDVA measured 20/26 at 3 months after surgery compared to 20/16 before surgery, but improved on its own to 20/20 by 12 months, and to 20/17 by 24 months after surgery. The lens misalignment in this case did not require intervention.

All other ocular AEs that occurred during this clinical trial were previously anticipated in nature, severity, and frequency based on prior clinical studies as well as the published literature regarding the ICL family of lenses and refractive surgery.

Table 79: Ocular Adverse Events

Cumulative Ocular AEs	All Eyes (N=232)	
	Eyes n(%) ¹	Events n
Eye experienced any ocular AE ²	83 (35.8%)	141
Intraocular pressure increased	40 (17.2%)	49 ³
Dry eye	10 (4.3%)	10
Residual refractive error	9 (3.9%)	9
Narrow anterior chamber angle	6 (2.6%)	6
Pupillary block	5 (2.2%)	5
Residual cylindrical refractive error (TMICL misalignment)	5 (2.2%)	5
Iris transillumination defect	4 (1.7%)	4
Punctate keratitis	4 (1.7%)	4
Corneal erosion	3 (1.3%)	3
Ocular Allergies	3 (1.3%)	3
Conjunctivitis viral	3 (1.3%)	3
Corneal abrasion	3 (1.3%)	3
Non-patent YAG PI ⁴	3 (1.3%)	3
Chalazion	1 (0.4%)	2
Eye pain	1 (0.4%)	2
Increase in cylindrical refractive error (TMICL rotation)	1 (0.4%)	2
Conjunctivitis allergic	2 (0.9%)	2
Giant papillary conjunctivitis	2 (0.9%)	2
Glare/Halo	2 (0.9%)	2
Increased TM pigment	2 (0.9%)	2
Iritis ⁴	2 (0.9%)	2
Meibomian gland dysfunction	2 (0.9%)	2
Visual acuity reduced	2 (0.9%)	2
Anterior subcapsular opacity	2 (0.9%)	2
Eye irritation	1 (0.4%)	1
Posterior subcapsular cataract	1 (0.4%)	1
Residual Astigmatism	1 (0.4%)	1
Subjective visual disturbance	1 (0.4%)	1
Vitreous detachment	1 (0.4%)	1
Intraocular lens dislocation	1 (0.4%)	1
Anterior chamber cell/flare	1 (0.4%)	1
Episcleritis	1 (0.4%)	1
Ocular discomfort	1 (0.4%)	1
Pigment on ICL	1 (0.4%)	1
Pupil fixed	1 (0.4%)	1
Drug hypersensitivity	1 (0.4%)	1

¹ Percentage calculated as (n/N)*100.

² Only the first incidence of an event is counted for any given eye.

³ One AE of increased IOP in one eye occurred prior to ICL surgery.

⁴ These AEs occurred prior to ICL surgery.

Ocular SAEs and SSIs

A total of 40 ocular SAEs were reported in this study (**Table 80**). The most common reported event was intraocular pressure increased with 11 (11/40, 27.5%) followed by residual refractive error 10 (10/40, 25%) (SAEs of residual refractive error were treated with PRK or LASIK), narrowing of the anterior chamber with 6 (6/40, 15%), pupillary block with 5 (5/40, 12.5%), and 4 (4/40, 10%) events of residual cylindrical refractive error (TMICL misalignment).

Table 80: Ocular Serious Adverse Events

Cumulative Ocular SAEs	All Eyes (N=232)	
	Eyes ¹ n (n/N%)	Events n
Eye experienced any ocular SAE	29 (12.5%)	40
Intraocular pressure increased	9 (3.9%)	11
Residual refractive error	10 (4.3%)	10
Narrow anterior chamber angle	6 (2.6%)	6
Pupillary block	5 (2.2%)	5
Residual cylindrical refractive error (due to TMICL misalignment)	4 (1.7%)	4
Eye pain	1 (0.4%)	1
Increase in cylindrical refractive error (due to TMICL rotation)	1 (0.4%)	1
Posterior subcapsular cataract	1 (0.4%)	1
Intraocular lens dislocation	1 (0.4%)	1

Percentage calculated as (n/N)*100.

¹ Only the first incidence of an event is counted for any given eye.

Secondary surgical interventions were reported for 28 (28/232, 12.1%) eyes of 21 subjects (**Table 81**). In 10 eyes (10/232, 4.3%) of 7 subjects, LASIK enhancement procedures were performed to address complaints of blurred vision (UDVA ranged from 20/16 to 20/32) after TICL surgery. Surgery to remove the lens was performed in 7 eyes (7/232, 3.0%) of 5 subjects due to: narrowing of the anterior chamber in 3 eyes of 2 subjects; increased IOP, eye pain, and TICL dislocation in 1 eye of each of 3 subjects; and in 1 eye of 1 subject prior to subsequent cataract surgery. Surgery to reposition the lens was performed in 5 eyes (5/232, 2.2%) of 4 subjects: to address residual refractive error in 3 eyes of 2 subjects; and to correct lens misalignment in 2 eyes of 2 subjects. Iridotomy/iridectomy was performed in 4 eyes (4/232, 1.7%) of 4 subjects: to treat pupillary block in 3 eyes of 3 subjects; and as treatment for increased IOP due to non-patent YAG PI in 1 eye of 1 subject. In 3 eyes (3/232, 1.3%) of 2 subjects, PRK enhancement procedures were performed to address residual refractive error (UDVA ranged from 20/20 to 20/40) after TICL surgery. In 2 eyes (2/232, 0.9%) of 1 subject a PI enlargement was performed to treat increased IOP and in 2 eyes (2/232, 0.9%) of another subject, lens exchange procedure was performed to treat narrowing of the anterior chamber angle. Cataract surgery was performed in 1 eye (1/232, 0.4%) of 1 subject to treat clinically significant cataract. All events reported in the study that required SSIs resolved without sequelae.

Table 81: Secondary Surgical Interventions

Secondary Surgical Interventions	All Eyes (N=232)	
	Eyes ¹ n (n/N%)	Events n
Eye experienced any AE requiring SSI²	28 (12.1)	34
LASIK	10 (4.3)	10
Surgery to remove lens	7 (3.0)	7
Surgery to reposition lens	5 (2.2)	5
Iridotomy/Iridectomy for pupillary block	4 (1.7)	4
PRK	3 (1.3)	3
Peripheral Iridotomy Enlargement	2 (0.9)	2
Surgery to exchange lens	2 (0.9)	2
Cataract Surgery	1 (0.4)	1

Percentage calculated as (n/N)*100.

¹ Only the first incidence of an event is counted for any given eye.

² Six eyes each experienced 2 events: ICL repositioning and LASIK in one eye, lens removal and cataract surgery in one eye, lens repositioning and PRK in 2 eyes of one subject and PI enlargement and lens removal in 2 eyes of one subject.

Secondary Effectiveness Endpoints

Absolute Lens Rotation

The absolute lens rotation between visits in all eyes was a secondary effectiveness endpoint in the study. Mean (SD) absolute lens rotation of all eyes between the 18- and 24-Month postoperative visits was 0.2° (0.19°, **Table 82**). No eye rotated ≥ 5° between 18-to-24 months postoperative.

Table 82: Absolute Lens Rotation at Each Visit – All Eyes

	PO V2 (Day 1-2)	PO V3 (Day 5-9)	PO V4 (Day 21-35)	PO V5 (Day 70-98)	PO V6 (Day 147-182)	PO V7 (Day 330-420)	PO V8 (Day 510-630)	PO V9 (Day 690-810)
Total # (N)	230	228	227	223	217	220	212	208
Rotation from Op Visit								
Available for Analysis	210	213	218	215	211	210	204	200
Mean (SD)	1.3 (1.90)	1.5 (2.64)	1.9 (6.21)	1.6 (2.80)	2.0 (6.33)	2.0 (6.29)	2.1 (6.39)	1.6 (2.86)
Median	0.9	1.1	1.1	1.0	1.1	1.1	1.1	1.0
Range	0.0 to 18.8	0.0 to 34.9	0.0 to 83.7	0.0 to 35.5	0.0 to 83.0	0.0 to 83.2	0.0 to 81.9	0.0 to 35.4
Rotation from PO V2								
Available for Analysis	-	207	208	205	202	201	198	194
Mean (SD)	-	1.0 (2.61)	1.4 (6.36)	1.0 (2.23)	1.5 (6.20)	1.4 (6.15)	1.5 (6.20)	1.0 (2.03)
Median	-	0.6	0.5	0.6	0.6	0.7	0.6	0.6
Range	-	0.0 to 27.7	0.0 to 83.9	0.0 to 25.1	0.0 to 83.2	0.0 to 83.5	0.0 to 82.1	0.0 to 25.0
Rotation from PO V3								
Available for Analysis	-	-	218	215	209	208	201	199
Mean (SD)	-	-	1.0 (5.81)	0.7 (1.31)	1.2 (5.90)	1.1 (5.85)	1.2 (5.94)	0.7 (0.90)
Median	-	-	0.4	0.4	0.4	0.5	0.4	0.5
Range	-	-	0.0 to 84.6	0.0 to 15.0	0.0 to 83.9	0.0 to 84.1	0.0 to 82.8	0.0 to 8.1
Rotation from PO V4								
Available for Analysis	-	-	-	220	215	213	206	203
Mean (SD)	-	-	-	0.5 (0.82)	0.6 (0.85)	0.7 (0.87)	0.6 (0.90)	0.6 (0.89)
Median	-	-	-	0.4	0.4	0.4	0.4	0.4
Range	-	-	-	0.0 to 9.1	0.0 to 8.7	0.0 to 9.1	0.0 to 9.3	0.0 to 9.4
Rotation from PO V5								
Available for Analysis	-	-	-	-	214	211	204	201
Mean (SD)	-	-	-	-	0.5 (0.42)	0.5 (0.47)	0.5 (0.46)	0.5 (0.47)
Median	-	-	-	-	0.3	0.4	0.4	0.4
Range	-	-	-	-	0.0 to 3.1	0.0 to 2.4	0.0 to 3.2	0.0 to 2.9
Rotation from PO V6								
Available for Analysis	-	-	-	-	-	211	204	201
Mean (SD)	-	-	-	-	-	0.4 (0.39)	0.5 (0.39)	0.4 (0.40)
Median	-	-	-	-	-	0.3	0.3	0.3
Range	-	-	-	-	-	0.0 to 2.2	0.0 to 2.3	0.0 to 2.2
Rotation from PO V7								
Available for Analysis	-	-	-	-	-	-	203	203
Mean (SD)	-	-	-	-	-	-	0.3 (0.28)	0.3 (0.24)
Median	-	-	-	-	-	-	0.2	0.2
Range	-	-	-	-	-	-	0.0 to 1.7	0.0 to 1.4
Rotation from PO V8								
Available for Analysis	-	-	-	-	-	-	-	197
Mean (SD)	-	-	-	-	-	-	-	0.2 (0.19)
Median	-	-	-	-	-	-	-	0.2
Range	-	-	-	-	-	-	-	0.0 to 0.9

The absolute lens rotation from intended orientation and the absolute rotation <5 degrees, <10 degrees, <20 degrees, and <30 degrees from the intended orientation at each visit for all eyes were additional secondary effectiveness endpoints in the study. Lens rotation from the intended axis of greater than 5° occurred in 8 primary eyes (6.4%) of 8 subjects (one subject experienced 2 events of significant rotation in one eye), and 7 fellow eyes (6.5%) of 7 subjects (**Table 83**).

Table 83: Absolute Lens Rotation from Intended Orientation at Each Visit – All Eyes

Absolute Rotation from Intended Orientation	Op Visit (Day 0)	PO V2 (Day 1-2) (N=230)	PO V3 (Day 5-9) (N=228)	PO V4 (Day 21-35) (N=227)	PO V5 (Day 70-98) (N=223)	PO V6 (Day 147-182) (N=217)	PO V7 (Day 330-420) (N=220)	PO V8 (Day 510-630) (N=212)	PO V9 (Day 690-810) (N=208)
n (%)	232 (100%)	210 (91.3%)	213 (93.4%)	218 (96%)	215 (96.4%)	211 (97.2%)	210 (95.5%)	204 (96.2%)	200 (96.2%)
Mean (SD)	0.3 (2.8)	1.4 (2.46)	1.6 (3.06)	2.2 (6.77)	1.8 (3.90)	2.3 (6.89)	2.3 (6.84)	2.4 (6.94)	1.8 (3.36)
Median	0.0	0.9	1.1	1.1	1.0	1.1	1.1	1.1	1.0
Range	0.0 to 34.0	0.0 to 24.1	0.0 to 34.9	0.0 to 83.7	0.0 to 35.5	0.0 to 83.0	0.0 to 83.2	0.0 to 81.9	0.0 to 35.4
	n (%)	n (%)	n (%)	n (%)	n (%)	n (%)	n (%)	n (%)	n (%)
<5°	230 (99.1%)	205 (89.1%)	207 (90.8%)	210 (92.5%)	206 (92.4%)	199 (91.7%)	200 (90.9%)	193 (91.0%)	192 (92.3%)
<10°	230 (99.1%)	206 (89.6%)	211 (92.5%)	213 (93.8%)	211 (94.6%)	205 (94.5%)	205 (93.2%)	198 (93.4%)	197 (94.7%)
<20°	230 (99.1%)	209 (90.9%)	211 (92.5%)	214 (94.3%)	212 (95.1%)	207 (95.4%)	206 (93.6%)	200 (94.3%)	198 (95.2%)
<30°	231 (99.6%)	210 (91.3%)	212 (93.0%)	215 (94.7%)	213 (95.5%)	208 (95.9%)	207 (94.1%)	201 (94.8%)	199 (95.8%)
≥30°	1 (0.4%)	0 (0.0)	1 (0.4%)	3 (1.3%)	2 (0.9%)	3 (1.4%)	3 (1.4%)	3 (1.4%)	1 (0.5%)
Missing	0 (0.0)	20 (8.7%)	15 (6.6%)	9 (4.0%)	8 (3.6%)	6 (2.8%)	10 (4.5%)	8 (3.8%)	8 (3.8%)

¹ Percentage calculated as (n/N)*100.

Postoperative MRSE and Cylinder

Postoperative MRSE and cylinder at each visit was a secondary effectiveness endpoint in this study. The MRSE and cylinder improved from -8.79 D and -2.06 D preoperative to -0.20 D and -0.54 D, respectively, at the Month 1 visit and remained stable throughout the study (**Table 84**).

Table 84: Manifest Refraction by Visit – All Eyes

	Preop (Day -120 to -1)	PO V4 (Day 21 – 35)	PO V5 (Day 70 – 98)	PO V6 (Day 147 – 182)	PO V7 (Day 330-420)	PO V8 (Day 510-630)	PO V9 (Day 690-810)
All Eyes (N)	232	227	223	217	220	212	208
Cylinder (D)							
n	232	226	223	217	220	212	208
Mean (SD)	-2.06 (0.86)	-0.54 (0.42)	-0.55 (0.43)	-0.50 (0.44)	-0.47 (0.45)	-0.43 (0.44)	-0.41 (0.44)
Median	-1.88	-0.50	-0.50	-0.50	-0.50	-0.25	-0.25
Min, Max	-4.00, -1.00	-3.00, 0.00	-2.25, 0.00	-2.75, 0.00	-2.75, 0.00	-3.00, 0.00	-3.00, 0.00
Missing	0	1	0	0	0	0	0
MRSE (D)							
n	232	226	223	217	220	212	208
Mean (SD)	-8.79 (2.41)	-0.20 (0.38)	-0.19 (0.36)	-0.21 (0.36)	-0.23 (0.36)	-0.21 (0.32)	-0.19 (0.31)
Median	-8.750	-0.125	-0.250	-0.125	-0.250	-0.125	-0.125
Min, Max	-15.00, -3.50	-1.25, 1.00	-1.38, 0.88	-1.50, 0.63	-1.63, 0.63	-1.38, 0.75	-1.38, 0.75
Missing	0	1	0	0	0	0	0

ADDITIONAL CLINICAL OUTCOMES

Accuracy to Target - MRSE

Table 85: MRSE within ± 0.50 D and ± 1.00 D of Target through Month 24 – All Eyes

	Month 6 (PO V6, Day 147 – 182) All Eyes (N=217)		Month 12 (PO V7, Day 330 – 420) All Eyes (N=220)		Month 24 (PO V9, Day 690 – 810) All Eyes (N=208)	
	Eyes n (n/N%)	Proportion (95% CI) ¹	Eyes n (n/N%)	Proportion (95% CI)	Eyes n (n/N%)	Proportion (95% CI)
MRSE within ± 1.00 D of target	210 (96.8%)	96.774 (93.47-98.69)	214 (97.3%)	97.273 (94.16-98.99)	205 (98.6%)	98.558 (95.84-99.70)
MRSE within ± 0.50 D of target	170 (78.3%)	78.341 (72.26-83.63)	165 (75%)	75.000 (68.74-80.58)	165 (79.3%)	79.327 (73.18-84.62)

¹ 95% Clopper-Pearson confidence interval for proportion of eyes available for analysis.

Percentage calculated as (n/N)*100.

Visual Acuity

At Month 24, 99.3% of eyes achieved postoperative UDVA 20/40 or better and 72.9% achieved 20/20 or better (**Table 86**). Similar results are shown at the Month 6 and Month 12 postoperative visits.

Table 86: Postoperative UDVA through Month 24 among Subjects with CDVA 20/20 or better at Preoperative Visit - All Eyes

	Month 6 (PO V6, Day 147 – 182) All Eyes (N=217)		Month 12 (PO V7, Day 330 – 420) All Eyes (N=220)		Month 24 (PO V9, Day 690 – 810) All Eyes (N=208)	
	Eyes n (n/N%)	Proportion (95% CI) ¹	Eyes n (n/N%)	Proportion (95% CI)	Eyes n (n/N%)	Proportion (95% CI)
UDVA 20/20 or better	110/148 (74.3%)	74.324 (66.50–81.15)	119/152 (78.3%)	78.289 (70.88–84.56)	105/144 (72.9%)	72.917 (64.89–79.98)
UDVA 20/40 or better	147/148 (99.3%)	99.324 (96.29–99.98)	150/152 (98.7%)	98.684 (95.33–99.84)	143/144 (99.3%)	99.306 (96.19–99.98)

¹ 95% Clopper-Pearson confidence interval for proportion of eyes available for analysis.

Percentage calculated as (n/N)*100.

At Month 24, 100% of eyes achieved postoperative CDVA 20/40 or better and 99.3% achieved 20/20 or better (**Table 87**). Similar results are shown at the Month 6 and Month 12 postoperative visits.

Table 87: Postoperative CDVA through Month 24 among Subjects with CDVA 20/20 or better at Preoperative Visit - All Eyes

	Month 6 (PO V6, Day 147 – 182) All Eyes (N=217)		Month 12 (PO V7, Day 330 – 420) All Eyes (N=220)		Month 24 (PO V9, Day 690 – 810) All Eyes (N=208)	
	Eyes n (n/N%)	Proportion (95% CI) ¹	Eyes n (n/N%)	Proportion (95% CI)	Eyes n (n/N%)	Proportion (95% CI)
CDVA 20/20 or better	146/148 (98.6%)	98.649 (95.20–99.84)	149/152 (98.0%)	98.026 (94.34–99.59)	143/144 (99.3%)	99.306 (96.19–99.98)
CDVA 20/40 or better	148/148 (100.0%)	100.00 (97.54–100.00)	152/152 (100.0%)	100.00 (97.60–100.00)	144/144 (100.00%)	100.00 (97.47–100.00)

¹ 95% Clopper-Pearson confidence interval for proportion of eyes available for analysis.

Percentage calculated as (n/N)*100.

Study Strengths and Limitations

Strengths of this study include the following observations:

1. A central independent reading center was used for analysis of retroillumination photographs using objective methods.
2. The rotational stability of the TICL between 18 and 24 months postoperative and correction of manifest cylinder confirm the effectiveness of the TICL established in the PMA study of the TICL.
3. This study confirmed the results of prior clinical studies of the MICL and TICL showing well-preserved CDVA, and high levels of improved UDVA.
4. No new concerns were raised regarding the performance and established safety of the toric ICL platform.

Limitations of this study include the following observations:

1. Endothelial cell density assessments may have added to the robust safety profile of the lens.
2. This study did not include the EVO Toric ICL, the TICL model that includes the central hole.

PRE-APPROVAL PMA CLINICAL TRIAL AND RESULTS - VISIAN MICL LENS FOR MYOPIA

The Visian MICL lens was evaluated in a prospective nonrandomized study of 526 eyes of 294 subjects, 470 of which were followed for 1 year and 369 followed for 3 years. Demographics for the Study Cohort are presented in the following table:

Table 88: Demographics – Visian MICL Study

N=526 Eyes (294 Subjects)

Age	
Average	36.55 ±5.8 years
Range	22 to 45 years
Race	
N (294)	n, % ¹
Black	6, 2.0%
Caucasian	249, 84.7%
Hispanic	23, 7.8%
Other	16, 5.4%
Gender	
Female	178, 60.5%
Male	116, 39.5%

¹ Percentage calculated as (n/N)*100.

In the study, surgeons supplied the following parameters to STAAR: manifest refraction – sphere, cylinder, axis; back vertex distance in millimeters; ACD in millimeters (posterior surface of the cornea to the anterior surface of the crystalline lens); and corneal thickness in millimeters. STAAR calculated the appropriate Visian MICL lens power using proprietary software.

Adverse Events

A total of 526 eyes of 294 subjects were evaluated in the clinical trial to determine the safety of the Visian MICL lens. Anterior subcapsular opacities, not all clinically significant, were observed postoperatively in 14/526 eyes (2.7%). An increase in postoperative cylinder >2 D at 3 years from surgery was present in 0.4% (2/256) of eyes. Loss of best corrected visual acuity (CDVA) >2 lines occurred in 4/526 eyes (0.8%) and a 2 line loss in 6/526 eyes (1.2%).

The AEs/complications experienced during the clinical study of the Visian MICL lens (between 1 and 36 months) included 3 retinal detachments (3/526, 0.6%), 2 cases of glaucoma (2/526, 0.4%), clinically significant cataract (2 anterior (2/526, 0.4%); 5 nuclear (5/526, 1%)), 1 case of elevated IOP >25 mmHg / >10 mmHg change from baseline at last visit (1/526, 0.2%), 1 macular hemorrhage (1/526, 0.2%) and 1 subretinal hemorrhage (1/526, 0.2%).

There were 20 cases of raised IOP requiring secondary surgical intervention in the early postoperative time period. Of these, 17 eyes were treated with YAG laser iridotomy for pupillary block, and 3 eyes were treated with repeat irrigation and aspiration for removal of retained viscoelastic. There were 16 cases of secondary surgical intervention for Visian MICL lens removal, replacement, or repositioning. In addition, most patients experienced some degree of endothelial cell loss after Visian MICL lens implantation.

Incidence of AEs/complications (compared with the FDA Grid for cataract extraction and posterior chamber IOL implantation) and incidence of surgical reinterventions are shown in **Table 89** and **Table 90**.

Table 89: Adverse Events – Visian MICL Study

Adverse Event	Cumulative %* (n/N)	FDA Grid %	Persistent (36 Mo) %* (n/N)	FDA Grid %
Endophthalmitis	0% (0/526)	0.1%	0% (0/526)	---
Hyphema	0% (0/526)	2.2%	0% (0/526)	---
Hypopyon	0% (0/526)	0.3%	0% (0/526)	---
IOL Dislocation	0% (0/526)	0.1%	0% (0/526)	---
Cystoid Macular Edema	0% (0/526)	3.0%	0% (0/526)	0.5%
Pupillary Block	3.2% (17/526)	0.1%	0% (0/526)	---
Retinal Detachment ¹	0.6% (3/526)	0.3%	0% (0/526)	---
Surgical Reintervention ²	6.8% (36/526)	0.8%	0% (0/526)	---
Corneal Edema (after 1 week)	0% (0/526)	---	0% (0/526)	0.3%
Iritis ³ (after 1 week)	0% (0/526)	---	0% (0/526)	0.3%
Surgical Treatments Not Monitored in FDA Grid				
Refractive Procedures	3.9% (20/526)	---	---	---
Iris Prolapse Repair	0.2% (1/526)	---	0% (0/526)	---

* Study percentage calculated as (n/N)*100.

¹ Comparison should be made to literature for retinal detachment rates for high myopia. Retinal detachment rates increase with increasing myopia. The risk of retinal detachment within one year of implantation of this device is 0.2%. The risk of retinal detachment for high myopes following implantation is more than 10 times the risk without surgery, i.e., greater than 10 fold the background rate of retinal detachment for high myopes (>-3 D) 5.0% in myopes >-6 D and 0.8% to 7.5% in pseudophakic eyes with high axial myopia.

Ogawa A, Tanaka, M. The relationship between refractive errors and retinal detachment, Jpn J Ophthalmol 32;310:1988.

Dellone-Larkin G, Dellone CA. Retinal detachment. Available at: <http://www.emedicine.com/emerg/topic504.htm>

Jacobi F, Hessemer V. Pseudophakic retinal detachment in high axial myopia. J Cat Ref Surg 23; 1095:1997. Refractive procedures include: AK and LASIK

² Refer to table below for details on Surgical Reinterventions.

³ There is no FDA Grid Rate for cumulative iritis.

Surgical reinterventions (see **Table 90** below) were not shown to have an impact on efficacy. Surgical reinterventions occurred in 6.8% (36/526) of cases.

Table 90: Visian MICL Lens Related Additional Surgery

	n	% ^{1*}
Visian MICL Lens Repositioning	4	0.8%
Visian MICL Lens Replacement, then Removal	1	0.2%
Visian MICL Lens Replacement	8	1.5%
Visian MICL Lens Removal	3	0.6%
Raised IOP Requiring Surgery	20	3.8%
TOTAL	36	6.8%

¹ Total Study Cohort (N = 526)

* Percentage calculated as (n/N)*100.

Refer to the Section “POST-APPROVAL CONTINUATION OF THE VISIAN MICL LENS CLINICAL STUDY” for a detailed discussion of AEs and complications that occurred in the PMA study cohort from day of surgery throughout the long-term post-approval phase of the study.

Visual Acuity

The postoperative results demonstrated that the Visian MICL lens can provide full correction for high myopia up to -15 D and only partial correction up to -20 D. The visual acuities at 12 and 36 months are described in the following tables:

Table 91: UDVA – Visian MICL Study (Where emmetropia was the goal (±0.50 D) and Preoperative CDVA was 20/20 or better)

N	12 Months 240 n, % ¹	36 Months 189 n, % ¹
20/20 or better	157, 65.4%	112, 59.3%
20/40 or better	232, 96.7%	179, 94.7%
20/80 or better	239, 99.6%	187, 98.9%
Worse than 20/80	1, 0.4%	2, 1.1%

¹ Percentage calculated as (n/N)*100.

Table 92: CDVA – Visian MICL Study (Eyes with Preoperative CDVA 20/20 or better)

N	12 Months 321 n, % ¹	36 Months 253 n, % ¹
20/20 or better	307, 95.6%	244, 96.4%
20/25 or better	320, 99.7%	253, 100%
20/40 or better	321, 100%	253, 100%

¹ Percentage calculated as (n/N)*100.

Predictability of Refraction

The refraction was predictable with 91.6% (417/455) of subjects achieving ±1.0 D from target at the 12 month examination.

Table 93: Spherical Equivalent (Target Variance) Distribution – Visian MICL Study

N	12 Months 455 n, % ¹	36 Months 363 n, % ¹
±0.50 D	314, 69%	248, 68.3%
±1.0 D	417, 91.6%	325, 89.5%

¹ Percentage calculated as (n/N)*100.

Stability

The refraction was stable with 97.6% (329/337) of eyes achieving less than or equal to ±1.0 D of shift at 36 months.

Table 94: MRSE Change between Visits – Visian MICL Study

N	6-12 Months 424 n, % ¹	12-24 Months 413 n, % ¹	24-36 Months 337 n, % ¹
±0.25 D	320, 75.5%	317, 76.8%	253, 75.1%
±0.5 D	386, 91.0%	371, 89.8%	304, 90.2%
±1.0 D	414, 97.6%	403, 97.6%	329, 97.6%
>1.0 D	10, 2.4%	10, 2.4%	8, 2.4%

¹ Percentage calculated as (n/N)*100.

Optical Visual Symptoms

Table 95 reports the subjective optical visual symptoms reported by subjects during this clinical study after Visian MICL lens implantation compared to before the Visian MICL surgery:

Table 95: Subjective Symptoms – Visian MICL Study

Symptom	Improved/No Change at 36 Months n/N (%) ¹
Glare	317/351 (90.4%)
Halos	310/350 (88.5%)
Double Vision	345/351 (98.3%)
Night Vision	308/350 (88.0%)
Night Driving Difficulties	301/335 (89.8%)

¹ Percentage calculated as (n/N)*100.

Additional Clinical Outcomes

Table 96 provides predictability of intended refraction (± 0.50 D and ± 1.0 D) for all eyes and by the level of preoperative refraction.

Table 96: MRSE vs. Intended Target¹ by Preoperative MRSE – Visian MICL Study

Lens Group	Exam Interval	N	± 0.5 D n, % ¹	± 1.0 D n, % ¹	± 2.0 D n, % ¹
Study Cohort	1 Week	501	324, 64.7%	438, 87.4%	487, 97.2%
	1 Month	506	344, 68.0%	445, 87.9%	495, 97.8%
	3 Months	485	310, 63.9%	430, 88.7%	475, 97.9%
	6 Months	479	320, 66.8%	426, 88.9%	470, 98.1%
	12 Months	455	308, 67.7%	411, 90.3%	447, 98.2%
	24 Months	443	293, 66.1%	399, 90.1%	434, 98.0%
	36 Months	363	245, 67.5%	320, 88.2%	356, 98.1%
New Calculation Method ³	36 Months	363	254, 70.0%	324, 89.3%	357, 98.3%
≤ -7 D Cohort	36 Months	72	61, 84.7%	70, 97.2%	72, 100%
New Calculation Method ³	36 Months	72	62, 86.1%	70, 97.2%	72, 100%
> -7 D to -10 D Cohort ²	36 Months	131	93, 71.0%	122, 93.1%	131, 100%
New Calculation Method ³	36 Months	131	92, 70.2%	121, 92.4%	131, 100%
> -10 D to -15 D Cohort	36 Months	130	84, 64.6%	112, 86.2%	128, 98.5%
New Calculation Method ³	36 Months	130	91, 70.0%	115, 88.5%	129, 99.2%
> -15 D Cohort	36 Months	30	7, 23.3%	16, 53.3%	25, 83.3%
New Calculation Method ³	36 Months	30	9, 30.0%	18, 60.0%	25, 83.3%

¹ All Study Cohort Eyes

² Note % lower with new Power Calculation Method

³ The new calculation method was used to correct for a change in power labeling to allow standard phakic IOL power formulas to be used without modification. It is a theoretical calculation only.

Percentage calculated as $(n/N) \times 100$.

Table 97 shows the UDVA for all eyes and by the level of preoperative refraction for all eyes implanted that were targeted for emmetropia and had a CDVA of 20/20 or better preoperatively.

Table 97: UDVA¹ by Preoperative MRSE – Visian MICL Study

Lens Group	Exam Interval	n	20/20 or Better n, % ¹	20/40 or Better n, % ¹
Study Cohort	1 Week	259	129, 49.8%	238, 91.9%
	1 Month	262	148, 56.5%	249, 95.0%
	3 Months	251	160, 63.7%	242, 96.4%
	6 Months	248	171, 60.9%	242, 96.4%
	12 Months	240	171, 65.4%	232, 96.7%
	24 Months	228	136, 59.6%	213, 93.4%
	36 Months	189	112, 59.3%	179, 94.7%
≤ -7 D	36 Months	58	42, 72.4%	57, 98.3%
> -7 D to -10 D	36 Months	83	52, 62.7%	77, 92.8%
> -10 D to -15 D	36 Months	48	18, 37.5%	45, 93.8%
> -15 D	36 Months	0	NA ²	NA ²

¹ Eyes with preoperative CDVA 20/20 or better and emmetropia targeted correction

² No Eyes > -15 D group with this preoperative status

Percentage calculated as $(n/N) \times 100$.

Table 98: Subjective Quality of Vision – Visian MICL Study (All Eyes)

Quality of Vision Grading	Preoperative N(524) n (%)	36 Months N(346) n (%)
Very Good/Excellent	288 (55.0%)	267 (77.0%)
Poor/Very Poor	61 (11.6%)	20 (5.8%)

¹ Percentage calculated as $(n/N) \times 100$.

Subjective Symptoms Stratified by Optic Diameter

Subjective symptoms reported by subjects were stratified into 4 groups based on the optic diameter: 4.9 mm, 5.2 mm, 5.5 mm and 5.8 mm. Glare was absent/mild in 82.4% (75/91) of subjects in the 4.9 mm, 90.3% (65/72) in the 5.2 mm, 91.8% (45/49) in the 5.5 mm and 89.9% (125/139) in the 5.8 mm groups. Marked/severe glare occurred in 3.3% (3/91) of eyes with the 4.9 mm, 2.8% (2/72) with the 5.2 mm, 4.1% (2/49) with the 5.5 mm and 1.4% (2/139) with the 5.8 mm optic at 36 months postoperatively.

The smaller the optic diameter, the greater the incidence of halos. Halos were absent/mild in 80.2% (73/91) of subjects in the 4.9 mm, 87.3% (62/71) in the 5.2 mm, 89.8% (44/49) in the 5.5 mm and 87.8% (122/139) in the 5.8 mm. Marked/severe halo was dependent upon the Visian MICL lens optic diameter and was 9.9% (9/91) with the 4.9 mm, 2.8% (2/71) with the 5.2 mm, 4.1% (2/49) with the 5.5 mm and 1.4% (2/139) with the 5.8 mm.

Double vision was absent in all eyes with the 5.8 mm optic diameter. Double vision was reported as absent in 95.6% (87/91) of the subjects with the 4.9 mm, 98.6% (71/72) with the 5.2 mm, and 98.0% (48/49) with the 5.5 mm at 36 months. The incidence of marked/severe night driving difficulties negatively correlated with the optic diameter. Marked/severe night driving difficulties was reported in 16.7% (15/90) of eyes in the 4.9 mm group compared to 0% (0/135) with the 5.8 mm. Night driving difficulties were absent / mild in 71.1% (64/90) of eyes using the 4.9 mm, 83.8% (57/68) with the 5.2 mm, 85.4% (41/48) with the 5.5 mm, and 91.9% (124/135) with the 5.8 mm.

A similar trend between the subjective symptom and the 36-month follow-up shows a negative correlation between the incidence/severity of night vision difficulties and the optic diameter. No cases of marked/severe night vision difficulties occurred with the 5.8 mm. Subjective night vision difficulties 36 months after Visian MICL lens insertion were absent/mild in 73.6% (67/91) of eyes with 4.9 mm, 84.7% (61/72) with the 5.2 mm, 83.7% (41/49) with the 5.5 mm, and 90.6% (126/139) with the 5.8 mm.

POST-APPROVAL CONTINUATION OF THE VISIAN MICL PMA CLINICAL TRIAL

As a condition of approval, STAAR Surgical was required to follow subjects of the original Visian MICL PMA study cohort through 60 months (5 years) after lens implantation with the specific objective of collecting long-term data on endothelial cell loss and on AEs/complications.

Study Objective

The objective of this post-approval study was to collect new long-term data on endothelial cell loss and on AEs/complications in order to assess long-term safety of the Visian MICL lens. Only data on these safety parameters are updated in this section.

Study Design

This post-approval study consisted of the extended follow-up of the original Visian MICL PMA study cohort. It was a single-arm study with follow-up visits scheduled at 48 and 60 months (4 and 5 years) post-implantation.

Total Number of Enrolled Study Sites and Subjects

Of the 526 eyes (294 subjects) enrolled at 14 sites in the United States in the original Visian MICL PMA study, 335 eyes of 192 subjects were seen at 60 months (5 years) or later, as shown in **Table 97**. However, since this post-approval study was initiated a number of years after the first implants of the Visian MICL lens in the original Visian MICL PMA study, some subjects were more than 60 months postoperative at the time of initiation of the post-approval study. These subjects were seen for a final visit and are included in the "≥ 60 Months" columns.

Table 99: Accountability – Post-Approval Continuation of Visian MICL Study

Accountability (all implanted eyes, N=526)	12 Months	24 Months	36 Months	48 Months	60 Months	≥ 60 Months
Available for Analysis	472	459	384	248	225	335
Discontinued (ICL Lens Removals) ¹	0	1	5	5	10	11
Missed Visit/CRF not Received	40	44	84	192	176	NA
Missing	0	0	0	1	4	NA
Lost to Follow-up	14	22	53	80	111	180
% Accountability ²	89.7%	87.4%	73.7%	47.7%	43.9%	65.6%

¹ Cumulative total number of eyes discontinued is 11

² % Accountability is equal to (Available for analysis)/(All Implanted Eyes-Discontinued-Missing)

Adverse Events and Complications

The incidence of AEs, complications and surgical reinterventions reported from time of surgery through the end of the post-approval study period (≥60 months), are shown in **Table 100** through **Table 102**.

Table 100: Adverse Events Through ≥60 months – Post-Approval Continuation of Visian MICL Study

Adverse Event ¹	Cumulative n/526 (%*)	≤ 12 Mo n/526 (%*)	>12-24 Mo n/462 (%*)	>24-36 Mo n/426 (%*)	>36-48 Mo n/276 (%*)	>48-60 Mo n/346 (%*)	≥60 Mo n/348(%*)
Endophthalmitis	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Hyphema	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Hypopyon	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
IOL Dislocation	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Cystoid Macular Edema	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Pupillary Block	17, (3.2%)	17, (3.2%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Retinal Detachment	3, (0.6%)	1, (0.2%)	1, (0.2%)	1, (0.2%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Surgical Reintervention ²	43, (8.2%)	28, (5.3%)	4, (0.9%)	4, (0.9%)	2, (0.7%)	4, (1.2%)	1, (0.3%)
Corneal Edema (after 1 week)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Iritis (after 1 week)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)
Iris Prolapse Repair	1, (0.2%)	1, (0.2%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)	0, (0.0%)

¹ Loss of VA, cataract development, raised IOP requiring pharmacologic intervention, endothelial cell loss and other unclassified complications are not included in the table but are discussed in the sections below.

² Refer to section below for details on Surgical Reinterventions.

* Percentage calculated as (n/N)*100.

Surgical Reinterventions

A total of 43/526 eyes (8.2%) underwent surgical reintervention during the study (**Table 101**). Of these, 23/526 (4.4%) eyes had repositioning (4 eyes), removal (10 eyes) or replacement (8 eyes) of the Visian MICL lens, and 1 eye had Visian MICL lens replacement and then removal. Each case of Visian MICL lens removal during the study was performed in conjunction with cataract surgery. An additional 20/526 eyes (3.8%) underwent repeat YAG iridotomy or additional irrigation/aspiration during the early postoperative time period. Of these, 17 eyes were treated with YAG laser iridotomy for pupillary block, and 3 eyes were treated with repeat irrigation and aspiration for removal of retained viscoelastic.

Table 101: ICL Lens Related Additional Surgery Through ≥60 months – Post-Approval Continuation of Visian MICL Study

ICL Lens Related Additional Surgery	Cumulative n/526 (%*)	≤ 12 Mo n/526 (%*)	>12-24 Mo n/462 (%*)	>24-36 Mo n/426 (%*)	>36-48 Mo n/276 (%*)	>48-60 Mo n/346 (%*)	≥60 Mo n/348(%*)
Visian MICL Lens Repositioning	4 (0.8%)	4 (0.8%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Visian MICL Lens Replacement, then Removal	1 (0.2%)	0 (0.0%)	0 (0.0%)	1 (0.2%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Visian MICL Lens Replacement	8 (1.5%)	4 (0.8%)	2 (0.4%)	2 (0.4%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
Visian MICL Lens Removal	10 (1.9%)	0 (0.0%)	2 (0.4%)	1 (0.2%)	2 (0.7%)	4 (1.2%)	1 (0.3%)
Raised IOP Requiring Surgery ¹	20 (3.8%)	20 (3.8%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)	0 (0.0%)
TOTAL	43 (8.2%)	28 (5.3%)	4 (0.9%)	4 (0.8%)	2 (0.7%)	4 (1.2%)	1 (0.3%)

¹ Refer to section on IOP for details.

* Percentage calculated as (n/N)*100.

Refractive Procedures

A total of 22/526 eyes (4.2%) underwent refractive procedures during the study; this consisted of 17/526 LASIK (3.2%) procedures and 5/526 Arcuate Keratotomy (AK) (1.0%) procedures, as seen in **Table 102**.

Table 102: Refractive Procedures Through ≥60 months – Post-Approval Continuation of Visian MICL Study

Refractive Procedure	≤ 12 Mo n/N (%) ¹	>12-24 Mo n/N (%)	>24-36 Mo n/N (%)	>36-48 Mo n/N (%)	>48-60 Mo n/N (%)	≥60 mo n/N (%)	Total n/N (%)
LASIK	15/526 (2.9%)	1/524 (0.2%)	0/448 (0.0%)	0/256 (0.0%)	0/231 (0.0%)	1/117 (0.9%)	17/526 (3.2%)
AK	3/526 (0.6%)	2/524 (0.4%)	0/448 (0.0%)	0/256 (0.0%)	0/231 (0.0%)	0/117 (0.0%)	5/526 (1.0%)

¹ Percentage calculated as (n/N)*100.

Best Corrected Visual Acuity (CDVA) Loss

Eighteen eyes of 16 subjects reported a significant vision loss of ≥2 lines in CDVA between 12 months and ≥ 60 months. Reasons for significant vision loss included cataract development (9 eyes), myopic degeneration (1 eye), retinal detachment (1 eye) and unknown etiology was reported for 4 eyes. For 3 eyes, decrease in CDVA was transient without intervention.

At the final study visit (which ranged from 18 to 62 months), 11 of these 18 eyes reported an improvement in CDVA of 2 to 10 lines compared to preoperative CDVA, attributed to cataract surgery, refractive surgery or reversal of transient vision loss. In the remaining 7 eyes, vision loss of ≥2 lines was persistent at the final study visit (which ranged from 36 to 60 months).

The number of eyes reporting a decrease in either 2 lines or > 2 lines is reported in **Table 103**.

Table 103: CDVA Loss Through ≥60 months – Post-Approval Continuation of Visian MICL Study

Decrease in CDVA	12 Mo n/N (%)	24 Mo n/N (%)	36 Mo n/N (%)	48 Mo n/N (%)	60 Mo n/N (%)	≥ 60 Mo n/N (%)
Decrease >2 Lines	1/469 (0.2%)	2/456 (0.4%)	3/384 (0.8%)	1/242 (0.4%)	1/222 (0.4%)	2/331 (0.6%)
Decrease =2 Lines	2/469 (0.4%)	3/456 (0.6%)	1/384 (0.3%)	1/242 (0.4%)	2/222 (0.8%)	2/331 (0.6%)

¹ Percentage calculated as (n/N)*100.

Lens Opacity and Visually Significant Cataract Formation

Table 104 provides the type of cataracts of grade trace or greater that developed over time for the PMA Study cohort. The long-term incidence of anterior subcapsular opacity secondary to implantation of the Visian MICL lens has been studied in 526 eyes of 294 subjects followed for up to 7.5 years, with 334 eyes available for analysis at 5 or more years. A total of 31 eyes developed an anterior subcapsular opacity.

Table 104: Cataract Through ≥60 months – Post-Approval Continuation of Visian MICL Study

Cataract Type	Preop n/N (%)	<12 Mo n/N (%)	12 Mo n/N (%)	24 Mo n/N (%)	36 Mo n/N (%)	48 Mo n/N (%)	60 Mo n/N (%)	≥60 Mo n/N (%)	Cumulative Number of Eyes
Nuclear	4/526 (0.8%)	4/526 (0.8%)	2/472 (0.4%)	1/457 (0.2%)	3/381 (0.8%)	0/245 (0.0%)	0/225 (0.0%)	3/334 (0.9%)	13
Cortical	2/526 (0.4%)	2/526 (0.4%)	0/472 (0.0%)	1/457 (0.2%)	4/380 (1.1%)	1/245 (0.4%)	0/225 (0.0%)	0/334 (0.0%)	8
Posterior Subcapsular	0/526 (0.0%)	0/526 (0.0%)	0/472 (0.0%)	0/457 (0.0%)	2/381 (0.5%)	0/245 (0.0%)	0/225 (0.0%)	2/334 (0.6%)	4
Anterior Subcapsular	0/526 (0.0%)	8/526 (1.5%)	3/472 (0.6%)	4/457 (0.9%)	2/381 (0.5%)	8/245 (3.3%)	2/225 (0.9%)	4/335 (1.2%)	31
Total Number of Eyes*	6	13	3	5	8	9	2	5	45

¹ Percentage calculated as (n/N)*100.

* Final row may not sum to number of nuclear, cortical or subcapsular cataracts, as some eyes had multiple types of cataracts.

Visually significant cataracts of all types, involving a vision loss of ≥2 lines in CDVA, were reported in 9/526 eyes (1.7%) through the extended follow-up study period: 1 anterior subcapsular cataract (ASC) at 18 months, 3 ASC at 48 months and 1 surgically induced ASC which was reported to have a 2 line loss of CDVA at 24 months after Visian MICL lens implantation; 1 nuclear cataract (NC) at 12 months, 1 at 30 months, 2 at 36 months.

Per eye, the risks of developing any anterior subcapsular opacity, developing a visually significant anterior subcapsular opacity, or of having cataract surgery for any type of cataract were calculated using Kaplan-Meier analyses. As provided in **Table 105**, these risks were 6.1%, 1.2% and 3.1% at 60 months and 12.4%, 1.2% and 3.1% at 84 months, respectively.

Table 105: Lens Opacification Risk Analysis – Post-Approval Continuation of Visian MICL Study

	≤12 Mo	>12 -24 Mo	>24 -36 Mo	>36-48 Mo	>48-60 Mo	≥60 Mo
Any Anterior Subcapsular Opacity (ASC)						
Number at risk at period start	526	499	477	441	366	251
Events during period	9	4	3	4	7	4
Survival estimate at period end	98.3%	97.5%	96.9%	95.9%	93.9%	87.6%
1-survial estimate (risk)	1.7%	2.5%	3.2%	4.1%	6.1%	12.4%
Visually Significant ASC						
Number at risk at period start	526	507	487	450	379	261
Events during period	1	1	0	0	3	0
Survival estimate at period end	99.8%	99.6%	99.6%	99.6%	98.8%	98.8%
1-survial estimate (risk)	0.2%	0.4%	0.4%	0.4%	1.2%	1.2%
Cataract Surgery for Any Type of Cataract						
Number at risk at period start	526	505	484	448	376	258
Events during period	3	3	2	3	3	0
Survival estimate at period end	99.4%	98.8%	98.4%	97.7%	96.9%	96.9%
1-survial estimate (risk)	0.6%	1.2%	1.6%	2.3%	3.1%	3.1%

Intraocular Pressure (IOP)

a) Changes in IOP from Baseline

Postoperatively, IOP >25 mmHg or an increase of >10 mmHg over preoperative was reported in 62/526 eyes (11.8%) of the Visian MICL Lens PMA cohort through ≥ 60 months.

Table 106: Changes in IOP from Baseline Through ≥60 months – Post-Approval Continuation of Visian MICL Study

IOP (mmHg)	Preop n/N (%)	1 D n/N (%)	14 D n/N (%)	1 Mo n/N (%)	3 Mo n/N (%)	6 Mo n/N (%)	12 Mo n/N (%)	24 Mo n/N (%)	36 Mo n/N (%)	48 Mo n/N (%)	60 Mo n/N (%)	≥60 Mo n/N (%)
> 10 mmHg over Baseline	NA	23/526 (4.4%)	20/526 (3.8%)	12/524 (2.3%)	2/522 (0.4%)	0/511 (0%)	2/501 (0.4%)	2/469 (0.4%)	1/410 (0.2%)	2/348 (0.6%)	6/262 (2.3%)	4/263 (0.4%)
> 25 mmHg	0/526 (0%)	23/526 (4.4%)	16/526 (3.0%)	11/524 (2.1%)	2/522 (0.4%)	0/511 (0%)	2/501 (0.4%)	2/469 (0.4%)	1/410 (0.2%)	3/348 (0.9%)	7/262 (2.7%)	4/263 (0.4%)

¹ Percentage calculated as (n/N)*100.

b) Raised IOP Requiring Surgery

A total of 20/526 eyes (3.8%) experienced raised IOP requiring intervention. An additional YAG iridotomy was performed on 17 of the eyes for pupillary block and 3 eyes had repeat irrigation and aspiration at 1 day postoperative to remove retained viscoelastic. All of these events occurred in the early postoperative period, most frequently at 1 to 2 days postoperative.

c) Raised IOP Requiring Pharmacologic Intervention

A total of 7 eyes of 4 subjects in the PMA cohort developed glaucoma during the clinical trial. Open angle glaucoma was diagnosed for 4 eyes (2 subjects) and the remaining 3 eyes of 2 subjects the type of glaucoma was not specified. None of these eyes required secondary surgical intervention for treatment of IOP during the study.

Upon gonioscopic examination, no anterior synechiae, transillumination defects, or abnormal angle depth was observed in any of these 7 eyes. However, abnormal pigmentation was observed in 6 eyes, with 2 eyes of a single subject diagnosed with open angle glaucoma and possible secondary pigment dispersion at 6 years postoperatively.

Table 107: Glaucoma – Post-Approval Continuation of Visian MICL Study

No. of Eyes	Type of Glaucoma	Onset	Abnormal Pigmentation	Anterior Synechiae	Transillumination Defects	Angle Depth
1 (1 subject)	Unspecified	62 Mo	None	None	None	Normal
2 (1 subject)	Unspecified	5 Mo, 12 Mo	Yes	None	None	Normal
2 (1 subject)	Open Angle	37 Mo, 53 Mo	Yes	None	None	Normal
2 (1 subject)	Open Angle, possibly 2° pigment dispersion	71 Mo, 73 Mo	Yes	None	None	Normal

Gonioscopic Findings

In the post-approval study, investigators were asked to perform gonioscopy at the 48 Month (Form 9) and/or > 60 Month (Form 10) study visits. Specifically, investigators were to report on the absence or presence of peripheral anterior synechiae, the absence or presence of abnormal pigment suggestive of pigment dispersion and normal or abnormal angle depth.

Table 108: Gonioscopic Findings – Post-Approval Continuation of Visian MICL Study

Finding	48 Month (Form 9)			≥ 60 Months (Form 10)		
	Absent % ¹ (n)	Present % ¹ (n)	Total*	Absent % ¹ (n)	Present % ¹ (n)	Total*
Peripheral Anterior Synechiae	99.05% (104)	0.95% (1)	105	99.66% (293)	0.34% (1)	294
Abnormal Pigment Suggestive of Pigmentary Dispersion	95.28% (101)	4.72% (5)	106	94.00% (282)	6.00% (18)	300
	Normal	Abnormal	Total	Normal	Abnormal	Total
Angle Depth	99.05% (104)	0.95% (1)	105	100% (298)	0.00% (0)	298

* Total number of eyes with gonioscopy was performed at that visit. (N)

¹ Percentage calculated as (n/N)*100.

Other Findings

At the 48 month visit, no “other findings” were reported. At the > 60 month visit, there were a total of 24 comments reported under “other findings”. They were: “Heavy Pigment” (n = 8); “Moderate Pigment” (n = 8); “Light Pigment” (n = 4); “Transillumination defects” (n = 2) and “Myopic Degeneration and Pigment Changes in Macula” (n = 2).

Slit Lamp Findings

Table 109 summarizes the incidence of pigment on cornea, pigment on Visian MICL lens and transillumination defects that occurred at different time points reported throughout the study follow-up period:

Table 109: Slit Lamp Findings – Post-Approval Continuation of Visian MICL Study

Finding/Onset	<12 Mo % ¹ (n/N)	12 Mo % ¹ (n/N)	24 Mo % ¹ (n/N)	36 Mo % ¹ (n/N)	48 Mo % ¹ (n/N)	≥ 60 Mo % ¹ (n/N)
Pigment on cornea	0.0% (0/526)	0.0% (0/472)	0.0% (0/459)	0.0% (0/384)	0.0% (0/248)	1.2% (4/335)
Pigment on Visian MICL lens	2.5% (13/526)	0.4% (2/472)	2.0% (9/459)	1.8% (7/384)	2.0% (5/248)	5.1% (17/335)
Transillumination defects	0.6% (3/526)	0.0% (0/472)	0.2% (1/459)	0.3% (1/384)	0.4% (1/248)	0.9% (3/335)

¹ Percentage calculated as (n/N)*100.

Endothelial Cell Density (ECD)

Specular microscopy was performed on a subgroup of the original PMA study cohort with data available through ≥60 months postoperatively. A central reading center was used to minimize the inherent variability associated with endothelial cell counts.

Table 110 provides detail on the number of readable specular microscopy images captured at each time point in the study.

Table 110: Specular Microscopy – Post-Approval Continuation of Visian MICL Study

	Preop	3 Mo	12 Mo	24 Mo	36 Mo	48 Mo	60 Mo	72 Mo	84 Mo
Total Cohort (N)	526		472	459	384	248	225	86	44
Eyes with readable ECD n (%)	192 (36.5%)	209	246 (52.1%)	220 (47.9%)	174 (45.3%)	146 (58.9%)	113 (50.2%)	37 (43.0%)	27 (61.4%)
Eyes with both Preop and Postop readable ECD n (%)	NA	162	175 (37.1%)	151 (32.9%)	132 (34.4%)	109 (44.0%)	85 (37.8%)	15 (17.4%)	19 (43.2%)

¹ Percentage calculated as (n/N)*100.

The analysis of ECD over time was conducted on eyes with both pre and postoperative ECD counts. Mean ECD results from clinical trial subjects are shown in **Table 109**.

Table 111: ECD Analysis Through ≥60 months – Post-Approval Continuation of Visian MICL Study

Visit	Mean	SD	90% Confidence Limits
Preop	2657	290	2622 to 2692
3 Mo	2570	340	2532 to 2609
12 Mo	2548	349	2511 to 2584
24 Mo	2479	357	2439 to 2518
36 Mo	2454	348	2411 to 2498
48 Mo	2396	367	2346 to 2447
≥60 Mo	2298	354	2252 to 2345

During the PMA trial and subsequent long-term follow-up of the PMA cohort, 13 eyes of 10 subjects (11.5% 13/113 of those available for evaluation ≥ 60 months after surgery) reported significant endothelial cell loss (> 30% loss of central ECD). Of these 13 eyes, 3 eyes of 3 subjects experienced this level of endothelial cell loss (30.8 – 45.6%) between baseline and the first 12 months of follow-up, and it was presumed to be the result of surgery; the remaining 10 eyes of 7 subjects had this level of endothelial cell loss (30.9 – 42.6%) at the final study visit (≥ 60 months, between 5.0 and 6.7 years).

Table 112: ECD loss from Preoperative Values – Post-Approval Continuation of Visian MICL Study

ECD loss from Preop (%)	12 Mo N=175 n (%)	36 Mo N=132 n (%)	≥ 60 Mo N=115 n (%)
≥10%	22 (12.6%)	44 (33.3%)	77 (67.0%)
≥15%	8 (4.6%)	22 (16.7%)	50 (43.5%)
≥20%	4 (2.3%)	12 (9.1%)	30 (26.1%)
≥30%	3 (1.7%)	2 (1.5%)	13 (11.3%)

¹ Percentage calculated as (n/N)*100.

The available data from the clinical study demonstrate a mean percentage change from baseline to 60 months of 12.3% (SD 9.4%), based on subjects with data at both baseline and ≥ 60 months.

Table 113 provides the mean, standard deviation, median, interquartile range, and range of percent change in ECD. These data represent changes in ECD between:

- The preoperative visit and the 12 month visit (for all eyes with ECD data at both visits);
- The 1 year visit to the 3 year visit (for all eyes with ECD data at both visits); and
- The 3 year visit to the final visit at 5 years or later (for all eyes with ECD data at both visits)

Table 113: Change in ECD over Time – Post-Approval Continuation of Visian MICL Study

Endothelial Cell Density	Preoperative Visit to 12 month visit	For all eyes with ECD data at both visits:	
		1 year visit to 3 year visit	3 year visit to Final Visit at 5 years or later
N (ECD observations with data at both visits)	175	150	108
Mean (SD) % Change in ECD	-3.19 (7.59)	-5.04 (8.09)	-6.74 (5.15)
Median % Change in ECD	-2.45	-4.27	-6.24
Interquartile Range % Change in ECD (Q1 to Q3)	0.97 to -2.45	-1.41 to -4.27	-3.04 to -9.93
Range (Min, Max) % Change in ECD	16.22,-42.94	11.62, -23.15	4.27, -22.52

The following table provides the predicted percent endothelial cell loss, by year, for a hypothetical patient with preoperative ECD equal to the mean level in the clinical study. For this hypothetical patient, there is 90% confidence that the endothelial cell loss will be between the lower and upper prediction interval bounds at each point in time. The entries in this table are calculated assuming a bi-exponential loss in ECD, i.e., a rapid initial phase of cell loss in the early postoperative period related to surgical trauma, followed by a slow, chronic phase of cell loss thereafter. Rates of predicted long term loss are derived from clinical data collected through 5 to 7 years postoperatively. The calculated chronic rate of loss from this post-approval data is approximately 1.8% per year.

Table 114: Predicted Percent Endothelial Cell Loss

Time from procedure	Predicted Percent Cell Loss	90% prediction interval*	
		Lower	Upper
3 months	1%	-20%	23%
1 year	4%	-18%	25%
2 years	5%	-16%	27%
3 years	8%	-14%	29%
4 years	9%	-12%	31%
5 years	11%	-10%	33%
10 years	20%	-2%	42%
15 years	28%	6%	50%
20 years	35%	13%	57%
25 years	42%	19%	64%
30 years	47%	25%	70%
35 years	53%	30%	75%
40 years	57%	35%	80%
45 years	62%	39%	84%
50 years	66%	43%	88%
55 years	69%	46%	92%

* Note: Positive values represent levels of % ECD loss; negative values represent levels of % ECD gain.

Other Complications

No cases of endophthalmitis, hyphema, hypopyon, cystoid macular edema or corneal ulcer were reported during the study. Corneal haze, corneal edema or iritis were not reported after the 1 week visit. One case each of iris prolapse (1/526, 0.2%), macular hemorrhage (1/526, 0.2%) and subretinal hemorrhage (1/526, 0.2%) were reported at 1 day, 1 week and 3 months postoperative, respectively. Retinal detachment was reported in 3 eyes (3/526, 0.6%) at 4, 22 and 31 months after Visian MICL lens implantation.

A case of anisocoria (unequal pupil size) has been reported for a subject implanted with an ICL lens in another clinical trial.

Study Strengths and Limitations

This post-approval study uses the original Visian MICL IDE study cohort, following patients who had already completed 36 months of follow-up; therefore, long-term data (60-months or later) is available sooner as opposed to a new-enrollment study. Additionally, this is the only post-approval sub-study that collected ECD data. However, the 60 month follow-up rate of 65.3% (335/515) is less than optimal. Biases could have been introduced into the study results because of the loss to follow-up, which could limit the generalizability of the study results.

POST-APPROVAL STUDY OF THE EFFECT OF THE VISIAN MICL LENS ON AXIAL LENGTH MEASUREMENT

The Visian MICL lens was evaluated in a prospective, non-randomized study of 30 eyes of 30 subjects to assess the effect of the lens on the measurement of the eye's axial length, and to determine whether the Visian MICL lens affects this measurement. Study inclusion criteria were:

- Moderate to high myopia (-3 D to -20 D measured as spherical equivalent of the manifest refraction) scheduled to undergo implantation of the commercially available Visian MICL lens.
- Subject meets all of the Indications for Use criteria for the commercially available Visian MICL lens.
- Ability to be measured with the IOL Master Axial Length measurement device.
- Willingness to comply with the sub-study preoperative and postoperative visit requirements.

There were no study exclusion criteria.

The subjects underwent implantation of the commercially available Visian MICL lens. The axial length was measured preoperatively and between one week and one month postoperatively. All axial length measurements were obtained using a Carl Zeiss IOL Master, a non-contact partial coherence laser interferometer. The difference in the pre and postoperative axial length was calculated individually for each eye.

Of the 30 subjects, 11 were male, 19 female, 29 Caucasian and 1 Asian. The Visian MICL lens power of the lens implanted averaged -10.68 D (range -3.50 D to -16.00 D). The preoperative axial length averaged 27.28 mm (range 23.69 mm to 34.32 mm) and the postoperative axial length averaged 27.28 mm (range 23.72 mm to 34.51 mm). The average difference in preoperative and postoperative axial lengths is -0.02 mm (range -0.23 mm to + 0.19 mm).

The correlation coefficient was calculated based on a regression analysis on the pre and postoperative data. The results of the analysis show that the variance preoperative is statistically equivalent to the variance postoperative at 95% confidence. The average difference of -0.02 mm in axial length measurement pre and postoperative would change IOL power prediction by 0.05 D, which is well below the measurement of error of IOL power manufacturers.

The data in this study suggests that the Visian MICL lens has a negligible influence on axial length measurements for IOL power calculations, when measurements are based on partial coherence laser interferometry. The accuracy of ultrasound-based measurement of axial length is unknown.

Study strengths include its representative sample (no exclusion criteria) and relevance to clinical questions surrounding axial measurement. Study limitations include its applicability only to laser interferometry-based measurement and not to ultrasound measurement of axial length, and the use of only 2 investigational sites.

POST APPROVAL ADVERSE EVENT STUDY - VISIAN MICL LENS FOR MYOPIA

A survey study was conducted in the US after the Visian MICL lens was approved by the FDA. The goal of this study was to collect safety information from patients who had Visian MICL surgery in the general population. All patients who consented to participate were asked to complete surveys at scheduled times up to 5 years after their Visian MICL surgery. The surveys asked patients to report any complications or additional eye surgeries because of the Visian MICL lens.

Description of the Study Patient Group:

- 2999 eyes of 1547 patients implanted with the Visian MICL lens participated;
- Most patients were white (Caucasian) and over half of the patients were female;
- Patients ranged from 17 to 77 years of age at time of surgery.

The surveys asked for information about the following adverse events:

- Problems with endothelial cells;
- Cataract formation;
- Medical treatment for inflammation inside the eye;
- Medical treatment for intraocular pressure and damage to the optic nerve caused by glaucoma;
- Surgery because of retinal detachment;
- Surgery to remove, replace or reposition the Visian MICL lens;
- Other complications in the eye.

The cumulative incidence per eye for each of the events assessed in the survey in addition to the cumulative incidence of the same events from the Visian MICL Lens PMA clinical study for comparison are presented in **Table 115**.

Table 115: Visian MICL Cumulative AEs – Post-Approval AE Study, Comparison to PMA Clinical Study

Survey Questionnaire	60 months - Cumulative %, (n/N)	PMA Study > 60 months - Cumulative %, (n/N)
1-Corneal problems	0.3%, (5/2999)	0%, (0/526) Corneal Edema (after 1 week)
2-Cataract development	5.1%, (154/2999)	8.6%, (45/526)
3-Treated intraocular inflammation	0.5%, (14/2999)	0.0%, (0/526)
4-Treated IOP or glaucoma	1.6%, (47/2999)	1.3%, (7/526)
5-Retinal Detachment Surgery	0.4%, (13/2999)	0.6%, (3/526)
6-Remove, replace or reposition Visian MICL lens	4.2%, (126/2999)	4.3%, (23/526)

¹ Percentage calculated as (n/N)*100.

The Visian MICL Lens PMA clinical study only enrolled subjects ≤ 45 years of age. A comparison of the cumulative incidence of the events between the PMA Clinical Study and the survey questions for patients ≤ 45 years of age at the time of Visian MICL surgery are provided in **Table 116**.

Table 116: Cumulative AEs, Comparison to Visian MICL PMA Clinical Study (Ages ≤ 45 yrs old at time of Surgery)

Survey Questionnaire	60 months - Cumulative %, (n/N)	PMA Study ≥ 60 months - Cumulative %, (n/N)
1-Corneal problems	0.0%, (0/2527)	0%, (0/526) Corneal Edema (after 1 week)
2-Cataract development	3.0%, (75/2527)	8.6%, (45/526)
3-Treated intraocular inflammation	0.5%, (13/2527)	0.0%, (0/526)
4-Treated IOP or glaucoma	1.5%, (38/2527)	1.3%, (7/526)
5-Retinal Detachment Surgery	0.3%, (7/2527)	0.6%, (3/526)
6-Remove, replace or reposition Visian MICL lens	2.9%, (74/2527)	4.3%, (23/526)

¹ Percentage calculated as $(n/N) \times 100$.

Glare was reported for 2.8% (85/2999) and halos were reported for 5.2% (156/2999) of all implanted eyes in the survey study. The cumulative per eye incidence of glare and halo at 36 months after surgery from the survey questionnaire was compared to the PMA clinical study data on worsening of glare and halo at 36 months compared to baseline. The comparison between the studies is made for patients ≤ 45 years of age at the time of Visian MICL surgery and is provided in **Table 117**.

Table 117: Cumulative reports of Glare and Halos at 36 Months, Comparison to PMA Study, (Ages ≤ 45 yrs old at time of Surgery)

Survey Questionnaire	%, (n/N)	PMA Study %, (n/N)
Glare	2.6%, (66/2527)	9.6%, (34/351)
Halos	5.6%, (142/2527)	11.5%, (40/350)

¹ Percentage calculated as $(n/N) \times 100$.

Overall, patient responses to surveys provided similar information to what was found in the FDA safety and effectiveness clinical study of 526 eyes of 294 patients. This study included patients over 45 years of age. This age group was not included in the FDA safety and effectiveness study of the Visian MICL lens. These older patients reported a higher rate of cataracts and need for a second surgery than patients who were 45 or younger at the time of initial Visian MICL surgery.

INSTRUCTIONS FOR USE

CAUTION: Implantation of an EVO/EVO+ ICL/TICL lens should only be attempted by a surgeon who is highly skilled in the required surgical technique and has completed the EVO/EVO+ ICL/TICL Certification Program.

CAUTION: Do not use EVO/EVO+ ICL/TICL lens if package has been opened or damaged. The sterility of the lens may be compromised.

ICL Lens Handling Precautions

- Choice of the proper EVO/EVO+ ICL/TICL lens size should be carefully considered prior to surgery.
- Check the label of the EVO/EVO+ ICL/TICL lens package for proper lens model and power.
- Open the package to verify the dioptric power of the lens.
- Handle the EVO/EVO+ ICL/TICL lens by the haptic portion. Do not grasp the optic with forceps as this could potentially lead to damage to the smooth anterior and posterior optical surfaces.
- Never touch the center of the optic with instruments once the EVO/EVO+ ICL/TICL lens is placed inside the eye. Inadvertent pressure through the optic could potentially damage the central crystalline lens resulting in a lens opacity.
- STAAR Surgical recommends using only the Accuject Refra-AR2900, LIOLI-24, or MicroSTAAR™ Injector systems (Models MSI-TF and MSI-PF with SFC-45 Cartridge), to insert the EVO/EVO+ ICL/TICL lens in the folded state.
- The EVO/EVO+ ICL/TICL lens should be carefully examined in the operating room prior to implantation.
- The EVO/EVO+ ICL/TICL lens should not be exposed to any solutions other than the normally used intraocular irrigating solutions (e.g., isotonic saline, BSS, viscoelastic, etc.)
- Keep the EVO/EVO+ ICL/TICL lens moist. It is recommended that the EVO/EVO+ ICL/TICL lens be held in sterile BSS solution prior to implantation.
- The EVO/EVO+ ICL/TICL lens should be handled carefully. No attempt should be made to reshape or cut any portion of the lens. Do not apply undue pressure to the EVO/EVO+ ICL/TICL lens optical portion with a sharp object since this could perforate the optic.
- The intended location of the EVO/EVO+ ICL/TICL lens is behind the iris within the posterior chamber and in front of the anterior capsule of the crystalline lens.
- The EVO/EVO+ ICL/TICL lens is manufactured so that rotation of no more than 22.5 degrees (2/3 clock hours) is necessary.
- It is recommended that the surgeon not rotate the EVO/EVO+ ICL/TICL lens more than 22.5 degrees from horizontal.
- Complete irrigation and aspiration of viscoelastic from the eye after completion of the surgical procedure is essential. Viscoelastic products that may be difficult to aspirate should not be used.

NOTE: The long term effects of phakic IOL implantation have not been determined. Physicians should continue to monitor implant patients postoperatively on a regular basis.

SURGICAL PRECAUTIONS/INFORMATION

Preoperative Information

Preoperative ECD Measurements

An ECD measurement should be performed preoperatively to determine if candidates meet the minimum ECD requirements based upon age and true ACD. The true ACD is defined as the distance from the apex of the **posterior** corneal surface to the apex of the **anterior** crystalline lens surface. Many measuring devices provide an ACD measurement defined as the distance from the apex of the anterior corneal surface to the apex of the anterior crystalline lens surface. If the surgeon is using an instrument that measures from the anterior corneal surface, the thickness of the cornea must be subtracted to get the true ACD.

EVO/EVO+ ICL/TICL Lens Length Determination

During the original US PMA clinical study, sizing of the ICL myopic lenses (12.1 mm to 13.7 mm) was determined by the horizontal white-to-white and the ACD measurements (true ACD, defined as the distance from the apex of the posterior corneal surface to the apex of the anterior crystalline lens surface). For eyes with ACD measurements of ≤ 3.5 mm, the lens size was calculated by adding 1.1 mm to the horizontal white-to-white measurement. Eyes exhibiting an ACD greater than 3.5 mm required the addition of up to 1.6 mm to the white-to-white measurement, up to a maximum length of 13.7 mm. Calculated lens sizes between the available lens diameters (in 0.5 mm steps) were generally rounded down if the ACD was ≤ 3.5 mm and rounded up if the ACD was >3.5 mm.

Analyses of all of the collected clinical data resulted in slightly modified recommendations for sizing of the EVO/EVO+ ICL/TICL lens as compared to those used in the clinical trial. A table of recommended EVO/EVO+ ICL/TICL lens lengths based upon white-to-white and ACD measurements is given below.

Table 118: Recommended EVO/EVO+ ICL/TICL Lens Overall Diameter by White to White and ACD Measurements

White-to-White (mm)	True ACD (mm)		
	All	≤ 3.5	>3.5
<10.5	Not Recommended	—	—
10.5–10.6	—	Not Recommended	12.1
10.7–11.0	12.1	—	—
11.1	—	12.1	12.6
11.2–11.4	12.6	—	—
11.5–11.6	—	12.6	13.2
11.7–12.1	13.2	—	—
12.2	—	13.2	13.7
12.3–12.9	13.7	—	—
≥ 13	Not Recommended	—	—

White-to-White Measurements

The white-to-white measurement is an indirect measurement and does not correlate with sulcus-to-sulcus measurements. Newer advancements in the direct measurement of the ciliary sulcus such as ultrasonic biomicroscopy (UBM) should be considered as alternative methods for the determination of the desired EVO/EVO+ ICL/TICL lens overall diameter. At present there is no large series study demonstrating the effectiveness of UBM in EVO/EVO+ ICL/TICL lens sizing.

Learning Curve/Individual Surgeon Variability Issues

A learning curve and individual surgeon variability was seen in the clinical trial in terms of early anterior subcapsular lens opacities, removals and reinsertions of the lens at the time of surgery, and lens replacements due to sizing.

Refraction

A cycloplegic refraction is recommended to confirm the accuracy of the manifest refraction.

EVO/EVO+ ICL/TICL Lens Power Calculation

Implantation of the EVO/EVO+ ICL/TICL lens requires that a preoperative determination of the dioptric power of the implanted lens be calculated. Achievement of emmetropia is not necessarily a desirable postoperative goal and factors such as visual status of the fellow eye and patient lifestyle should be considered when determining the lens power to be used.


In order to achieve refractive results similar to those found in the PMA study, EVO/EVO+ ICL/TICL lens power and size calculation should be performed using the STAAR ICL Calculation Software.

The ICL calculator will recommend a range of spherical powers along with their expected postoperative values (i.e. residual sphere); or a cylinder power and a range of spherical powers along with their expected postoperative values (i.e. residual sphere, cylinder, axis and spherical equivalent). Selection of lens power is based on the treatment plan of the surgeon for a given eye.

In all cases it is recommended the EVO/EVO+ ICL/TICL lens be implanted horizontally in the eye through a temporal incision.

EVO/EVO+ TICL Lens Implantation Orientation

As part of the implantation procedure, the EVO/EVO+ TICL lens may need to be rotated up to 22.5 degrees clockwise or counterclockwise from the 0°-180° meridian in order to align the lens axis at the preoperative plus cylinder axis. The surgeon should mark the horizontal axis (0°-180°) of the eye at the slit lamp prior to surgery. These horizontal axis marks will be used as reference points to mark the desired orientation of the lens under the operating microscope, using a suitable corneal axis marking device. For example, if the preoperative plus cylinder axis is at 136° and the lens selected has the cylinder axis at 115°, the lens will need to be rotated 21° counterclockwise from the temporal meridian in eye. In this case the desired axis marked on the cornea would be 21° counterclockwise from the 0°-180° meridian. The online ordering software for the EVO/EVO+ TICL lens is designed to generate an Implantation Orientation Diagram (IOD) to guide the surgeon in determining the amount and direction of rotation for the specific lens selected. See example below:



Intraoperative Information

Preparation of the lens for use

CAUTION: Perform the following steps in a sterile field.

- Inspect the lens vial. Ensure that it is not damaged.
- While keeping the vial in a vertical position, remove the aluminum seal and remove the cap.
- Carefully remove the lens from the vial.
- Examine the lens carefully under the microscope for damage or particulate matter.

CAUTION: Do not allow the EVO/EVO+ ICL/TICL lens to dry after removal from the glass vial.

Delivery System

STAAR Surgical recommends using only the Accuject Refra-AR2900, LIOLI-24, or MicroSTAAR Injectors, Model MSI-TF or MSI-PF with SFC-45 Cartridge. For detailed loading instructions, see information provided with the MSI injection system or with the Accuject Refra-AR2900, or with the lens for the LIOLI-24 delivery system.

CAUTION: The EVO/EVO+ ICL/TICL lens should be injected within 1-2 minutes after loading. Viscoelastic materials tend to lose their lubricity if exposed to air too long.

Viscoelastic Usage

Complete removal of viscoelastic from the eye after completion of the surgical procedure is essential. Irrigation for a minimum of one minute with at least 10 – 20 cc of solution is recommended. STAAR Surgical recommends a low molecular weight 2% hydroxypropyl methylcellulose (HPMC) or dispersive, low viscosity ophthalmic viscosurgical device. Do not use short chain sodium hyaluronate acids (viscoelastics) due to increased risk of cataract formation related to trapped viscoelastic.

Inadequate flushing of the viscoelastic from the eye may lead to IOP spikes. IOP should be checked 1 – 6 hours postoperatively so that elevated IOP may be treated in a timely manner.

Postoperative Information

Postoperative EVO/EVO+ ICL/TICL Lens Vault

Lens vault (the distance between the anterior surface of the crystalline lens and the posterior surface of the EVO/EVO+ ICL/TICL lens) should be assessed 24 hours postoperatively at a slit lamp. Although the postoperative vault of the EVO/EVO+ ICL/TICL lens is intended to be approximately equal to the central corneal thickness, we believe that the optimal vault should be between 50% and 150% of central corneal thickness, this being equivalent to a range of 250 to 900 microns. However, in the absence of symptoms, lens vault outside of this range may not necessarily require exchange or removal.

EVO/EVO+ ICL/TICL Lens Removal

It is recommended that the EVO/EVO+ ICL/TICL lens be removed in cases where the vault is insufficient and the patient exhibits early anterior subcapsular cataract. Removal of the EVO/EVO+ ICL/TICL lens may be necessary in cases where the vault is excessive causing narrowing of the anterior chamber angle, thus decreasing aqueous flow. EVO/EVO+ ICL/TICL lens removal may also be necessary for other reasons on an individual basis. The risks involved in EVO/EVO+ ICL/TICL lens replacement have not been studied and are unknown.

Axial Length Measurement Correction for Intraocular Lens (IOL) Power Calculation

The accuracy of ultra-sound based measurement of axial length in an eye with an EVO/EVO+ ICL/TICL lens is unknown. Axial length measurements based upon partial coherence laser interferometry appear to not be significantly affected by implantation of the lens. See section on “Post-Approval Study of the Effect of the Visian MICL on Axial Length Measurement.”

NOTE: More detailed information regarding the recommended Surgical Technique is provided in conjunction with STAAR’s EVO/EVO+ ICL/TICL Physician Certification Program.

MRI Safety Information

The EVO/EVO+ ICL/TICL lens is MR Safe.

SURGICAL PROCEDURE

All physicians must complete the STAAR Surgical EVO/EVO+ ICL/TICL Physician Certification Program prior to using the EVO/EVO+ ICL/TICL lens in a clinical setting.

PATIENT IMPLANT CARD

Each patient who receives an EVO/EVO+ ICL/TICL lens must be provided with an Implant Identification Card. An Implant Identification Card is supplied in the unit package. This card should be given to the patient with instructions to keep it as a permanent record of the implant and to show the card to any eye care practitioner seen in the future.

REPORTING

Adverse Reactions and/or potentially sight-threatening complications that may reasonably be regarded as lens related and that were not previously expected in nature, severity or degree of incidence should be diligently reported to STAAR Surgical immediately at:

USA Phone: (800) 352-7842

Fax: (800) 952-4923

This information is being requested from all surgeons in order to document potential long-term effects of EVO/EVO+ ICL/TICL lens implantation, especially in younger patients. Physicians must report these events in order to aid in identifying emerging or potential problems with the EVO/EVO+ ICL/TICL lens.

HOW SUPPLIED

Each EVO/EVO+ ICL/TICL lens is provided sterile and non-pyrogenic in sealed vials within a sterile thermoform tray placed in a box with labels and product information. The tray and vial containing the EVO/EVO+ ICL/TICL lens are sterilized with steam and should be opened only under sterile conditions.

EXPIRATION DATE

The expiration date on the device package and unit box is the sterility expiration date. If the tray seal and vial seal are not punctured or damaged, sterility is assured until the expiration date indicated on the package label. This device should not be used past the indicated sterility expiration date.

RETURN POLICY FOR STAAR EVO/EVO+ ICL/TICL LENSES

Contact STAAR Surgical. The EVO/EVO+ ICL/TICL lens should be returned dry. Do not attempt to rehydrate.

LENS SPECIFIC RECOMMENDATION

The physician must use the STAAR recommended Injector and Cartridge delivery system for implanting the EVO/EVO+ ICL/TICL lens in the folded state.

WARRANTY AND LIMITATION OF LIABILITY

STAAR Surgical Company warrants that reasonable care was taken in making this product. STAAR Surgical Company shall not be responsible for any incidental or consequential loss, damage, or expense which arises directly or indirectly from the use of this product. Any liability shall be limited to the replacement of any STAAR EVO/EVO+ ICL/TICL lens which is returned to and found to be defective by STAAR Surgical Company.

































This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including but not limited to, any implied merchantability or fitness for use.

STORAGE

Store the EVO/EVO+ ICL/TICL lens at room/ambient temperature.

WARNING: Do not autoclave the EVO/EVO+ ICL/TICL lens. Do not expose to temperature greater than 40 °C. Do not freeze. If temperature requirements are not met, return the EVO/EVO+ ICL/TICL lens to STAAR Surgical.

SYMBOL GLOSSARY

	Medical device		Authorized representative in the European Community
	Do not re-use		CE conformity marking per European Council Directive 93/42/EEC or European Council Regulation (EU) 2017/745
	Do not re-sterilize		Manufacturer
	Do not use if the product sterile barrier system or its packaging is compromised		Date of manufacture
	Body diameter (Optic diameter)		Country of manufacture–United States
	Overall diameter		Country of manufacture–Switzerland
	Single sterile barrier system with protective packaging outside		Unique Device Identifier
	Use-by date		Catalogue number
	Diopter		Right eye
	Date		Left eye
	Caution		Serial number
	Contains biological material of animal origin		Spherical power
	U.S. (Federal) law restricts this device to sale by or on the order of a physician		Cylindrical power
	Store at room/ambient temperature. Do not freeze. Do not expose to temperature greater than 40°C		Axis
	Health care center or Doctor		Spherical equivalent power
	Sterilized using steam		Consult electronic instructions for use

edfu.staar.com
+1-800-352-7842
+41 32 332 8888