# ICL/ TICL Calculator/ IOD Software

eDFU-00481

## TABLE OF CONTENTS

|                      |    |                       |     |
|----------------------|----|-----------------------|-----|
| EN (English).....    | 1  | IT (Italiano).....    | 80  |
| BG (Български) ..... | 7  | LT (Lietuvių).....    | 86  |
| CS (Česky) .....     | 13 | LV (Latviski) .....   | 92  |
| DA (Dansk).....      | 19 | NL (Nederlands).....  | 98  |
| DE (Deutsch).....    | 25 | NO (Norsk).....       | 104 |
| EL (Ελληνικά) .....  | 31 | PL (Polski).....      | 110 |
| ES (Español).....    | 38 | PT (Português) .....  | 116 |
| ET (Eesti).....      | 44 | RO (Română).....      | 122 |
| FI (Suomi) .....     | 50 | SK (Slovenčina) ..... | 128 |
| FR (Français).....   | 56 | SL (Slovenščina)..... | 134 |
| HR (Hrvatski) .....  | 62 | SV (Svenska) .....    | 140 |
| HU (Magyar).....     | 68 | TR (Türkçe) .....     | 146 |
| IS (Íslenska).....   | 74 |                       |     |

EN

## ICL/TICL Calculator /IOD Software

FOR USE WITH STAAR IMPLANTABLE COLLAGEN LENS (ICL AND TICL)

## DIRECTIONS FOR USE

### PRODUCT INFORMATION

Please review this product information completely before performing your initial clinical procedure. All physicians must complete the STAAR Surgical ICL Physician Certification Program.

### DEVICE DESCRIPTION

#### ICL/TICL Calculator /IOD Software

The ICL/TICL Calculator /IOD Software consists of both ICL/TICL Calculator Software (Calculator) and Toric ICL Implantation Orientation Diagram Software (IOD Software). The Calculator and IOD Software resides on STAAR Surgical's e-commerce website in STAAR ICL planning support (Stella™): <https://stella.staar.com/>. Prior to implantation of the Implantable Collamer™ Lens (ICL) physicians use the Online Calculator Software as an aid in the calculation of the size and diopter power (with residual refraction) for physician selection of the lens. For Toric ICLs, an Implantation Orientation Diagram (IOD) is also generated to provide the physician with pictorial representation of axis rotation and alignment.

#### Intended Purpose

The ICL/TICL Calculator /IOD Software is designed to automate the calculation of the ICL/TICL lens power and size based on specific patient biometrics.

#### Intended Patient Population

The ICL/TICL Calculator /IOD Software does not directly interact with patients but is a tool for ophthalmic physicians and staff to calculate ICL/TICL lens power and size based on specific patient biometrics and calculate rotational positioning of the ICL/TICL in the eye.

#### Intended Use Environment

The Calculator /IOD software is a web-based program located on secure servers. The software is available to physicians and may be accessed on a Personal Computer (PC) with Windows 10 or higher with internet access and a secure browser. It is recommended that the user maintain an up-to-date operating system and a secure browser with up-to-date anti-virus software. The calculator may not function properly if used outside of the intended use environment.

Use of the ICL/TICL Calculator /IOD Software is restricted to ICL certified ophthalmic physicians, other refractive practice staff, and authorized intermediaries who have been trained to use the software. Each user of the software is granted access, authorized, and authenticated by STAAR Surgical using a unique ID and password as login credentials.

The STAAR authorized user (surgeon or other healthcare personnel) is responsible of proper management of the ID and password. The authorized user is also responsible for maintenance of PC security. Report any ID, password, or PC security breaches to STAAR Surgical immediately via the phone numbers provided in the reporting section. If there is a problem during data entry or calculation, the user can close the browser and exit the website and all information will be deleted. The user should reload the website and enter the biometric data again to begin calculations.

#### Software Version

ICL/TICL Calculator/IOD Software: Version 7.00

### INDICATIONS FOR USE

The ICL/TICL Calculator Software is designed as an aid in the calculation of the size and diopter power (with residual refraction) for physician selection of the lens. For Toric ICLs, an IOD is also produced to provide the physician with pictorial representation of the axis of rotation and alignment.

### CONTRAINDICATIONS

The use of ICL/TICL Calculator /IOD Software for any other phakic intraocular implants has not been tested or approved for use by STAAR Surgical Company.

### OPERATION

#### Accessibility

The ICL/TICL Calculator /IOD Software resides on STAAR Surgical's e-commerce website in the ICL planning support (Stella™): <https://stella.staar.com/>. Use of the ICL/TICL Calculator /IOD Software is restricted to ICL certified ophthalmic physicians and other refractive practice staff and intermediaries who have been authorized to use the software and have been granted access by STAAR Surgical.

## New Calculations

The calculator can be accessed from multiple entry points in Stella™.

+ Add new patient

Option 1: In the Dashboard, Patients or Surgeons tabs, click on **+Add New Patient** button to create a new patient and enter the patient information, then follow the prompt to access the calculator.

## PRE-OP DATA / CALCULATIONS

+ Add pre-op data

2024-06-25  
10:00:34

Dist/Surg Stella (75659) ^

- OD** OD Pre-Op Saved  
Calculation (MANIFEST) v7.00 OUS ✓
- OS** OS Pre-Op Saved  
Calculation (MANIFEST) v7.00 OUS ✓

## PRE-OP DATA / CALCULATIONS

+ Add pre-op data

2024-06-25  
10:00:34

Dist/Surg Stella (75659) ^

- OD** OD Pre-Op Saved  
Calculation (MANIFEST) v7.00 OUS ✓
- OS** OS Pre-Op Saved  
Calculation (MANIFEST) v7.00 OUS ✓

Option 2: Under the Dashboard, Patients and Reservations tabs, click on an existing patient under Patient ID, then under the PRE-OP DATA / CALCULATIONS select saved calculations by clicking on "Pre-Op Saved" or the associated date to reach the calculator.

## PRE-OP DATA / CALCULATIONS

+ Add pre-op data

2024-06-25  
10:00:34

Dist/Surg Stella (75659) ^

- OD** OD Pre-Op Saved  
Calculation (MANIFEST) v7.00 OUS ✓
- OS** OS Pre-Op Saved  
Calculation (MANIFEST) v7.00 OUS ✓

Option 3: Click on the "Copy this pre-op data to a new set" icon to copy the pre-op data and enter the calculator.

The screenshot shows the ICL/TICL Calculator interface for patient SP12345. The top navigation bar includes links for Dashboard, Patients, Orders, Reservations, Surgeons, Admin, and a Surgeon-specific dropdown for Dist/Surg Stella (75659). The main area is divided into two sections: OD (left) and OS (right). Each section has input fields for Refraction (MANIFEST CL OVER REFRACTION), Biometry, and Calculate For (ICL, TICL, EDOF, MONOCULAR). Below these are detailed eye diagrams for each eye, showing corneal topography, central corneal thickness, anterior chamber depth, and axis measurements (K1, K2, CYL). The OD diagram is labeled 'S' at the top and shows 'White to White' markers. The OS diagram is labeled 'N' at the top and shows 'EDIT OS' and 'White to White' markers. A legend at the bottom left identifies the axis colors: green for K1 Axis, blue for K2 Axis, and white for CYL Axis (refractive).

Figure 1: The ICL/TICL Calculator in Stella™

The user must choose to calculate for an **ICL** or a **Toric ICL**. If **ICL** is selected, EDOF certified users can also select for an **EDOF ICL**. The user then enters the Patient ID, Patient Name (optional), Date of Birth, gender, surgeon, operative eye (OD or OS), and pre-operative data. After entering and saving the pre-operative data, the user should select the **Calculate** button to calculate the residual refraction.

Error messages will appear in red if no value has been entered in a required field, or if the value entered is invalid/outside the required range. The calculation cannot proceed if the errors are not corrected.

Pop-up messages notifying the user of unexpected values will appear for reasons including:

- STAAR Surgical does not manufacture a lens that would treat the refractive error for the patient information entered.
- The data entered includes a cylindrical power, but the user selected the calculation for a lens without cylindrical power.
- The refractive cylinder values and corneal cylinder values do not agree.

Please ensure that the input data is correct.

The user will be presented with a list of lens powers and the expected residual refraction for the patient for each of those lenses based upon the data entered.

### Alternative Lens Length Selection

Users have the ability to select an alternative length from the length recommended by the software at the time of calculation. The user should enter the pre-operative data, select the lens that will achieve the desired outcome and then click on the **Different Length Selected** toggle switch to access the drop-down list of available lengths and click on the length desired. Click **Save** to continue. Alternative lengths should be chosen only after careful consideration by the surgeon.

To select a different length, click the toggle switch for "Different Length Selected" and select desired length

| Sel Sphere | Sel Cylinder | Exp Sphere | Exp Cylinder | Exp Axis |
|------------|--------------|------------|--------------|----------|
| -06.50     | +2.00        | +00.40     | +00.20       | 070      |
| -06.00     | +2.00        | +00.01     | +00.21       | 070      |
| -05.50     | +2.00        | -00.38     | +00.21       | 070      |
| -05.00     | +2.00        | -00.77     | +00.22       | 070      |
| -04.50     | +2.00        | -01.17     | +00.22       | 070      |

**Target Lens**

Different Length Selected    12.6 mm ▾  
12.1 mm  
12.6 mm  
13.2 mm  
13.7 mm

Figure 2: Selecting Lens Length

### Saving and Printing Calculations

The pre-operative data and calculation results will be saved to the STAAR server and can be printed by the user by selecting **Print Report** button on the top navigation bar. After a Toric lens has been ordered, an IOD can be generated. See the IOD section below.

### IOD (Implantation Orientation Diagram) – For Toric Lenses Only

Once a Toric lens has been ordered or reserved, the user can generate an IOD. The IOD is produced to provide the physician with pictorial representation of the axis of rotation and alignment. The user can reach the IOD screen in several ways:

1. Under Patients tab, click on the **IOD** icon next to the target lens
2. Under Dashboard or Patients, click on the Patient ID, then click on ORDERS, then click on the **IOD** icon under Actions
3. Under Orders, click on the order Number to access the Order Details page, click on the **IOD** icon
4. On the Order Confirmation page, click on the **IOD** icon

Dashboard Patients Orders Reservations Surgeons Admin Checkout - 0 ? + Add new patient

## Patients

Patient Activity  
You are currently viewing all activity for all currently active patients  
Patients are sorted by most recent activity at the top.

| Activity<br>YYYY-MM-DD | Patient ID<br>Name (Last, First)       | DOB<br>YYYY-MM-DD | Surgeon<br>Last, First   | OD  | OS  |
|------------------------|--|-------------------|--------------------------|---|---|
| 2024-06-25             | <a href="#">BTS007</a>                 | 1999-01-01        | 75659: Stella, Dist/Surg | Calculations<br> Lens OD<br>Toric Myopic 12.1 mm<br>$-8.50 +1.50 \times 060$ | Calculations<br> Lens OS<br>Toric Myopic 12.1 mm<br>$-8.50 +1.50 \times 052$ |
| 2024-06-25             | <a href="#">SP12345</a>                | 2000-01-01        | 75659: Stella, Dist/Surg | Target lens pending   | Target lens pending   |
| 2024-06-25             | <a href="#">T001</a><br>Test           | 1996-12-12        | 75659: Stella, Dist/Surg | Myopic 12.6 mm<br>$-8.00$   | Myopic 12.6 mm<br>$-8.00$   |
| 2024-06-25             | <a href="#">01234</a><br>Sample, Belle | 1999-01-01        | 75659: Stella, Dist/Surg | Toric Myopic 12.6 mm<br>$-6.50 +1.50 \times 080$  | Toric Myopic 12.6 mm<br>$-6.50 +1.50 \times 080$  |
| 2024-06-13             | <a href="#">DS001</a><br>Testt, Meghan | 1993-08-14        | 75659: Stella, Dist/Surg | Toric Myopic 12.6 mm<br>$-6.50 +1.00 \times 090$  | Toric Myopic 12.6 mm<br>$-6.50 +1.00 \times 090$  |

Click on the IOD button to access the IOD preview window

Figure 3: Accessing Implantation Orientation Diagram (IOD) from the Patients tab

LENS / SURGERY INFO | RESERVATIONS | ORDERS

**Orders** You are currently viewing all orders for this patient.

| Date       | Number                 | Target Lens  | Ordered Lens   | Status  | Actions   |
|------------|------------------------|--|--|---------|---|
| 2024-06-25 | <a href="#">235941</a> |  Toric Myopic 12.6 mm<br>$-5.50 +2.00 \times 070$ | <b>VTICM5_12.6</b><br>$-5.50 +2.00 \times 079$ #T1666358 | Ordered |  |


Click on the IOD button to access the IOD preview window

Figure 4: Accessing Implantation Orientation Diagram (IOD) from the Patient Details page

**Items in this order**

---

Patient: [SP12345](#)

|   |  |   |
|---|--|---|
| <br>Primary - OD | <b>EVO+ Visian ICL</b><br>\$0.00<br>T1666358 VTICM5_12.6<br>$-5.50 +2.00 \times 079$<br>Calculation: <a href="#">SP12345 - OD</a><br>Surgeon: <a href="#">75659 - Dist/Surg Stella</a> |  |
|---|--|---|

Click on the IOD button to access the IOD preview window

Figure 5: Accessing Implantation Orientation Diagram (IOD) from the Order Details and Order Confirmation pages

Once the IOD preview window opens, review the information on the screen is correct, then click on Generate IOD.

The users have the option to open the IOD into full screen by clicking the **Fullscreen** button, or print the IOD by clicking the **Print** button.

# Implantation Orientation Diagram (IOD)

SP12345

(F)

2000-01-01  
24 yrs

SURGEON

Dist/Surg Stella

TARGET LENS

Toric Myopic 12.6 mm

Please confirm that Patient and Lens information is correct. Then press Generate IOD.

Lens Ordered

| Serial Number | Model       | Version         | Exp. SEQ | Exp. Ref         | Sphere | Cylinder | Axis |
|---------------|-------------|-----------------|----------|------------------|--------|----------|------|
| T1666358      | VTICM5_12.6 | EVO+ Visian ICL | -0.27    | -0.38 +0.21 x070 | -5.50  | +2.00    | 079  |

Override

Click on the generate IOD button to create the diagram

Cancel

Generate IOD

Figure 6: Generating the IOD

**Stella.**

SP12345 F

Patient ID

Birthdate: 2000-01-01 (24 years) Surgeon: Dist/Surg Stella (75659)  
Customer: [REDACTED] Surgery Date: —

**OD** Target Lens: Toric Myopic 12.6 mm -5.50 +2.00 x070

**Lens Ordered**  
Please confirm serial number.

| Serial Number | Model       | Version         | Exp. SEQ | Exp. Ref         | Sphere | Cylinder | Axis |
|---------------|-------------|-----------------|----------|------------------|--------|----------|------|
| T1666358      | VTICM5_12.6 | EVO+ Visian ICL | -0.27    | -0.38 +0.21 x070 | -5.50  | +2.00    | 079  |

PLACE LENS LABEL HERE

Rotate lens **clockwise 9°** after horizontal implantation to position lens at **171°**.  
No PIs are required with this lens.

**Implantation Orientation Diagram**

**IMPLANTATION ORIENTATION DIAGRAM**

Show Mendez ring  
 Show eye outline

**Fullscreen**

**Print**

**Cancel**

Information on which axis to position the lens and if peripheral iridotomy is required

**STAAR SURGICAL**  
2024-06-25 12:33 | [View Patient Report](#)

**EVO ICL™**  
Calculated with 700 OUS

Figure 7: Implantation Orientation Diagram (IOD)

## REPORTING

Adverse events and/or potentially sight-threatening complications that may reasonably be regarded as ICL/TICL Calculator /IOD Software related and that were not previously expected in nature, severity and degree of incidence should be reported to STAAR Surgical. For surgeons/patients located in the EU, the competent authority should also be notified in the EU member state where the surgeons/patients are established.

International Phone: +(41) 32 332 88 88

USA/Canada Phone: +1 (800) 352-7842

## SYMBOL GLOSSARY

Medical device

Manufacturer

Authorized representative in the European Community

Unique Device Identifier

CE conformity marking per European Council Directive 93/42/EEC or European Council Regulation (EU) 2017/745

Consult electronic instructions for use

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The screenshot shows the Stella software interface for ICL/TICL calculations. It is divided into two main sections: OD (Left) and OS (Right).  
**Top Section:**

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  
**Middle Section:**

- OD Eye View:** Shows corneal topography with axis markings (K1 Axis, K2 Axis, CYL Axis (refractive)) and a central thickness measurement of 12.00 mm.
- OS Eye View:** Shows corneal topography with axis markings (K1 Axis, K2 Axis, CYL Axis (refractive)) and a central thickness measurement of 12.00 mm.

  
**Bottom Section:**

- OD Eye View:** Shows corneal topography with axis markings (K1 Axis, K2 Axis, CYL Axis (refractive)) and a central thickness measurement of 12.00 mm.
- OS Eye View:** Shows corneal topography with axis markings (K1 Axis, K2 Axis, CYL Axis (refractive)) and a central thickness measurement of 12.00 mm.

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The screenshot displays the Stella software's ICL/TICL calculator. At the top, it shows a navigation bar with links to Instrumentpanel, Patienter, Bestillinger, Reservationer, Kirurger, Admin, and a patient ID (SP12345). Below the navigation is a header with date (2024-06-25), time (10:00:34), and location (Dist/Surg Stella, 75659). The main area is divided into two sections: **OD** (left) and **OS** (right). Each section contains:

- REFRAKTION: MANIFEST | CL-OVERREFRAKTION:** Includes fields for D, CYLINDER, AKSE, and BVD, with a value of 12.00 mm highlighted.
- BEREGN FOR:** Options for ICL, TICL, and EDOF (selected) or MONOKULÆR.
- BIOMETRI:** Displays K1, K2, ACD, CCT, and WW measurements.
- Eye Diagrams:** Detailed diagrams of the eye showing corneal curvature (S, N, T, S), central corneal thickness, and anterior chamber depth (For kammer dybde).
- Axis Markings:** Shows K1 Axis, K2 Axis, and CYL Axis (refractive) with corresponding angle markers (90°, 180°, 0°, 270°).

In the OS section, there is a red box around the text "REDIGER OS" over the eye diagram.

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The screenshot shows the STAAR surgical software interface for calculating lens parameters. On the left, there are sections for OD (Left Eye) and OS (Right Eye). Each section displays manifest refraction (-5.00 D, +2.00 D, 70°), biometry (44.00 D, 45.00 D), and calculated parameters (12.00 mm BVD). Below these are tables for OD and OS showing various lens options. A green box highlights the 'Beregning afsluttet' (Calculation completed) message and the instruction 'Vælg en linse fra listen nedenfor' (Select a lens from the list below). An orange callout box points to the 'Anden længde valgt' (Other length selected) button in the OD section, which is checked. A dropdown menu shows lens lengths: 12.6 mm, 12.1 mm, 12.6 mm (selected), 13.2 mm, and 13.7 mm. The right side of the screen shows a table of calculated values for different lens lengths, with the first row highlighted. At the bottom, there are 'Annuler' (Cancel) and 'Gem' (Save) buttons.

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The screenshot shows the STAAR Surgical software interface for calculating lens parameters. At the top, patient information is displayed: Name: NACHNAME, VORNAME, Date of Birth: 2000-01-01, Age: 24 Jahre, Location: CHIRURG, ID: 75659: Dist/Surg Stalla. Below this, the patient ID SP12345 is shown. The interface is split into two main sections: OD (Left) and OS (Right). Each section shows refractive parameters (Manifest refraction: -5.00 D, +2.00 D, 70°, 12.00 mm; BIOMETRY: 44.00 D, 70°, K1: K1-ACHSE, K2: K2-ACHSE; 45.00 D, 160°, ACD: 3.40 mm, CCT: 500 µm, WW: 11.40 mm), and calculation fields for ICL (OD) and KL (OS). A green box on the left indicates 'Berechnung abgeschlossen.' (Calculation completed). A callout box with orange text and border points to the 'Abweichende Länge gewählt' (Off-axis length selected) button and the dropdown menu for selecting an alternative lens length. The dropdown menu lists options: 12.6 mm, 12.1 mm, 12.6 mm (selected), 13.2 mm, and 13.7 mm. To the right, a list of lens options is shown, with the first item highlighted: Erw. Sphäre: +00.56, Erw. Zylinder: +00.00, Erw. Achse: 090, Erw. SEQ: +00.56. The status bar at the bottom right shows the date: 2024-06-25 and the message: 'geschlossen. Bitte wählen Sie die gewünschte Linse aus der untenstehenden Liste. Berechnet mit 7.00 OUS'.

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The main area is divided into two sections: **OD** (Left) and **OS** (Right). Each section has fields for REFRACTION (MANIFESTAR), BIOMETRÍA, and CALCULAR PARA (ICL, TICL, EDOF, MONOCULAR).

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**Eye Diagrams:** Large circular diagrams for OD and OS show the eye's anatomy with labels for S (Superior), N (Nasal), T (Temporal), and I (Inferior). The OD diagram includes a legend for K1 Axis (green dashed line), K2 Axis (blue dashed line), and CYL Axis (refractive) (white line). The OS diagram has an 'EDITAR OS' button overlaid on the eye diagram.

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The screenshot shows the Stella software interface for calculating ICL/TICL. It displays two main sections for the right eye (OD) and left eye (OS). Each section includes a summary table for refractive data and a detailed eye diagram. The eye diagrams show the cornea, lens, and retina, with various measurements labeled such as K1 and K2 axes, and CYL axis. The software interface also includes navigation menus at the top and bottom, and a search bar.

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The screenshot shows the software interface for surgical planning. It displays two sets of data for the left eye (OD) and right eye (OS). Each set includes manifest refraction, cylinder, axis, and biometric parameters. Below this, there are tables for selecting alternative lens thicknesses. A callout box highlights the 'Elterő hossz kiválasztva' (Alternative thickness selected) button and a dropdown menu showing options like 12.6 mm, 12.1 mm, 12.6 mm, 13.2 mm, and 13.7 mm. Another callout box points to the 'Másik hossz kiválasztásához kattintson a „Elterő hossz kiválasztva” kapcsolóra, majd válassza ki a kívánt hosszúságot' (Click the switch for selecting alternative thickness, then choose the desired length) instruction. At the bottom, there are buttons for 'Mégse' (Cancel) and 'Mentés' (Save).

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The screenshot shows the Stella software interface with two main sections for eyes OD and OS. Each section includes a lens diagram, a table of lens parameters, and a detailed eye diagram with axis measurements.

**OD Section:**

- Lens Parameters:**
  - LJÓSBROT: BIRTA | CL YFIR LJÓSBROT
  - D, D, ÅS, BVD
  - KÜLA, SÍVALNINGUR, ÅS, BVD
  - 12.00 mm
- Reikna fyrir:**
  - ICL, TOSIC
  - EDOF, EINS/EI
- Lífmælingar:**

|   |   |    |    |     |     |    |
|---|---|----|----|-----|-----|----|
| D | ° | K1 | ÅS | mm  | μm  | mm |
| D | ° | K2 | ÅS | ACD | CCT | HH |
- Eye Diagram:** Shows a circular eye with axis markings: S (90°), T (180°), N (0°), and E (270°). It also shows K1 Axis, K2 Axis, and CYL Axis (refractive).

**OS Section:**

- Lens Parameters:**
  - LJÓSBROT: BIRTA | CL YFIR LJÓSBROT
  - D, D, ÅS, BVD
  - KÜLA, SÍVALNINGUR, ÅS, BVD
  - 12.00 mm
- Reikna fyrir:**
  - ICL, TOSIC
  - EDOF, EINS/EI
- Lífmælingar:**

|   |   |    |    |     |     |    |
|---|---|----|----|-----|-----|----|
| D | ° | K1 | ÅS | mm  | μm  | mm |
| D | ° | K2 | ÅS | ACD | CCT | HH |
- Eye Diagram:** Shows a circular eye with axis markings: S (90°), T (180°), N (0°), and E (270°). It also shows K1 Axis, K2 Axis, and CYL Axis (refractive).

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The screenshot displays the Stella™ software interface for ICL/TICL calculations. At the top, it shows the date (2024-06-25), time (10:00:34), and patient information (KIRURGS 75659: Dist/Surg Stella). Below this, there are two main sections for 'OD' (Left Eye) and 'OS' (Right Eye).

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The screenshot shows the software interface for selecting alternative lens thickness. On the left, there are two columns for OD and OS eyes. Each column has a section for 'REFRAKCIJA: IZTEIKTA' (Programmed) and 'KAM APRĒKINĀT' (Calculate). The OD section shows values: -5.00 D, +2.00 D, 70°, 12.00 mm. The OS section shows values: -4.00 D, 0.00 D, 90°, 12.00 mm. Below these are 'BIOMETRIJA' tables. A large orange callout box highlights the 'Izvēlēts cits garums' (Selected alternative thickness) option in the dropdown menu for the OD eye, which is set to 12.6 mm. The dropdown also lists 12.1 mm, 12.5 mm, 13.2 mm, and 13.7 mm. An orange arrow points from the callout to the dropdown menu. At the bottom right, there is a note in Latvian: 'Lai atlasiitu citu garumu, noklikšķiniet uz nolaižamās izvēlnes, lai atlasiit opciju "Izvēlēts cits garums" un atlasiit vēlamo garumu'. The date '2024-06-25' and the note 'Lūdzu, atlasiit mērķa lēcu turpmāk esošajā sarakstā.' are also visible.

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The main area is divided into two sections: OD (Left) and OS (Right). Each section has a 'REFRACTIE: TONEN' field with input fields for D, CILINDER, AS, and AVA, and a 'BIOMETRIE' field with input fields for D, K1, K2, K1-AS, K2-AS, ACD, CHD, and WW. There are also 'UITREKKEN' and 'BEREKENEN VOOR' buttons.

Below these fields are detailed eye diagrams for OD and OS. The OD diagram shows a circular eye with axis markings S (90°), T (180°), N (0°), and E (270°). It includes labels for 'Centrale hoornvlies dikte', 'Diepte voorste oogkamer', and 'White-to-white'. The OS diagram is similar but labeled 'BEWERKEN OS'.

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A screenshot of the STAAR Surgical software interface. The top navigation bar includes 'Painel de controlo', 'Doentes', 'Encomendas', 'Reservas', 'Cirurgiões', 'Admin', 'Imprimir Relatório', 'Checkout - 0', and a user icon. The main area shows patient information: ID: SP12345, Name: F, Date: 2000-01-01, Age: 24 years, Surgeon: 75659: Dist/Surg Stell. There are sections for OD and OS eyes, each with refraction details and biometry tables. A callout box highlights the 'Distância Focal Diferente Selecionada' button (green checkmark) and the dropdown menu showing options: 12.6 mm, 12.1 mm, 12.5 mm, 13.2 mm, and 13.7 mm. An arrow points from the text in the callout box to the '12.6 mm' option in the dropdown.

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The screenshot shows the Stella software interface. At the top, there's a navigation bar with icons for Nadzorna plošča, Pacienti, Naročila, Rezervacije, Kirurgi, and Administrator. Below that, it shows ID PACIENTA (SP12345), DOB (2000-01-01), Age (24 leta), and KIRURG (75659: Dist/Surg Stella). On the right, there are buttons for Tiskanje poročila, Blokajna -0, and Izračunaj. The main area is divided into two sections: 'OD' and 'OS'. Each section has a 'REFRAKCIJA' dropdown set to 'SUBJEKTIVNA', a 'BIOMETRIJA' section with various measurements like D, CILINDER, OS, BVD, and CCT, and an 'IZRAČUNAJTE ZA' section with options for ICL, TORIC, EDOF, and MONOKULARNA. Below these are detailed eye diagrams for 'OD' and 'OS' eyes, showing corneal thickness, pupil position, and axis measurements (K1, K2, CYL). A large 'UREDI OS' button is visible in the OS section.

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The screenshot shows the software interface for calculating lens parameters. At the top, it displays patient information: ID PACIENTA (NAZIV: IME, PRIMEK) 2000-01-01 24 leta KIRURG 75659: Dist/Surg Stola. Below this, the patient ID SP12345 is shown. The interface is divided into sections for OD (Left Eye) and OS (Right Eye). Each section has fields for Refrakcija (Subjective and Objective), Biometrija (Biometry), and a summary table. A green bar at the bottom of each section indicates the calculation is complete ('Izračun je dokončan'). A callout box with orange border and text 'Če želite izbrati drugo dolžino, kliknite preklopno stikalo za »Izbrana drugačna dolžina« in izberite želeno dolžino.' points to the 'Izbrana drugačna dolžina' button in the dropdown menu for the OS section. The dropdown menu shows several options: 12.6 mm, 12.1 mm, 12.6 mm (selected), 13.2 mm, and 13.7 mm. The right side of the interface shows a table of calculated values for both eyes, with the date 2024-06-25 and time 7:00 OUS.

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© STAAR SURGICAL
STAAR Surgical AG